
Kit Cãnula Discografia d.View

Reg. ANVISA: 80195520032.

BRAMSYS

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USE INSTRUCTIONS

D.VIEW DISCOGRAPHY CANNULA KIT

The d.View Discography Cannula Kit is indicated for diagnostic procedures for discogenic origin pain, including provocative discography and/or disc analgesia, in order to administer fluids into the intervertebral disc and control its pressure during the procedure.

INDICATIONS FOR USE

The guide cannula leads the introducer cannula to the center of the possibly damaged disc nucleus, and, with the assistance of the inflator, the initial disc pressure is measured.

Subsequently, the contrast is delivered from the Inflator to the disc nucleus through the introducer cannula. During this process, the physician monitors the injected volume and disc pressure using the inflator.

COMPOSITION/COMMERCIAL PRESENTATION

Primary packaging: blister or surgical grade paper envelope containing 2 or 6 cannulas wrapped in a plastic protective tube.

Secondary packaging: blister or surgical grade paper envelope containing ALL the other units that make up the kit, as presented below.

Tertiary packaging: individual packaging box.

Quaternary packaging: transport box.

Code	Guide Cannula			Introducer Cannula		
	Ø	Length	Qty.	Ø	Length	Qty.
DVS-1822-07	18G	70mm	1	22G	150mm	1
DVS-1822	18G	100mm	1	22G	150mm	1
DVS-1822-20	18G	100mm	1	22G	200mm	1
DVS-1822-1520	18G	150mm	1	22G	200mm	1
DVS-2025-07	20G	70mm	1	25G	150mm	1
DVS-2025	20G	100mm	1	25G	150mm	1
DVS-2025-20	20G	100mm	1	25G	200mm	1
DVS-2025-1520	20G	150mm	1	25G	200mm	1
DVS-1822-07/3	18G	70mm	3	22G	150mm	3
DVS-1822/3	18G	100mm	3	22G	150mm	3

DVS-1822-20/3	18G	100mm	3	22G	200mm	3
DVS-1822-1520/3	18G	150mm	3	22G	200mm	3
DVS-2025-07/3	20G	70mm	3	25G	150mm	3
DVS-2025/3	20G	100mm	3	25G	150mm	3
DVS-2025-20/3	20G	100mm	3	25G	200mm	3
DVS-2025-1520/3	20G	150mm	3	25G	200mm	3

INSTRUCTIONS FOR USE

1. Ensure that the packaging is not damaged.
2. Remove the cannulas, the inflation device, the three-way adapter, the surgical planning scale, and the surgical marker from the sterile packaging in a clean environment using aseptic techniques.
3. Verify that the cannulas, the inflation device, the three-way adapter, the surgical planning scale, and the surgical marker are intact.
4. Position the patient according to the surgical procedure to be performed.
5. If necessary, use the surgical planning scale and the surgical marker to define the puncture site.
6. Perform the percutaneous insertion of the guide cannula into the patient with the assistance of fluoroscopy or X-ray.
7. Remove the mandrel from the guide cannula.
8. Insert the introducer cannula into the guide cannula and position it at the desired target with the assistance of fluoroscopy or X-ray.
9. Connect a needle to the luer lock inlet of the extension tubing of the inflation device.
10. Aspirate the radiopaque/anesthetic/medication solution to be injected.
11. Disconnect the needle from the extension tubing of the inflation device and discard the needle.
12. Remove the mandrel from the introducer cannula.
13. Attach the extension tubing of the inflation device to the introducer cannula.
14. Inject the radiopaque/anesthetic/medication solution.
15. Disconnect the extension tubing of the inflation device from the introducer cannula.
16. Remove the introducer cannula.
17. Remove the guide cannula.
18. After completing the procedure, the material should be destroyed and disposed of in an appropriate location, following the applicable regulations for the disposal of sharp and contaminated materials adopted by the hospital.

CONTRAINDICATIONS

- Infection or hematoma near the chosen puncture site
- Pregnancy;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Changes in coagulation, for example during anticoagulant therapy are a relative contraindication;
- Failure to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast, antibiotics and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical or psychological condition that increases the risk of reduced performance of the procedure.

STORAGE AND TRANSPORT CONDITIONS

To maintain the product's integrity, it must:

Be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging should not be perforated or torn and should be kept away from sharp or pointed instruments.

Do not place under weight or volume, as they could damage its technical characteristics.

Do not place near or with different categories, as they could damage its technical characteristics.

They should be stored and transported following good practices for the transport of sterilized products.

Transportation of the product should be done in its original packaging. Storage and transportation temperature: 5°C to 60°C.

Storage and transportation humidity: 10% to 80%.

OPERATING CONDITIONS

The user, before using the product, must carefully read the instructions for use that will assist in its correct and safe use, in addition to protecting against possible risks.

Do not use the product past its expiration date.

Ensure that the sterilization indication is correct.

The product should be used by qualified medical professionals with the necessary knowledge and experience for the proper use of the product and handling of sterile products.

Verify the integrity of the packaging before using the product; in case of violation, DO NOT USE it, as the guarantee of its sterility will be compromised.

Temperature: 10°C to 40°C.

Humidity: 20% to 90%.

PRECAUTIONS, RESTRICTIONS, STORAGE AND WARNINGS

- The correct and safe use of this product for its intended purpose is governed by the conditions described in this user manual.
- The product must not be reused and/or reprocessed, as it is for SINGLE USE only and must be destroyed and discarded as MEDICAL WASTE.
- This product is sterilized using Ethylene Oxide (ETO).
- No modifications to this product are allowed.
- Disposal must comply with the requirements of RDC No. 222, dated March 28, 2018, and other applicable laws and regulations.

No. Batch Number, Sterilization Date and Expiration Date: see product packaging.

ANVISA: 80195520032

Technical Manager: Marcos Rojas | **CREA:** 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance under number UY44WX94W79M, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br