
# Kit Cãnula All One Block

Reg. ANVISA: 80195520034.

## **BRAMSYS**

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## **USE INSTRUCTIONS**

### **PERCUTANEOUS DISK ACCESS CANNULA KIT L**

#### **PURPOSE**

The Percutaneous Disc Access Cannula Kit aims to access the disc with the minimally invasive percutaneous technique. In this way, the surgeon can decompress the protruding disc.

#### **INDICATIONS OF USE**

The Percutaneous Disc Access Cannula Kit is indicated for accessing the disc in the spine in minimally invasive surgical procedures for decompression of intervertebral disc herniation, improving functional capacity and pain caused by compression.

#### **COMPOSITION**

Primary packaging: blister or surgical grade paper envelope, containing: (1) guide wire 18G, 330mm, (1) guide wire 20G, 330mm, (1) introducer cannula 16G, 200mm, (1) trocar 10G, 165mm, (1) faceted cannula 10G, 165mm, (1) trephine 12G, 195mm, (1) obturator 12G, 195mm, (1) depth stop, (1) surgical planning scale, (1) surgical marker.

Secondary packaging: individual box.

Tertiary packaging: transport box.

### **USE INSTRUCTIONS**

- Define puncture site for access using fluoroscopy equipment.
- Insert introducer cannula to the desired target.
- Remove mandrel from the introducer cannula.
- Inject contrast to visualize the disc.
- Insert guide wire to the desired target.
- Carefully remove the introducer cannula, keeping the guide wire in the desired target.
- Select trocar or faceted cannula to perform the procedure.
- Fit obturator in the selected option.
- Position the depth stop on the selected option
- Introduce the set of the selected option, obturator, and depth stop into the patient through the guide wire.
- Position the depth stop on the selected option until it touches the patient's skin.
- Lock depth stop to stabilize the selected option on the patient.
- Remove obturator from selected option.
- Introduce trephine in the selected option.
- Press trephine against the disc wall making rotating movements until the disc is penetrated.
- Remove trephine from selected option.
- Perform the procedure on the disk through the selected option.
- After performing the procedure, remove instruments from the selected option.
- Remove selected option from patient.
- After completing the procedure, the material must be destroyed and disposed of in an appropriate place, following the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.

### **CONTRAINDICATIONS**

- Infection or hematoma near the site chosen for the puncture;
- Pregnancy;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Changes in coagulation, eg during anticoagulant therapy are a relative contraindication;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast, antibiotics and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical or psychological condition that increases the risk of reduced performance of the procedure.

### **STORAGE CONDITIONS**

The Percutaneous Disc Access Cannula Kit must always be stored protected against all types of humidity, dirt and contact with flammable, corrosive, toxic or explosive substances. The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.

Do not place it under weight or volumes, as they could damage its technical characteristics and compromise its integrity.

Do not store the product where it can be exposed to direct sunlight, high temperatures, humidity or dust as this will compromise its integrity/quality.

Do not place it next to or together with different categories, as they could damage its technical characteristics and compromise its integrity/quality.

Storage temperature: 5°C to 60°C.

Storage humidity: 10% to 80%.

### **TRANSPORTATION CONDITIONS**

The Percutaneous Disc Access Cannula Kit must be transported and stored always protected against all types of humidity, dirt and contact with flammable, corrosive, toxic and explosive substances. The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.

**IMPORTANT:** they must be packaged and transported, following the rules of good practices for the transport of sterilized products. The transport of the product must be done in its original packaging, duly packed.

## **PRECAUTIONS, RESTRICTIONS AND WARNINGS**

- The user, before using the product, must carefully read the instructions for use that will help correct and safe use, as well as protecting against possible risks.
- Do not use the product after its expiration date, always check if the sterilization indication is correct.
- The product must only be used by qualified medical professionals, who have the knowledge and experience necessary for the proper use of the product, knowledge of handling sterile products and, mainly, of the importance of maintaining a totally sterile environment in the operating rooms, where care and precaution with every detail are essential so that there is no risk of contamination.
- The product must not be reused and/or reprocessed, being for SINGLE use, and must be destroyed and disposed of in HOSPITAL WASTE.
- This product is sterilized by Ethylene Oxide (ETO).
- The product must be transported and stored, always protected against all types of moisture, dirt and contact with flammable, corrosive, toxic or explosive substances.
- The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.
- Do not store or store this product next to products of different categories.
- Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
- This product is for single use, sterilized by Ethylene Oxide (ETO).

This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

**Lot number, Sterilization and Validity Date:** see the product packing.

**MS Registration:** 80195520034

**Technical Responsible:** Marcos Rojas

**CREA:** 5061214879

## **WARRANTY CERTIFICATE**

**BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

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## **MODELOS / MODELS**

LNBH-1805-05; LNBH-1805-10; LNBH-1810-05; LNBH-1810-10; LNBH-1815-05;  
LNBH-1815-10; LNBH-2005-05; LNBH-2005-10; LNBH-2010-05; LNBH-2010-10; LNBH-  
2015-05; LNBH-2015-10; LNBH-2105-05; LNBH-2105-10; LNBH-2110-05; LNBH-2110-10;  
LNBH-2115-05; LNBH-2115-10; LNBH-2205-05; LNBH-2205-10; LNBH-2210-05; LNBH-  
2210-10; LNBH-2215-05; LNBH-2215-10; LNBH-1805-05; LNBH-1805-10; LNBH-1810-05;  
LNBH-1810-10; LNBH-1815-05; LNBH-1815-10; LNBH-2005-05; LNBH-2005-10; LNBH-  
2010-05; LNBH-2010-10; LNBH-2015-05; LNBH-2015-10; LNBH-2105-05; LNBH-2105-10;  
LNBH-2110-05; LNBH-2110-10; LNBH-2115-05; LNBH-2115-10; LNBH-2205-05; LNBH-  
2205-10; LNBH-2210-05; LNBH-2210-10; LNBH-2215-05; LNBH-2215-10; LNBK-1805-05;  
LNBK-1805-10; LNBK-1810-05; LNBK-1810-10; LNBK-1815-05; LNBK-1815-10; LNBK-  
2005-05; LNBK-2005-10; LNBK-2010-05; LNBK-2010-10; LNBK-2015-05; LNBK-2015-10;  
LNBK-2105-05; LNBK-2105-10; LNBK-2110-05; LNBK-2110-10; LNBK-2115-05; LNBK-  
2115-10; LNBK-2205-05; LNBK-2205-10; LNBK-2210-05; LNBK-2210-10; LNBK-2215-05;  
LNBK-2215-10; LNBT-1805-05; LNBT-1805-10; LNBT-1810-05; LNBT-1810-10; LNBT-  
1815-05; LNBT-1815-10; LNBT-2005-05; LNBT-2005-10; LNBT-2010-05; LNBT-2010-10;  
LNBT-2015-05; LNBT-2015-10; LNBT-2105-05; LNBT-2105-10; LNBT-2110-05; LNBT-  
2110-10; LNBT-2115-05; LNBT-2115-10; LNBT-2205-05; LNBT-2205-10; LNBT-2210-05;  
LNBT-2210-10; LNBT-2215-05; LNBT-2215-10.