
Kit Cãnula CTC

Reg. ANVISA: 80195520036.

BRAMSYS

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USE INSTRUCTIONS

CTC CANNULA KIT

The CTC Cannula Kit is a disposable single-use product designed to isolate and guide the radiofrequency electrode. Its composition includes: cannula with canalized introducer, radiopaque sonovisible anti-adherent insulated surface, thermocouple mandrel, and active tip.

INDICATIONS FOR USE

The CTC Cannula Kit is indicated for the treatment of chronic pain through minimally invasive surgical procedures involving controlled temperature radiofrequency, aiming to alleviate pain and improve functional capacity. The cannula features a radiopaque anti-adherent insulated coating, thermocouple mandrel, and sonovisible active tip, enabling guidance through ultrasound or fluoroscopy.

The CTC Cannula Kit is a disposable single-use product designed to isolate and guide the radiofrequency electrode. Its composition includes: cannula with canalized introducer, radiopaque sonovisible anti-adherent insulated surface, thermocouple mandrel, and active tip.

COMPOSITION/COMMERCIAL PRESENTATION

Primary packaging: Surgical-grade blister or envelope containing 1 or 2 cannulas, wrapped in a plastic protection tube, 1 surgical planning scale, 1 surgical marker.

Secondary packaging: Individual packaging box, containing primary packaging.

Tertiary packaging: Transport box.

Code	Ø	Length	Active tip	Qty Cannulas
CTC-2050-5/1	20G	50mm	5	1
CTC-2050-5/2	20G	50mm	5	2
CTC-21100-10/1	21G	100mm	10	1
CTC-21100-10/2	21G	100mm	10	2
CTC-21100-5/1	21G	100mm	5	1
CTC-21100-5/2	21G	100mm	5	2
CTC-21150-10/1	21G	150mm	10	1
CTC-21150-10/2	21G	150mm	10	2
CTC-21150-5/1	21G	150mm	5	1
CTC-21150-5/2	21G	150mm	5	2
CTC-2150-5/1	21G	50mm	5	1
CTC-2150-5/2	21G	50mm	5	2
CTC-18100-10/1	18G	100mm	10	1
CTC-18100-10/2	18G	100mm	10	2
CTC-18100-5/1	18G	100mm	5	1
CTC-18100-5/2	18G	100mm	5	2
CTC-18150-10/1	18G	150mm	10	1
CTC-18150-10/2	18G	150mm	10	2
CTC-18150-5/1	18G	150mm	5	1
CTC-18150-5/2	18G	150mm	5	2
CTC-1850-5/1	18G	50mm	5	1
CTC-20100-10/1	20G	100mm	10	1
CTC-20100-10/2	20G	100mm	10	2
CTC-20100-5/1	20G	100mm	5	1
CTC-20100-5/2	20G	100mm	5	2
CTC-20150-10/1	20G	150mm	10	1
CTC-20150-10/2	20G	150mm	10	2
CTC-20150-5/1	20G	150mm	5	1
CTC-20150-5/2	20G	150mm	5	2

INSTRUCTIONS FOR USE

- Verify that the packaging is not damaged or wet.
- Ensure that the sterilization indication is correct.
- Check the sterilization expiration date.
- Open the package and remove the packaged and sealed product.
- Following procedures for handling sterile products, open the second package and remove the product.
- Prepare the patient before the procedure according to standard techniques.
- Perform fluoroscopy or ultrasound procedure to define the desired target and puncture location.
- Insert the CTC Cannula up to the desired target, maintaining the thermocouple mandrel.
- Perform the radiofrequency procedure at the defined location.
- After completing the radiofrequency process.
- Remove the CTC Cannula with the thermocouple mandrel.
- At the end of the procedure, the material must be destroyed and disposed of appropriately, following applicable regulations for disposal of sharp and contaminated materials as adopted by the hospital in HOSPITAL WASTE.

CONTRAINDICATIONS

- Infection or hematoma near the chosen puncture site
- Pregnancy;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Changes in coagulation, for example during anticoagulant therapy are a relative contraindication;
- Failure to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast, antibiotics and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical or psychological condition that increases the risk of reduced performance of the procedure.

STORAGE AND TRANSPORT CONDITIONS

To maintain the product's integrity, it must:

Be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging should not be perforated or torn and should be kept away from sharp or pointed instruments.

Do not place under weight or volume, as they could damage its technical characteristics.

Do not place near or with different categories, as they could damage its technical characteristics.

They should be stored and transported following good practices for the transport of sterilized products.

Transportation of the product should be done in its original packaging.

Storage and transportation temperature: 5°C to 60°C.

Storage and transportation humidity: 10% to 80%.

OPERATING CONDITIONS

The user, before using the product, must carefully read the instructions for use that will assist in its correct and safe use, in addition to protecting against possible risks.

Do not use the product past its expiration date.

Ensure that the sterilization indication is correct.

The product should be used by qualified medical professionals with the necessary knowledge and experience for the proper use of the product and handling of sterile products.

Verify the integrity of the packaging before using the product; in case of violation, DO NOT USE it, as the guarantee of its sterility will be compromised.

Temperature: 10°C to 40°C.

Humidity: 20% to 90%.

PRECAUTIONS, RESTRICTIONS, STORAGE AND WARNINGS

- The correct and safe use of this product for its intended purpose is governed by the conditions described in this user manual.
- The product must not be reused and/or reprocessed, as it is for SINGLE USE only and must be destroyed and discarded as MEDICAL WASTE.
- This product is sterilized using Ethylene Oxide (ETO).
- No modifications to this product are allowed.
- Disposal must comply with the requirements of RDC No. 222, dated March 28, 2018, and other applicable laws and regulations.

ANVISA: 80195520036

Technical Manager: Marcos Rojas | **CREA:** 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance under number UY44WX94W79M, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br