
Kit Cãnula Low Vert

Reg. ANVISA: 80195529012.

BRAMSYS

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LOW VERT CANNULA KIT

PURPOSE

The Low Vert Cannula Kit is indicated for access in minimally invasive vertebroplasty procedures for the treatment of vertebral fractures.

OPERATING PRINCIPLE

1. Faceted working cannula

Straight faceted tip of the Cannula to facilitate the insertion of the bone cement Cannula.

Faceted Mandrel tip to facilitate puncture to the desired target.

11G diameter to ease the insertion and exchange of other items in the product.

Centimeter markings to assist the surgeon during the puncture.

2. Drill

Creates a cavity in the vertebra for the application of bone cement.

3. Bone cement application cannula

Straight tip of the Cannula and mandrel (not sharpened) to facilitate the application of the bone cement.

Hub with luer lock connection to facilitate the transfer of the cement into the cannula.

COMPOSITION

(1 or 2) Faceted working cannulas, (1) Drill, (3 or 6) Bone cement application cannulas.

COMMERCIAL PRESENTATION

Primary Packaging: Blister containing (1 or 2) Faceted working cannulas, (1) Drill, and (3 or 6) Bone cement application cannulas.

Secondary Packaging: Individual box.

		Code	
		LVU12011	LVB12011
Faceted cannula	Length (mm)	120mm	120mm
	Diameter (G)	11G	11G
	Qty	1	2
Drill	Length (mm)	187mm	187mm
	Diameter (G)	13G	13G
	Qty	1	1
Application cannula	Length (mm)	187mm	187mm
	Diameter (G)	11G	11G
	Qty	3	6

USE INSTRUCTIONS

- Ensure that the packaging is not damaged.
- Remove all items from the sterile packaging in a clean environment using aseptic techniques.
- Ensure all items are intact.
- Position the patient according to the surgical procedure to be performed.
- Define the puncture site for vertebral access.
- Insert the faceted working cannula to the desired target.
- Remove the mandrel from the faceted working cannula.
- Insert the Drill into the faceted working cannula.
- Rotate the Drill clockwise until no more bone friction is felt.
- Remove the Drill from the faceted working cannula.
- Remove the mandrel from the bone cement application cannulas and fill them.
- Insert the bone cement application cannula into the faceted working cannula.
- Insert the mandrel into the bone cement application cannula to transfer the cement to the vertebra.

NOTE: Repeat steps 11 to 13 until the vertebra is filled as desired.

- Remove the bone cement application cannula from the faceted

working cannula.

- Remove the faceted working cannula.
- Dispose of all used materials properly after the procedure.

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging must not be punctured or torn, and it should be kept away from sharp or pointed instruments.

Do not place under weight or heavy volumes, as this may damage the product's technical characteristics.

Do not place it near or together with different categories, as they may damage its technical characteristics.

It must be stored and transported following good practices for the transport of sterilized products.

This product is sterilized by ethylene oxide (ETO), and the packaging must not be damaged before use.

The product must be transported in its original packaging.

Storage temperature: Between 15°C and 30°C.

CONTRAINDICATIONS

The use of this product is not recommended in patients with severe hemorrhages, bleeding disorders, or those undergoing anticoagulant treatment. Medical opinion is essential when considering surgery in patients under these conditions.

ADVERSE EFFECTS

Adverse effects include, but are not limited to: discomfort, hematoma formation, infection, bleeding, damage to adjacent tissue, and blood vessel perforation.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- Avoid using the product if the packaging is damaged;
- Do not use the product if it has passed its expiration date;
- The product is intended for single use and must be discarded after use;
- Reuse is strictly prohibited;
- Reuse or resterilization may compromise the structural integrity of the product and/or cause damage that could result in injury or illness to the patient;
- Reuse or resterilization may also pose a risk of product contamination and/or cause infection to the patient or cross-infection, including but not limited to the transmission of infectious diseases between patients;
- It is advisable to avoid using the product if especially hazardous maneuvers are required, or when it is subject to torsional or bending forces that could compromise its integrity and proper functioning;
- Post-surgical care may vary depending on the technique used and the physiological conditions of each patient;
- A thorough clinical judgment is necessary when considering surgery in patients undergoing anticoagulant therapy or with bleeding disorders;
- This product must be used by a physician fully familiar with the indications, contraindications, limitations, typical outcomes, and potential adverse effects of vertebral puncture procedures.

MS Registration: 80195529012

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

The company BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency, guarantees for 24 MONTHS, from the sterilization date, the quality and sterile condition of this product, except in cases of misuse or non-compliance with the guidelines contained in this document.

SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC)
SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC)
CUSTOMER SERVICE (CS)
Phone: +55 11 5013 0010
E-mail: sac@bramsys.com.br
Site: www.bramsys.com.br