
Probe Reutilizável

Reg. ANVISA: 80195529011.

BRAMSYS

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REUSABLE PROBE

PURPOSE

The Reusable Probe is intended to promote the modification of peripheral nerves for pain relief.

OPERATING PRINCIPLE

The probe must be connected to the generator unit to conduct energy in a controlled manner to the tip of the cannula.

COMPOSITION

(1) Reusable probe.

The product may be accompanied by its accessory (optional): a sterilization case.

Commercial Presentation

Primary Packaging: Ziplock: (1) Reusable probe.

Secondary Packaging: Transport box.

Code	Length	Quantity of cannulas
REPJ-050	50mm	1
REPJ-100	100mm	1
REPJ-150	150mm	1
REPK-050	50mm	1
REPK-100	100mm	1
REPK-150	150mm	1
REPT-050	50mm	1
REPT-100	100mm	1
REPT-150	150mm	1

USE INSTRUCTIONS

Prerequisite: The following instructions assume that the probe is clean and sterilized.

- Connect the probe to the adapter cable;
- Connect the adapter cable to the energy generator unit;
- Insert the probe into the cannula already positioned in the patient;
- Perform the procedure as directed by the qualified professional;
- Remove the probe from the cannula;
- Follow the reprocessing instructions.

REPROCESSING INSTRUCTIONS

Pre-cleaning

• Remove visible residues from the instruments immediately after use. This can be done with running water and a disposable cloth that does not leave residues, suitable for hospital environments (Wipers).

Manual cleaning

- Use soft brushes and enzymatic detergents to clean the instruments. Make sure to follow the manufacturer's detergent instructions.
- Follow the cleaning procedure according to RDC 15 - Provides for good practice requirements for the processing of health products and other provisions, or its successor, and the ANVISA Cleaning and Surface Disinfection Manual.
- The product should not be submerged or cleaned by ultrasonic cleaning equipment.

Sterilization

• The medical device should be sterilized in a steam autoclave at 135°C for 15 to 25 minutes. It is necessary to wait for the product to reach room temperature before use or new sterilization cycle, according to the sterilization manual and RDC 15 - Provides for good practice requirements for the processing of health products and other provisions, or its successor.

OPERATING CONDITIONS

- The product must be clean and sterilized before use;
- The product should only be used by a duly qualified physician;
- The product must be discarded at the end of its useful life, in compliance with current disposal laws and regulations.

STORAGE AND TRANSPORT CONDITIONS

To maintain product integrity, it should:

- Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- Not be placed under weight or volumes, as they may damage the product's technical characteristics.
- Not be placed near or with different categories, as they may damage its technical characteristics.
- Be stored and transported following good practices for the transport of sterilized products.
- Be transported in its original packaging.

CONTRAINDICATIONS

- Infection or hematoma near the procedure site;
- Pregnancy;
- Preoperative neurological changes are a relative contraindication;
- Coagulation disorders, such as during anticoagulant therapy, are a relative contraindication;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast agents, antibiotics, and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical, or psychological condition that increases the risk of reduced procedure performance.

ADVERSE EFFECTS

Adverse effects include, but are not limited to, discomfort during and/or after the procedure.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The product is supplied in a non-sterile manner and must be sterilized before use;
- Before using the product, the user must read the instructions for use carefully, which will assist in the correct and safe use and protect against possible risks;
- The product's lifespan is significantly influenced by the degree of care in handling;
- This product should not be used for any purposes other than those stated in the instructions for use;
- The electrodes must not be bent at any time, either during the asepsis/sterilization process or during the operational procedure;
- The protective tube should be inserted into the Probe before sterilization;
- Before use, the condition of the cable, connector, and rod should be inspected. If any irregularity is found, it should not be used;
- The use of the Bramsys sterilization case accessory for autoclave sterilization is recommended;
- The product should only be used by qualified medical professionals who have the necessary knowledge and experience for its proper use, know how to handle reusable products, and especially understand the importance of maintaining a completely sterile environment in operating rooms, where care and attention to every detail are essential to prevent risks and contamination;
- Cleaning should be carried out according to WHO-established processes for Reusable Probes;
- Creutzfeld-Jakob Disease (CJD): If this device comes into contact with or is used on a patient with suspected or confirmed CJD, the device must not be used again and must be incinerated due to the risk of cross-contamination.
- The product can be re-sterilized for up to 20 cycles of use. Exceeding this amount may compromise the product's functionality, as well as the safety of the procedure to be performed.

MS Registration: 80195529011

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency, guarantees the quality of this product for **12 MONTHS** or **20 CYCLES OF USE** from the date of manufacture, except in cases of misuse or non-compliance with the instructions contained in this document.

SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC)
SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC)
CUSTOMER SERVICE (CS)
Phone: +55 11 5013 0010
E-mail: sac@bramsys.com.br
Site: www.bramsys.com.br