
Cabo adaptador p.Connect

Reg. ANVISA: 80195529006.

BRAMSYS

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P.CONNECT ADAPTER CABLE

PURPOSE

The p.Connect Adapter Cable is designed to connect electrodes, cannulas and cannulated introducers with a thermocouple mandrel to Radiofrequency Generators available on the market for sensitive stimulation, motor stimulation, ablation and neuromodulation with controlled temperature.

INDICATIONS OF USE

The p.Connect Adapter Cable is connected to electrodes, cannulas and cannulated introducers with a thermocouple mandrel, where it is connected to the Radiofrequency Generator for functions of sensitive stimulation, motor stimulation, ablation and neuromodulation with controlled temperature.

COMPOSITION

1 adapter cable.

USE INSTRUCTIONS

- Make sure the packaging is not damaged;
- Remove the p.Connect Adapter Cable from its packaging;
- Make sure that the p.Connect Adapter Cable is intact;
- Connect one of the p.Connect Adapter Cable connectors to the product used (electrode, cannula or cannulated introducer with thermocouple mandrel) and another connector to the Radiofrequency Generating Unit;
- After completing the procedure, disconnect the p.Connect Adapter Cable from the Radiofrequency Generator Unit and from the product used;
- Clean the product with the method recommended by the manufacturer before storing it.

CLEANING METHOD

- To remove visible residue, clean the product with isopropyl alcohol and a soft, lint-free cloth (gauze wipes can be used);
- Do not immerse the connectors in any type of liquid;

"In case of disposal, it must meet the requirements of RDC No. 222, of March 28, 2018 and other legislation in force at the federal, state and municipal level."

COMPATIBILITY WITH OTHER MEDICAL PRODUCTS

The p.Connect Adapter Cable was developed to work with radiofrequency generators compatible with the following technical characteristics:

- Electric voltage supply range: 90 - 260VAC
- Power supply frequency: 47 - 63Hz
- Maximum electrical output power: 50Wrms
- Wave frequency: 100 - 500kHz
- Output voltage range: 0 - 150V
- Output electric current range: 0 - 999mA
- Output pulse frequency: 1 - 200Hz
- Output wavelength: 0.1 - 50ms
- Impedance measurement: 0 to 5000ohms

STORAGE AND TRANSPORT CONDITIONS

The p.Connect Adapter Cable should be stored and transported in a cool, dry place, avoiding exposure to extreme temperatures.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- DO NOT STERILIZE BY ETHYLENE OXIDE (ETO), AUTOCLAVE AND/OR RADIATION;
- The correct and safe use of this product for its intended purpose is based on the conditions described in this instruction for use.
- The product should only be used by qualified medical professionals who have the necessary knowledge and experience in the proper use of the product.

Note:

- Do not use the adapter cable if the packaging or the product is

damaged in any way.

- Do not use this product if it shows any visible damage.
- Only use with the electrode, cannula or cannulated introducer combination with indicated thermocouple chuck and generator.
- Do not proceed with the procedure unless the electrode, cannula or cannulated introducer with thermocouple mandrel/generator are properly connected.

Batch Number and Expiry Date: see product packaging.

MS Registration: 80195529006

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees the quality and condition of this product for an **INDETERMINED** period, except in cases of misuse or non-compliance with the guidelines contained in this document.

**SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC)
SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC)
CUSTOMER SERVICE (CS)**

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br

MODELOS / MODELS

PCNT-AK50CM; PCNT-AT50CM; PCNT-AV50CM;
PCNT-BN50CM; PCNT-BS50CM; PCNT-BY50CM;
PCNT-CI50CM; PCNT-DS50CM; PCNT-MC50CM;
PCNT-MR50CM; PCNT-NT50CM; PCNT-RS50CM;
PCNT-SN50CM; PCNT-SR50CM; PCNT-TP50CM.