
Kit Cãnula HealBone

Reg. ANVISA: 80195529010.

BRAMSYS

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HEALBONE CANNULA KIT

PURPOSE / INDICATION FOR USE

The HealBone Cannula Kit is designed to collect, filter, prepare, and apply bone marrow stem cells to a patient, aiming to provide pain relief and improve functionality through regeneration.

PRINCIPLE OF OPERATION / MECHANISM OF ACTION

The HealBone Cannula Kit is designed to collect, filter, prepare, and apply bone marrow stem cells to a patient, aiming to provide pain relief and improve functionality through regeneration.

Functionality by Item:

1. Aspiration Cannula

Diamond tip to facilitate puncture.

Fixed tip to prevent contamination from peripheral blood.

Alternating side holes to provide greater variety during collection.

11G gauge for efficient aspiration of bone marrow stem cells through a minimally invasive approach.

Centimeter markings to facilitate cannula retraction and collection at various levels.

Luer lock connection to ensure secure sealing during collection.

2. Depth Control Device

Assists in retracting the collection cannula for aspiration at various levels.

3. Filter

Removes bone fragments and other debris.

Filtration membrane with specifically sized pores to allow the passage of cells and growth factors.

Luer lock connections at both ends to facilitate connection.

Transparent casing for visualization of the filtration process.

4. 3-Way Adapter

Enables secure connection of devices for handling and collection.

Three ports with luer lock connections to ensure secure sealing.

5. Extension Tube

Enhances ergonomics during the use of the application cannula.

Flexible, transparent tube for visualization during application.

Luer lock connections at both ends to ensure secure sealing.

6. Injection Cannula

Application of bone marrow stem cells to the injury site or treatment area.

Sharp tip to facilitate tissue penetration.

18G or 20G gauge for application.

Length of 100mm or 150mm.

Luer lock connection to ensure secure sealing.

Echogenic tip for ultrasound visualization.

Centimeter markings.

7. Planning Scale and Surgical Marker

Assists in marking and planning the desired target during the procedure.

COMPOSITION

(1) Aspiration cannula, (1) Filter, (1) 3-way adapter, (1) Depth control device, (1) Extension tube, (1) Injection cannula, (1) Surgical marker, (1) Surgical planning scale.

COMMERCIAL PRESENTATION

Primary Packaging:

Sealed blister with Tyvek paper containing:

(1) Aspiration cannula, (1) Filter, (1) 3-way adapter, (1) Depth control device, (1) Extension tube, (1) Injection cannula, (1) Surgical marker, (1) Surgical planning scale.

Secondary Packaging:

Individual packaging box containing the primary packaging.

Code	Injection Cannula			Aspiration Cannula
	Length (mm)	Diameter (G)	Qty.	Qty.
KCHB-18100	100	18	1	1
KCHB-18150	150	18	1	1
KCHB-20100	100	20	1	1
KCHB-20150	150	20	1	1

USE INSTRUCTIONS

- Ensure that the packaging is not damaged.
- Remove all items from the sterile packaging in a clean environment using aseptic techniques.

• Ensure that all items are intact.

• Position the patient according to the surgical procedure to be performed.

Note: If necessary, use the Surgical Planning Scale and Surgical Marker to assist with puncture planning and locating the desired target.

• Prepare the skin puncture area using aseptic and antiseptic techniques.

• Perform the percutaneous puncture using the Aspiration cannula.

Note: If necessary, use fluoroscopy, ultrasound, or X-ray to assist with puncture and target monitoring.

• Reach the deepest point to begin the collection.

• Using the Depth Control Device, retract the cannula to perform additional collections as needed.

• Remove the cannula after completing the collection.

Note: If necessary, use the filter to remove bone fragments and other debris, and the 3-way adapter to prepare the collected substance before administration.

• Position the patient according to the desired application area.

• Prepare the skin puncture area using aseptic and antiseptic techniques.

• Perform the percutaneous puncture using the injection cannula.

Note: If necessary, use the extension tube to apply the collected material.

• Remove the cannula after application.

• After the procedure is complete, all materials should be decharacterized and disposed of in an appropriate location, following applicable regulations for the disposal of sharp and contaminated materials adopted by the hospital.

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it should:

Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging should not be punctured or torn and should be kept away from cutting or sharp instruments.

Do not place under weight or pressure, as it may damage the product's technical characteristics.

Avoid placing it next to or with different categories of products, as this could damage its technical characteristics.

Be handled and transported following best practices for transporting sterilized products.

This product is sterilized with ethylene oxide (ETO) and must not have damaged packaging before use.

The product should be transported in its original packaging.

Storage temperature: Between 15°C and 30°C.

CONTRAINDICATIONS

It is not recommended for use in patients receiving intensive anticoagulant treatment or those with severe hemorrhagic conditions. Medical advice is essential when considering performing a biopsy in patients with bleeding disorders or those on anticoagulant therapy.

ADVERSE EFFECTS

Adverse effects may include, but are not limited to: discomfort, bruising, infection, bleeding, damage to adjacent tissue, and blood vessel perforation.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Avoid using the product if the packaging shows damage.

Do not use the product if it has passed its expiration date.

This product is intended for single use only and must be discarded after use. Reuse is strictly prohibited.

Reuse or re-sterilization may compromise the product's structural integrity and/or cause damage that can result in injury or illness to the patient.

Reuse or re-sterilization may also present a risk of contamination and/or cause infection in the patient or cross-contamination, including but not limited to the transmission of infectious diseases between patients.

It is advisable to avoid using the product if particularly dangerous maneuvers are required, or when subjected to twisting or bending forces that could compromise its integrity and proper functioning.

Post-collection care may vary depending on the technique used and the physiological conditions of each patient.

A thorough clinical judgment is required when considering performing a collection in patients undergoing anticoagulant therapy or with bleeding disorders.

This product should be used by a physician who is fully familiar with the indications, contraindications, limitations, typical results, and possible adverse effects of using the collection puncture.

MS Registration: 80195529010

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

The company BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency, guarantees for 24 MONTHS, from the sterilization date, the quality and sterile condition of this product, except in cases of misuse or non-compliance with the guidelines contained in this document.

SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC)
SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC)
CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br