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WARNING: Please read the instruction manual carefully before using the equipment.

Responsável Técnico / Responsable Técnico / Technical Responsible:	Marcos Rojas CREA 5061214879
---	---------------------------------

Fabricado por / Manufactured by: Bramsys Indústria e Comércio Ltda.
Av. Dr. Alberto San Juan, 350
Jd. Baccili – Tietê – SP
CEP 18530-000 – Brasil / Brazil
☎ +55 11 5013-0010
✉ sac@bramsys.com.br
www.bramsys.com.br

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Bramsys reserves the right to make any changes without prior notice.

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TABLE OF CONTENTS 

1. INTRODUCTION	37
2. INSTRUCTION FOR USE	38
2.1 Initial precautions	38
2.2 Installation	39
2.3 Attaching the frame to the patient	39
2.4 Assembling the segments	41
2.5 Preparation for tomographic study	41
2.6 Assembly and verification of the Stereotactic Apparatus	45
2.7 Use of the cranial fixture adapter	43
2.8 Scales	46
2.9 Tolerance	46
2.10 Sterilization	47
3. PRECAUTIONS, RESTRICTIONS, WARNINGS AND SPECIAL CARE	47
3.1 Precautions	47
3.2 Restrictions	47
3.3 Warnings and cautions	47
4. STORAGE, CLEANING, CONSERVATION AND TRANSPORT	48
4.1 Storage	48
4.2 Cleaning and disinfection	48
4.3 Lubrication	48
4.4 Conservation	48
4.5 Transport	48
5. DESCRIPTION OF PRODUCT EFFECTIVENESS AND SAFETY	48
5.1 Undesirable side effects and contraindications	48
5.2 Preventive maintenance	49
5.3 Corrective maintenance	49
5.4 Description of safety and efficacy	49
6. SYMBOLS	49
7. WARRANTY TERMS	49

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1. INTRODUCTION

The BMS-125E stereotaxic apparatus system was designed to assist the neurosurgeon in locating targets for performing procedures such as biopsies, brain abscess drainage, guided surgeries, hematoma drainage and removal, isotope implantation for inside-out or intracavitary irradiation, applications in brachytherapy, electrodes implantation for stimulation (DBS), and surgeries involving encephalic/cervical structural injuries.

ATTENTION: this equipment should only be used by a specialized medical professional.

– The Bramsys BMS-125E stereotaxic apparatus consists of the following items:

- 1 open frame
- 4 extension bars
- 4 extension bar fixing screws
- 4 thread reducer adapters
- 4 allen keys (2 of 4mm/2 of 5mm)
- 2 millimeter rulers
- 6 millimeter ruler fixing screws
- 1 graduated arch
- 1 trolley
- 2 millimeter tubes
- 1 long “Z” vertical bar
- 1 long “Z” vertical bar guide
- 1 short “Z” vertical bar
- 1 short “Z” vertical bar guide
- 3 tomography benchmarks (2 lateral and 1 frontal, both large)
- 3 tomography benchmarks (2 lateral and 1 frontal, both small)
- 6 aluminum screws for fixing the benchmarks
- 1 straight segment
- 1 curved segment
- 2 aluminum screws for fixing the segments

– Accessories:

- 2 to 3 storage boxes
- 1 30mm stainless steel screw
- 4 50mm stainless steel screws
- 4 75mm stainless steel screws
- 4 75mm pediatric stainless steel screws
- 4 90mm pediatric stainless steel screws
- 1 set of isocentric retractor (2 isocentric retractors and 6 spatulas)
- 2 stainless steel universal adapter
- 2 titanium retractor arms
- 2 stainless steel retractor arms
- 1 pediatric fixture (2 lateral and 1 frontal)
- 2 ear aligners
- 1 laser guide
- 3 reusable bushing adapters for biopsy
- 3 reusable bushing adapters for guided surgery
- 3 reusable bushing adapters for drainage
- 3 reusable bushing adapters for electrode implantation
- 3 MRI benchmarks (1 frontal and 2 lateral)

- 1 open carbon fiber frame
- 1 carbon fiber straight segment
- 1 carbon fiber curved segment
- 2 segment fixing screws
- 4 carbon fiber extension bars
- 4 carbon fiber thread reducer adapters
- 1 neuro endoscope holder
- 4 reusable 30mm M8 aluminum screws - adult
- 4 reusable 50mm M8 aluminum screws - adult
- 4 reusable 75mm M8 aluminum screws - adult
- 4 reusable 50mm M8 aluminum screws - pediatric
- 4 reusable 75mm M8 aluminum screws - pediatric
- 4 reusable 90mm M8 aluminum screws - pediatric
- 4 50mm stainless steel screw - pediatric
- 4 30mm titanium M8 screws - adult
- 4 50mm titanium M8 screws - adult
- 4 75mm titanium M8 screws - adult
- 4 50mm titanium M8 screws - pediatric
- 4 75mm titanium M8 screws - pediatric
- 4 90mm titanium M8 screws - pediatric
- 4 30mm plastic M8 screws with titanium tip - adult
- 4 50mm plastic M8 screws with titanium tip - adult
- 4 75mm plastic M8 screws with titanium tip - adult
- 4 50mm plastic M8 screws with titanium tip - pediatric
- 4 75mm plastic M8 screws with titanium tip - pediatric
- 4 90mm plastic M8 screws with titanium tip - pediatric
- 4 30mm plastic M8 screws with ceramic tip
- 4 50mm plastic M8 screws with ceramic tip
- 4 75mm plastic M8 screws with ceramic tip
- 4 50mm plastic M12 screws with ceramic tip
- 4 75mm plastic M12 screws with ceramic tip
- 4 50mm plastic M12 screws with titanium tip
- 4 75mm plastic M12 screws with titanium tip
- 1 Bilateral arch set (2 half-graduated arches, 1 trolley, 1 millimeter tube, 1 long Z vertical bar, 1 long Z bar guide, 1 millimeter ruler)
- 1 adapter trolley for microregistration (according to the microregister model)
- 1 benchmark holder

– Optional parts:

- 1 cranial fixture adapter (as per the cranial fixture model)
- 1 tomography table holder (as per the model of the tomograph)
- MNPS MEVIS surgical planning software
- Transport case

2. INSTRUCTION FOR USE

2.1 Initial Precautions

- When unpacking the equipment, check that there is no apparent damage caused by dropping or improper handling during transport. If any damage is observed, please contact the manufacturer/representative immediately for the measures to be taken.

- Before starting the procedure, lubricate all parts that slide together, the ends of the locking screws of the millimeter side ruler holder and the vertical “Z” bar holder, the underside of the washers, and the contact surface of the arch on the millimeter tube using silicone or sterile Vaseline. This will allow less effort to loosen or tighten the screws and, consequently, the fixation will be improved, in addition to prolonging the good performance of the equipment.

2.2 Installation

- Examine all accessories before starting surgery. Make sure they are working correctly. Incorrect assembly can result in malfunctions and undesirable surgical effects.

ATTENTION: Bramsys will be responsible for the safety, reliability, and performance of this equipment as long as it is used in accordance with the instructions for use provided in this manual.

2.3 Attaching the frame to the patient


- **Target position**

1. Check the target position through a previous survey before fixing the frame to the patient.
2. Position the frame above or below the region of interest, never at the same height, as this will prevent visualization.
3. Thread the screws directly into the frame or extension bars using the thread adapters.

WARNING: select the screw length according to the size of the patient’s head, making sure that the screw does not exceed the frame parallels.


- **Fixing with 3 screws**

For the 3-screw fixation, use 1 small or medium screw at the front and 2 at the rear. When using the millimeter ruler on the front of the stereotactic frame, use the small front screw only so that it does not exceed the straight surface of the frame.

- **Fixing with 4 screws**

For the 4-screw fixation, use 2 front screws and 2 rear screws. When using the millimeter ruler on the front of the stereotactic frame, make sure that the front screws do not go beyond the straight surface of the frame


Note: screws are registered and sold separately.

- **Fixing with extension bars**


There are 4 interchangeable extension bars to be used when the fixing of the stereotactic frame directly to the patient's skull hampers the target area visualization.

Fix the extension bars to the stereotactic frame with Allen Key #4 using the corresponding screws.


- **Thread reducer adapters**

When using the extension bar, place the thread reducer adapters by threading them from the inside of the bars.

2.4 Assembling the segments

1. The need to assemble the curved or straight segment to increase the working area must be evaluated by the person in charge of surgery and assembly.
2. In order to replace the straight segment with the curved one, use the Allen Key #5 to loosen the fixing screws, replace the parts and tighten the screws again, making sure the assembly is properly rigid.


2.5 Preparation for tomographic study

- **Mounting the holder on the tomograph table**


1. Fix the holder on the tomograph table.
2. After unlocking the latch, fit the frame to the holder, moving it until you find the most comfortable position for the patient.
3. Lock the frame holder with the plastic lever and attach the frame to the holder with the two aluminum screws.

The screws that lock the frame into place must always be in the opposite position to the benchmarks, as shown in the figures below, thus preventing the artifact from the screw head from being produced in the image.

Upward Mounting of the Benchmarks


Downward Mounting of the Benchmarks

- **Placing the tomography benchmarks on the frame**

Place the benchmarks (2 lateral and 1 frontal) on the stereotactic frame, fitting them into the corresponding holes, with the frontal benchmark being distinguished by having a hole in the center and the lateral ones having the word "Occipital" in the back of the patient's head. According to the positioning of the frame on the patient's head, the benchmarks can be placed "UPWARD" or "DOWNWARD".

WARNING: when using the extension bars, the position of the benchmarks should always be "UPWARD."

WARNING: incorrect assembly of the benchmarks leads to errors in the calculation of stereotactic coordinates.

- **Calculation of coordinates with the MNPS software**

Converting DICOM images to IMG

1. Insert the CD with the tomograph stereotactic images into the PC.
2. Copy them to the MNPS/DICOM directory.
3. Open the TAC directory;
4. Drag any image over the CT MevDICOM icon.
5. Show 12 images with the fourth button.
6. Sort the images from the menu Options → Sort images → z (transversal).
7. Cut a frame in which the benchmarks are included (use the fifth button).
8. Double-click on the red frame.
9. Deactivate the Scout View images or those with linked benchmarks by using the Space key.
10. Go to File → Export to IMG or press the second button (Save).
11. Answer Yes to the prompt to open the MNPS.

Preparation of procedure planning

1. Select the Upward or Downward benchmark positions to be used by the software.
2. Select the marking position (right, left, frontal, occipital).
3. Select the type of frame to be used.
4. Perform a section-by-section inspection to check if the automatically found benchmarks are correct, approving the image with the Enter key.
5. In the mosaic window (where 9 images appear together), double-click on any one of them to go to the main MNPS window.
6. In the main window, go to a section that is approximately at the height of the tumor (biopsies) or at the height of the third ventricle junction (functional).
7. Adjust the window using the right mouse button on the main image.
8. For a finer tuning of the image contrast, click on the brightness and contrast button and, with the keyboard arrows, move the histogram levels on the right of the screen.

Definition of stereotactic coordinates in case of biopsy

1. Go to the level of a section approximately in the center of the tumor.
2. Zoom in on the lesion.
3. Start drawing the tumor perimeter. For such, click on the ROI Menu button (the one with a pencil and green paper) and then on the Draw ROI Segment option on that menu. After choosing the ROI number and writing its name, start drawing over the image. Remember that you can use the Delete button to delete part of an ROI.
4. Change sections and draw over other sections. Finish with the pencil button.
5. Save the planning by clicking on the button with the floppy disc icon.
6. Go to the ROI Menu again (pencil and green paper button) and select ROI center to POI to automatically set a point at the center of the tumor.
7. Using the POI Table button (or the Control+End keys), edit the coordinates of that point. To visualize the trajectory, write the word OUT in the FROM field of the central point of the tumor and adjust the desired angles for the biopsy cannula.
8. Save the planning and write down the coordinates to be used in the Stereotactic Apparatus.
9. To visualize the chosen trajectory, use the axial sections, visualize in 3D, and also in coronal and sagittal planes. You can also do an oblique plane visualization.
10. To switch from one view to another (coronal or sagittal), return to the axial view using the Go to button or the mosaic button.

Definition of stereotactic coordinates in case of functional surgery

1. Go to the AC level and check AC with the POI button. Do the same for PC and for IHP.
2. Rebuild the images parallel to the AC-PC plane with Alt-F6.
3. Modify the position of the AC and PC points at will or rebuild them with Alt-F6.
4. Measure the distance between the AC and PC points with the Control+D keys. It should be around 26 mm in normal anatomy.
5. Check the lateral position of the IHP by observing the anatomy symmetry at the ocular orbit height in AC-PC parallel planes.
6. Rebuild planes inclined at 7 degrees in relation to the AC-PC plane using Shift-F6.
7. Activate the atlas using the atlas button.
8. Check that the atlas matches the anatomy in the regions close to the third ventricle and the cerebral aqueduct.
9. If it does not, go back to the AC-PC plane and modify the AC or PC point.
10. Locate the target structure and place a POI within it.
11. Edit the POI table using the POI Table button, registering the coordinates of the target point (GPI etc.) that will be used in the Stereotactic Apparatus.
12. To capture the planning images, use Alt-C.

Image fusion

If the structures to be worked on are poorly defined in the tomography, the image fusion with the magnetic resonance can be used.

1. Insert the images in the DICOM filter for MRI and write down the patient's name, indicating the orientation of the images (coronal, sagittal, or axial).
2. Sort the images.
3. Save the images.
4. Exit the DICOM filter.
5. In MNPSShell, open Planning → CAT3D for Fusion.
6. In CAT3DShell, open CAT3D and load the patient information/images (MRI) with Load Anatomy – Start New Patient.
7. Double-click on a mosaic image.
8. Mark 4 points, for example, OD, OI, Fast, and T.
9. Save the planning.
10. Go to CAT3D's start menu using the ESC key and select Export Fusion.

11. Return to Windows with the Windows key and select the window.
12. In MNPS, with the TAC, mark the same 4 points with the same 4 names.
13. Press the Fusion button and choose the patient (axial, coronal, or sagittal).
14. Move the resonance window with the Control+F5 keys.
15. Start the mutual information maximization with the Alt+F5 keys.
16. Check target point position and base points.
17. Exit fusion using the same button.
18. Visualize the atlas in CT with Shift+F6 and reset AC-PC.
19. Locate the target to be injured and write down the coordinates to be used in the Stereotactic Apparatus.
20. Save the planning and close the system.
21. We recommend backing up patient images and information.

2.6 Assembly and verification of the stereotactic apparatus

WARNING: all parts should slide easily against each other.

If there is any resistance, check if the part is correctly mounted, fitted, or if any fixing screws are causing interference.


If the problem persists, please contact Bramsys' technical assistance.

- **Trolley and adapters assembly**

Mount the adapters on the trolley and on the adjustable guide as follows: place the two longer adapters on the graduated part of the trolley and the shorter one on the underside of the adjustable guide.

The longer adapter has an indication of the position of the hole for fixing the cannula through the trolley locking screw.


Note: adapters are registered and sold separately.

- **Instrument assembly**

Insert the cannula or electrode into the trolley until the end of its course, where its base touches the adapter.


Note: cannulas and electrodes are registered and sold separately.

- **Crosshairs for isocenter control**

Check if the instrument is properly calibrated at the isocenter, as follows:

1. Mount the millimeter tube, arch, and trolley assembly.
2. Adjust the coordinates from the surgical planning in the Stereotactic Apparatus.
3. Insert the cannula or the electrode with the trolley in its "0" (zero) position.
4. Align the two ends of the crosshair looking through the millimeter tube.
5. Adjust the area of interest on the cannula or electrode so that it is positioned in the center of the millimeter tube crosshairs.
6. Tighten all motion lock screws to ensure the adjusted coordinates remain correct.
7. Complete the arch and trolley assembly on the stereotaxic frame and carry out the procedure.

2.7 Use of the cranial fixture adapter

1. Mount the adapter on the stereotactic frame and secure it properly.
2. Adjust the cranial fixture, choosing the best position for the procedure.

2.8 Scales

All scales are in millimeters (mm) and angles in sexagesimal degrees.

2.9 Tolerance

The equipment coordinates accuracy is $\pm 0.5\text{mm}$.

2.10 Sterilization

The stereotactic apparatus can be sterilized with ethylene oxide, hydrogen peroxide, or steam autoclave (requires a special storage box), according to the information in the table below:

Method	Sterilization Parameters	Notes
Ethylene oxide (ETO)	2 to 3 hours, between 50 and 60°C.	After sterilization, the equipment must undergo mechanical aeration for 240 minutes or environmental aeration for 24 to 72 hours.
Hydrogen peroxide (STERRAD)	Approximately 75 minutes.	-
Steam autoclave (requires special storage box)	Approximately 1 hour and 15 minutes at 134°C.	Do not sterilize on rapid cycles or above 150°C.

Do not use the steam autoclave method on the following accessories:

- Plastic screws with titanium tip, plastic screws with ceramic tip, MRI benchmarks.

3. PRECAUTIONS, RESTRICTIONS, WARNINGS, AND SPECIAL CARE

3.1 Precautions

Do not use the equipment and its accessories after their validity period. They may malfunction, making the user’s work difficult and/or causing harm to the patient. Although not using disposable items or generating waste after use, the equipment must be segregated for possible revalidation or discarded in accordance with current local legislation at the end of its useful life.

3.2 Restrictions

- All parts and accessories accompanying the equipment are for exclusive use with such equipment.
- The following products are recommended for use with the BMS-125E:
 - Cannulas for Bramsys stereotactic apparatus
 - Bramsys adapters
 - Bramsys screws
 - Bramsys electrodes for functional surgery

ATTENTION: the use of any part, accessory, or material not specified herein is the user’s sole responsibility.

3.3 Warnings and cautions

- Avoid dropping, twisting, crushing, impacting, or submitting the equipment to any actions that may deform it.
- Deformations in the structure of the equipment will compromise its accuracy and may cause damage to the patient.
- Properly secure the equipment’s moving parts.

- Movements of the stereotaxic coordinates, resulting from the accidental displacement of these parts, will change the position of the predetermined target, which may be harmful to the patient.
- The equipment must be properly cleaned and sterilized according to the instructions contained in the user manual. Failure to do so may result in bacterial and/or viral contamination, which may cause infectious complications for the patient.
- Excessive force when fixing the stereotactic frame to the skull beyond that supported by the bone structure may cause injuries to the patient.
- Improper reading of the equipment scales due to user misinterpretation or wear of the paint on the scales can cause bad execution of the procedure, which may be harmful to the patient.

4. STORAGE, CLEANING, CONSERVATION, AND TRANSPORT

4.1 Storage

- Keep the equipment in its original packaging, in a place protected from excessive humidity, rain, and direct sunlight.

4.2 Cleaning and disinfection

- Use a neutral soap and warm water only to clean the equipment externally after each procedure. Use compressed air or a clean, dry cloth to dry it.
- Do not use organic solvents, such as thinners, or detergents containing solvents, such as ether, stain remover, gasoline, etc., to clean the equipment. Do not use aerosol or liquid spray applicators.
- Do not use micro-abrasive material or steel wool for cleaning.
- Wash the box using a soft sponge or wet cloth with water and neutral soap, and dry it after washing.

4.3 Lubrication

- Lubricate all parts that slide together, the ends of the locking screws of the millimeter side ruler holder and the vertical "Z" bar holder, the underside of the washers, and the contact surface of the arch on the millimeter tube using silicone or sterile Vaseline before starting the procedure.

4.4 Conservation

- Keep the equipment clean and dry for the next use.
- Keep it in a clean and dust-free place.

4.5 Transport

- Avoid vibrations and impact, and do not drop the equipment.
- Do not allow the equipment to be exposed to rain or excessive humidity.
- All parts must be packed inside the shipping box.
- Check that all parts are correctly positioned and secured in place.

5. DESCRIPTION OF PRODUCT EFFECTIVENESS AND SAFETY

5.1 Undesirable side effects and contraindications

- This equipment does not produce any undesirable side effects if all the recommendations described in this document are strictly complied with. The equipment has no contraindications. However, it can only be used or operated by a specialized medical professional or supervised by such professional.

5.2 Preventive maintenance

- Return the equipment to Bramsys every 12 months for preventative maintenance to ensure its good working order and safety.

5.3 Corrective maintenance

- If there is a problem with the equipment or any of its parts, return it to Bramsys with a description of the problem.

WARNING: request technical assistance through the Customer Service Center via +55 (11) 5013-0010 or email at sac@bramsys.com.br.

5.4 Description of safety and efficacy

- Safety and efficacy in the use of the equipment can be ensured as long as the recommendations prescribed in the User’s Manual are strictly complied with.

6. SYMBOLS

	Fragile, handle with care		Careful! Consult accompanying documents
	This side up		Maximum stacking limit
	Keep away from sunlight		Keep dry
	Temperature limits		Nonsterile

7. WARRANTY TERMS

Bramsys Indústria e Comércio Ltda. guarantees its original user customers that the equipment it manufactures is produced using state-of-the-art technology and strict quality control, thus ensuring its perfect functioning, provided that the conditions and terms below are complied with.

WARRANTY

The warranty, without any charge to the buyer, is limited to the replacement and/or repair of any defective parts or the correction of any production defect upon verification by our technical assistance.

The replacement and/or repair referred to in the previous item does not apply to parts with natural wear and tear, as well as damage due to malpractice or negligence in the use of the equipment, or to parts that have been repaired or modified by people not accredited by Bramsys.

Under no circumstances, should there be a need to replace any component covered by this term, the warranty period will be extended by adding any supplements to the replaced component.

EQUIPMENT

This warranty is valid for all Bramsys equipment produced by Bramsys Indústria e Comércio Ltda.

INSTALLATION AND USE

The installation and/or operation of the equipment, as well as its working conditions, must comply with the Bramsys standards contained in this Instruction manual. The warranty clauses will be void if the conditions are different from those indicated.

WARRANTY LOCATION

Equipment repair will be carried out at our facilities. The transport freight (round trip) will be paid by the buyer.

WARRANTY PERIOD

Starting from the date of issue of the sales invoice by Bramsys or its resellers, being valid for a period of 12 months. The warranty period for accessories is determined by their useful life or ninety (90) days, whatever comes first.

RESPONSIBILITY

This warranty is only valid for equipment in use and in the possession of the original buyer. Bramsys' liability is limited to the replacement and/or repair of components, with the exception of any damages for lost profits or compensation for any other indirect or immediate damages.

Request technical Assistance through the Consumer Assistance Service via +55 (11) 5013-0010, or email sac@bramsys.com.br.

Revisão: 5.0 – 06/09/23