
Cânulas para Estereotaxia

Reg. ANVISA: 80195520014.

BRAMSYS

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BRAMSYS STEREOTAXY CANNULAS

PURPOSE

The Bramsys Biopsy Cannulas are intended for stereotactic surgical procedures.

INDICATIONS OF USE

- The Biopsy Cannulas are specifically indicated for removing tissue samples from the brain region.
- The instructions for use that accompany the product must be read before opening the sterile packaging to ensure its correct and safe use.

COMPOSITION

1 Cannula, 1 Mandrel, 1 Stopper Cap, and 1 Extension Tube. All parts and components of the product are marketed for exclusive use with the same. The cannula may be accompanied by adapters and screws (registered separately under numbers 80195520007 and 80195520006, respectively).

USE INSTRUCTIONS

- Ensure that the stereotactic frame is properly mounted on the patient;
- Ensure that the surgical planning of the stereotactic coordinates is ready and reviewed;
- Verify that the packaging is not damaged or wet;
- Verify that the sterilization indication is correct;
- Check the sterilization expiration date;
- Open the packaging box and remove the sealed blister;
- Following sterile handling procedures, open the blister and remove the product;
- Locate the assembly location of the adapters on the equipment;
- Assemble the adapters in the designated locations, observing the correct order;
- Perform the pre-assembly and adjustment of the coordinates on the stereotactic equipment, ensuring that the cannula window is at the target;
- The biopsy cannula length to the center of the window is set at 243 mm. If using equipment that requires different lengths, position the stopper cap to meet the stereotactic coordinates to reach the desired target;
- Attach the male luer connector of the extension tube to the mandrel and turn it clockwise to secure it;
- Attach a syringe to the female luer connector of the extension tube, turning it clockwise until secured;
- Perform the final assembly of the stereotactic frame on the patient, insert and secure the mandrel into the cannula with the window closed, and gently introduce it through the adapter hole until it reaches the target point in the brain;
- Secure the locking screw of the carriage to fully lock the cannula movement;
- Turn the mandrel to open the window, ensuring that the marks are aligned;
- Pull the syringe plunger to create negative pressure and keep it pulled;
- Turn the mandrel 180° to close the window and collect the material that entered the cannula;
- Release the syringe plunger;
- Carefully remove the mandrel;
- Remove the collected material that will be in the window;
- Reinsert the mandrel and repeat the procedure in as many quadrants as necessary;
- After the procedure, insert the mandrel, close the window, unlock the cannula, remove it from the equipment, dispose of it, and discard it in MEDICAL WASTE.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will assist in its correct and safe use and protect against potential risks.
- Insertion and removal of the cannula in the patient's brain must be done carefully and always with the "window" closed, otherwise, it may cause serious brain damage to the patient, including the risk of death.
- Each type of collected material has different consistency, which may vary the quantity and difficulty of collection.
- Deformations in the stereotactic system structure can compromise the procedure's accuracy, posing a risk of death to the patient.
- Properly secure the movable parts of the stereotactic equipment.
- Excessive force in securing the stereotactic frame to the skull, beyond what the bone structure can support, may compromise the cranial cavity and cause damage to the patient.
- Incorrect placement of the adapters on the carriage can cause a discrepancy in the planned stereotactic coordinate, compromising the system's accuracy and posing a risk of death to the patient.
- Do not use the product past the expiration date; always verify the sterilization indication is correct.
- The product should only be used by qualified medical professionals who possess the necessary knowledge and experience for its proper use, knowledge of handling sterile products, and, most importantly, the importance of maintaining a completely sterile environment in operating rooms, where care and caution with every detail are fundamental to prevent contamination risks.
- The product should not be reused and/or re-sterilized; it is for SINGLE use only and should be destroyed and discarded in MEDICAL WASTE.
- The product must always be transported and stored protected from all types of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging should not be punctured or torn and must be kept away from sharp or pointed instruments.
- Do not store or keep this product near products of different categories.
- Verify the integrity of the packaging before using the product; if compromised, DO NOT USE, as its sterility guarantee will be compromised.
- Reprocessing this product is prohibited.
- This product is single-use, sterilized by Ethylene Oxide (ETO) or Gamma Radiation.
- This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

A completely safe product when used by qualified professionals, respecting the usage specifications described in this document.

Lot number, Sterilization and Expiration Date: refer to the product packing

MS Registration: 80195520014 / 80195520007

Technical Responsible: Marcos Rojas | CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

Commercial Models	Cannula				Mandrel		
	Length to the center of the window	Useful length	Window length	Qty	Length	Window length	Qty
BIOPSY CANNULA 14G-243/5	243mm	235mm	5mm	1	262mm	5mm	1
BIOPSY CANNULA 14G-243/10	243mm	235mm	10mm	1	262mm	10mm	1
BIOPSY CANNULA 14G-243/8 – BMS-D490	243mm	235mm	8mm	1	262mm	8mm	1