
Eletrodo Bipolar PLMe

Reg. ANVISA: 80195520025.

BRAMSYS

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ELECTRODE BIPOLAR PLMe

PURPOSE

The Bipolar PLMe is intended for producing soft tissue ablation and coagulation, as well as for hemostasis of blood vessels in arthroscopic procedures.

COMPOSITION

(1) Electrode

COMMERCIAL PRESENTATION

Primary Packaging: Sealed blister with Tyvek paper containing (1) electrode.

Secondary Packaging: Individual packing box containing the primary packaging

Code	Tip model	Drive
PLME-BTW	Bolt	Pedal
PLME-IFW	Infinity	Pedal
PLME-SBW	Super Bolt	Pedal
PLME-STW	Star	Pedal
PLME-XGW	Xgon	Pedal

USE INSTRUCTIONS

- Check that the packaging is not damaged or wet.
- Ensure the sterilization indication is correct.
- Verify the sterilization expiration date.
- Open the packaging and remove the sealed product.
- Following sterile product handling procedures, open the second package and remove the product.
- Prepare the patient for the procedure according to standard procedures.
- For suction functionality, attach the suction line adapter to normal hospital suction equipment.
- Connect the end of the Bipolar PLMe cable to the receptacle of the generator unit.
- Insert the Bipolar PLMe through the trocar and place it in contact with the tissue for ablation.
- Select the ablation function using the manual or foot pedal control, depending on the model. The electrodes should always be in contact with the target tissue, and the electrode must always be in motion when active. The Bipolar PLMe will emit a single-tone beep when activated. Some bubbling and an orange glow of the Bipolar PLMe are normal. If a different setting is desired, adjust the configuration using the foot pedal or manual control, or the control buttons on the generator unit's front panel.
- Periodically check the electrodes to ensure they are intact and that the PLMe is functioning as expected. Use the recommended predefined point to achieve the desired result.
- For coagulation, select the coagulation function using the manual or foot pedal control, depending on the model. If a different setting is desired, adjust it using the front panel buttons of the generator unit until the desired configuration is achieved.
- After the procedure, turn off the generator unit, disconnect the suction tubing, disconnect the electrode, and dispose of it in HOSPITAL WASTE.

CONTRAINDICATIONS

- Non-arthroscopic surgical procedures.
- Procedures without the use of a conductive irrigant.
- Patients for whom an arthroscopic procedure is contraindicated for any reason.
- Patients with cardiac pacemakers or other implanted electronic devices.

ADVERSE EFFECTS

- Failure to apply the suction adapter and provide adequate suction may cause thermal injuries to the physician or patient.
- During the procedure, ensure adequate fluid flow and saline solution circulation to avoid unnecessary heating, which can lead to tissue injury.
- Electrosurgery may cause damage to surrounding tissues due to iatrogenic injuries.

STORAGE AND TRANSPORT CONDITIONS

The product must be transported and stored protected from any type of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging should not be damaged.

Keep the product away from sharp or pointed objects.

Check the packaging integrity before using the product. If the package is compromised, DO NOT USE, as sterility may be compromised.

This product is for single-use only and is sterilized with Ethylene Oxide (ETO)

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The Bipolar PLMe is designed for use with generator units compatible with the following specifications:
 - Operating Frequency: 100-450 kHz
 - Operating Voltage: 0-320 Vrms
 - Operating Power: 400W max
- Before use, users must carefully read the instructions to ensure correct and safe usage and to protect against potential risks.
- Do not use the product if the packaging is damaged.
- The product should only be used by qualified medical professionals with knowledge and experience in handling sterile products, maintaining a fully sterile environment in operating rooms, where careful attention to each detail is essential to avoid contamination risks.
- Do not touch the active tip of the Bipolar PLMe when it is activated.
- Only use conductive media (e.g., normal saline solution or Ringer's lactate). DO

NOT use non-conductive media (e.g., sterile water, air, gas, glycine).

- Avoid fluid contact with the connection end of the Bipolar PLMe cable connected to the electrosurgical unit.
- The Bipolar PLMe is intended for ablation and/or coagulation purposes only and not for mechanical tissue displacement by applied force, which may bend or separate the electrodes, causing device damage or spacer cracks.
- DO NOT use the Bipolar PLMe as a lever to increase the surgical site or to access tissues, as this could bend or separate the electrodes, causing device damage or spacer cracks.
- The Bipolar PLMe with suction incorporates suction lumens intended to remove bubbles and/or small particles from the surgical field, not for large-volume suction or evacuation.
- Under arthroscopic guidance, ensure that the Bipolar PLMe tip (including the return electrode) is fully surrounded by irrigation solution during use.
- Prevent patient contact with grounded metal objects.
- Excessive electrode use against bony surfaces may cause excessive wear.
- Take fire precautions at all times. Spark and heat production associated with electrosurgery may be flammable.
- DO NOT use flammable anesthetics or oxidizing gases such as nitrous oxide (N₂O) and oxygen.
- As with any other electrosurgical unit, electrodes and cables may offer current paths. Position the cable to avoid patient or other wire contact. Other electrical equipment may experience interference when placed near the system.
- Do not use metal cannulas, as they may damage the Bipolar PLMe insulation and cause burns.
- Position monitoring electrodes as far away as possible from surgical electrodes when high-frequency surgical equipment and physiological monitoring equipment are used simultaneously on the same patient. Needle monitoring electrodes are not recommended.
- Monitoring equipment with high-frequency limiting devices is recommended.
- Ensure adequate flow and saline solution circulation to prevent unnecessary heating, which could cause tissue or thermal injuries.
- Do not use pre-warmed saline solution with the Bipolar PLMe, as it may lead to tissue or thermal injuries.
- For the Bipolar PLMe with suction, properly position the suction adapter to avoid thermal injuries to the physician or patient. Lack of adequate suction may result in thermal injuries.
- For the Bipolar PLMe with suction lumens, using the lumens as the sole drainage route for irrigation is not recommended. It may lead to thermal injuries for the physician and patient. The Bipolar PLMe with suction should be used with a secondary flow output source.
- The product should not be reused or re-sterilized, as it is SINGLE USE ONLY, and should be destroyed and disposed of in HOSPITAL WASTE
- Product is completely safe if used by qualified professionals and respecting the usage specifications described in this document.
- A reduction in the ablation rate may indicate that the PLME-RBW electrode is not functioning correctly. Replace the Bipolar PLMe if the problem persists.
- The Bipolar PLMe is not intended for indefinite soft tissue ablation. Electrode wear during use depends on variables that cannot be replicated in all cases. Factors influencing electrode wear include, but are not limited to, ablation rate, high power settings, prolonged use against bony surfaces, and minimal fluid and suction management.
- Reprocessing this product is prohibited.

MS Registration: 80195520025

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency, guarantees for 24 MONTHS, from the sterilization date, the quality and sterile condition of this product, except in cases of misuse or non-compliance with the guidelines contained in this document.

SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC)
SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC)
CUSTOMER SERVICE (CS)
Phone: +55 11 5013 0010
E-mail: sac@bramsys.com.br
Site: www.bramsys.com.br