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
Reg. ANVISA: 80195520009.

BRAMSYS

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USE INSTRUCTIONS

FACET ABLATION CANNULA KIT

Bramslys ablation cannulas are designed to guide electrodes that produce radiofrequency ablation in various types of surgical procedures.

INDICATIONS FOR USE

The cannulas are specifically indicated for facet ablation procedures in the cervical region at one lateral level. The instructions for use included with the product should be read before opening the sterile packaging to ensure correct and safe usage.

COMPOSITION

Electrode and cannula (quantity as per the table below). All parts and components of the product are sold exclusively for use with this system.

Code	Diameter	Length	Active Tip	Quantity
BMS-D501	18G	50mm	5mm	1
BMS-D533	18G	50mm	10mm	1
BMS-D502	18G	100mm	5mm	1
BMS-D534	18G	100mm	10mm	1
BMS-D503	18G	150mm	5mm	1
BMS-D535	18G	150mm	10mm	1
BMS-D517	21G	50mm	5mm	1
BMS-D549	21G	100mm	10mm	1
BMS-D518	21G	150mm	5mm	1
BMS-D550	21G	50mm	10mm	1
BMS-D519	21G	100mm	5mm	1
BMS-D551	21G	150mm	10mm	1
BMS-D525	22G	50mm	5mm	1
BMS-D526	22G	100mm	5mm	1
BMS-D527	22G	150mm	5mm	1
BMS-D505	18G	50mm	5mm	2
BMS-D537	18G	50mm	10mm	2
BMS-D506	18G	100mm	5mm	2
BMS-D538	18G	100mm	10mm	2
BMS-D507	18G	150mm	5mm	2
BMS-D539	18G	150mm	10mm	2
BMS-D521	21G	50mm	5mm	2
BMS-D553	21G	50mm	10mm	2
BMS-D522	21G	100mm	5mm	2
BMS-D554	21G	100mm	10mm	2
BMS-D523	21G	150mm	5mm	2
BMS-D555	21G	150mm	10mm	2
BMS-D529	22G	50mm	5mm	2
BMS-D530	22G	100mm	5mm	2
BMS-D531	22G	150mm	5mm	2
BMS-D571	21G	50mm	5mm	4
BMS-D572	21G	100mm	5mm	4
BMS-D573	21G	150mm	5mm	4

INSTRUCTIONS FOR USE

- Check if the packaging is undamaged or wet.
- Verify that the sterilization indicator is correct.
- Check the sterilization expiration date.
- Open the packaging and remove the sealed product.
- Following sterile product handling procedures, open the secondary packaging and remove the product.
- Connect the electrode's connection cable to the radiofrequency unit, ensuring it is properly connected to the corresponding power supply.
- Connect the return cable to the radiofrequency unit and the

disposable return pad.

- Apply the return pad to the patient, positioning it on a flat and clean area (thigh, buttocks, or upper back).
- Turn on the radiofrequency unit.
- Adjust the SETUP parameters (time and stimulation) according to the procedure.
- Using the impedance function, check whether the reading corresponds to the expected value for the electrode's location.
- Perform the appropriate stimulation process.
- Carry out the ablation procedure, carefully adjusting the power to achieve the desired temperature.
- If necessary, repeat the steps above according to the adopted protocol for this procedure.
- After the procedure, turn off the radiofrequency unit, disconnect the electrode and return pad, remove the cannula, discard these products, and dispose of them in HOSPITAL WASTE.

PRECAUTIONS, RESTRICTIONS, STORAGE, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use to ensure correct and safe usage and to prevent potential risks.
 - When using rigid return pads (steel plates), always apply conductive gel between the pad's contact surface and the patient, positioning it to maximize contact area.
 - Do not handle or position the electrode on the patient or user while the radiofrequency is activated.
 - Do not use in flammable or explosive atmospheres.
 - Do not use the product beyond its expiration date. Always verify that the sterilization indicator is correct.
 - The product should only be used by qualified medical professionals with the necessary knowledge and experience for proper use, including sterile product handling and, most importantly, maintaining a fully sterile environment in operating rooms. Attention to detail is essential to prevent contamination risks.
 - The product must not be reused or re-sterilized; it is for SINGLE USE • ONLY and must be destroyed and disposed of in HOSPITAL WASTE.
 - The product must always be transported and stored protected from moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
 - The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.
 - Do not store this product near other product categories.
- Check the integrity of the packaging before use. In case of any damage,
- DO NOT USE, as sterility cannot be guaranteed.
 - Reprocessing of this product is strictly prohibited.
 - This product is for single use, sterilized by Ethylene Oxide (ETO), Hydrogen Peroxide, or Gamma Radiation.

Batch No., Sterilization Date, and Expiration Date: See product packaging.

MS Registration No.: 80195520009

Responsible Technician: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for 24 MONTHS, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

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