
Kit Cãnula

Artroscopia ATM

Reg. ANVISA: 80195520037.

BRAMSYS

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ATM ARTHROSCOPY CANNULA KIT

PURPOSE

The ATM Arthroscopy Cannula Kit helps to create access for partial or total removal of joints, soft tissues and cartilage in minimally invasive procedures in the treatment of ATM.

INDICATIONS OF USE

The ATM Arthroscopy Cannula Kit is indicated for minimally invasive surgeries of the temporomandibular joint, allowing the manipulation of the soft and cartilaginous tissues present in the articulation.

COMPOSITION

Primary packaging: Blister or Envelope surgical grade paper containing: 1 ATM ablation electrode, 2 trocar cannulas, 1 hook probe, 1 curette probe, 1 planing probe, 1 sickle cutting probe, 1 U tip probe, 1 blunt plug, 1 cutting plug.

Secondary packaging: individual packaging box, containing primary packaging.

Tertiary packaging: Transport box.

USE INSTRUCTIONS

- Check that the packaging is not damaged or wet;
- Check that the sterilization indication is correct;
- Check the sterilization validity period;
- Open the packaging and remove the packed and sealed product;
- Following the procedures for handling sterile products, open the second package and remove the product;
- Prepare the patient before the procedure according to normal techniques;
- Insert the trocar cannula (s) at the point (s) chosen by the surgeon to access the joint cavity;
- Instrumentals that come with the ATM Arthroscopy Cannula Kit must be selected according to the surgeon's technique and assessment;

USING THE ATM ABLATION ELECTRODE

- Connect the electrode to the electro-surgical unit that must be connected to the corresponding electrical network;
- Turn on the electro-surgical unit;
- Adjust the ablation power (Ablate) to 10W and the coagulation power (coag) to 5W;
- Perform the applicable surgical procedure;
- If necessary, readjust the ablation (ablate) or coagulation (coag) powers to better suit the type of procedure being performed, knowing that the maximum powers to be used are 20W for ablation and 10W for coagulation ;
- After completing the procedure, remove the instruments, turn off the electro-surgical unit, disconnect the electrode, disable it and dispose of it in HOSPITAL WASTE.

The ATM Arthroscopy Cannula Kit is supplied sterile and must be used only once. Do not clean, re-sterilize or reuse the ATM Arthroscopy Cannula Kit, as it may cause malfunction, failure or injury to the patient and may expose the patient to the risk of transmitting infectious diseases.


Refer to the product's Instructions for Use for specific information associated with its use.

USE OF INSTRUMENTS

Electrode for ATM ablation: Performs controlled ablation and tissue modification during surgical procedures of ATM arthroscopy, ATM arthroplasty.

Trocar cannula: Allows the insertion and removal of instruments with 2mm diameters and prevents the flow of fluids, when necessary. It has a skin limiter and graduation.

Cutting plug: In conjunction with the Trocater Cannula, it is used to perforate and break tissue.

Blunt plug: In conjunction with the Trocater cannula, it is used to break the articular disc capsule.

Sickle cutting probe: In conjunction with the Trocater Cannula, it is used to remove adhesion lysis around the disc.

Hook probe: In conjunction with the Trocater Cannula, it pulls the suture threads and assists in the identification of nerves and tendons.

U tip probe: In conjunction with the Trocater Cannula, it holds the articular disc, assisting the Sickle Cutting Probe and the Curette Probe.

Planing probe: In conjunction with the Trocater Cannula, it helps to rough the disc.

Curette probe: In conjunction with the Trocater Cannula, it helps in scraping the cartilage.

CONTRAINDICATIONS

- Non-arthroscopic surgical procedures;
- Patients for whom an arthroscopic procedure is contraindicated for any reason;
- Patients with allergies to contrast, antibiotics or local anesthetics used in the procedure;
- Pregnant women;
- Patients with neurological disorders, preoperative knowledge are a relative contraindication
- Changes in coagulation, for example, during anticoagulant therapy are a relative contraindication;
- Known hypersensitivity to the materials used;

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The user, before using the product, must carefully read the instructions for use that will assist in the correct and safe use in addition to protecting against possible risks;
- If the packaging is damaged, do not use the product;
- Do not use the product after the expiration date, always check if the indication of sterilization is correct;
- The product should only be used by qualified medical professionals, who have the necessary knowledge and experience for the proper use of the product, knowledge of handling sterile products and, mainly, the importance of maintaining a totally sterile environment in the operating rooms.
- Where care and precaution with every detail are essential so that there is no risk of contamination;
- The product must not be reused and / or reprocessed being for SINGLE use, and must be destroyed and disposed of in HOSPITAL WASTE;
- The product must be transported and packed always protected against all types of humidity, dirt and contact with flammable, corrosive, toxic or explosive substances;
- The package must not be punctured or torn, it must be far from sharp or pointed instruments;
- Do not store or store this product near products of different categories;

Totally safe product if used by a qualified professional and respecting the usage specifications described in this document.

Storage temperature: 5°C to 60°C

Storage humidity: 10% to 80%

Attention: It is forbidden to reprocess this product.

This product is for single use, sterilized by Ethylene Oxide (ETO).

This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

Lot number, Sterilization and Validity Date: see the product packing

MS Registration: 80195520037

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillanc, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

www.bramsys.com.br