
# Kit Cãnula Vertebroplastia Full Vert

Reg. ANVISA: 80195520033.

## **BRAMSYS**

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## VERTEBROPLASTY FULL VERT CANNULA KIT

### PURPOSE

The Cannula Vertebroplasty Full Vert Kit is indicated for the application of bone cement in vertebral fractures in order to stabilize them, improving the functional capacity and pain caused by compression.

### INDICATIONS OF USE

The Vertebroplasty Full Vert Cannula Kit aims at minimally invasive percutaneous access to the fractured or damaged vertebra through the working cannula. After accessing the vertebra, it must be filled by applying bone cement. The Kit also allows for the collection of samples for biopsy.

### COMPOSITION

Primary packaging: Blister or Envelope surgical grade paper containing: 4 faceted working cannulas, 1 injector device, 1 biopsy cannula, 1 Guide wire, 1 bone cement mixer.

Secondary packaging: individual packaging box, containing primary packaging.

Tertiary packaging: Transport box.

Code	Diameter	Length	Qty.
VHU12011	11G	120mm	1
VHU15011	11G	150mm	1
VHB12011	11G	120mm	2
VHB15011	11G	150mm	2
VHT12011	11G	120mm	3
VHT15011	11G	150mm	3
VHQ12011	11G	120mm	4
VHQ15011	11G	150mm	4
VHC12011	11G	120mm	5
VHC15011	11G	150mm	5
VHS12011	11G	120mm	6
VHS15011	11G	150mm	6

### USE INSTRUCTIONS

- Define the location for puncture to access the vertebra
- Insert an expansive working cannula to the desired target
- Remove mandrel from faceted work cannula
- Insert the Biopsy Cannula into the faceted cannula.
- Perform circular movements to collect the samples
- Remove the biopsy cannula from the faceted cannula.
- Remove mandrel from the biopsy cannula
- Insert the Guide Wire into the biopsy cannula to remove collected material and deposit it in an appropriate container
- Prepare bone cement with the aid of bone cement mixer.
- Attach the injection device to the bone cement mixer
- Collect bone cement in the injection device
- Uncouple the injector from the bone cement mixer
- Connect the external tube to the injection device
- Fill the outer tube with bone cement through the injector device
- Connect the external tube to the faceted work cannula
- Remove the faceted cannula
- Cover the incision with a sterile razor and suture, if necessary
- After the procedure is finished, the material must be destroyed and disposed of in an appropriate place, following the regulations applicable to the disposal of sharps and contaminated material, adopted by the hospital.

## CONTRAINDICATIONS

- Infection or hematoma in the vicinity of the site chosen for a puncture;
- Pregnancy;
- Known neurological changes in the preoperative period are a relative contraindication;
- Changes in coagulation for example, during anticoagulant therapy are a relative contraindication;
- Inability to obtain the patient's consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast antibiotics and local anesthetics, used in the procedure;
- Any medical, post-surgical, congenital, anatomic or psychological condition that causes the risk of reduced performance of the procedure.

## PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The user, before using the product, must carefully read the instructions for use that will help in the correct and safe use, in addition to protecting against possible risks.
- Do not use the product after the expiration date, always check that the indication of sterilization is correct.
- The product should only be used by qualified medical professionals, with the necessary knowledge and experience for the proper use of the product, knowledge of handling sterile products and, above all, the importance of maintaining a totally sterile environment in the operating rooms. , where care and precautions with all details are essential so that there is no risk of contamination.
- Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
- This product is sterilized by ethylene oxide (ETO).
- The product must be transported and packed always protected against all types of humidity, dirt and contact with flammable, corrosive, toxic or explosive substances.
- The packaging must not be perforated or torn, it must be kept away from sharp or pointed instruments.
- Do not store or store this product near products of different categories.
- Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
- It is prohibited to reprocess this product.
- The product must not be reused and / or re-sterilized, being for SINGLE use only, and must be destroyed and discarded in the HOSPITAL WASTE.

This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

**Lot number, Sterilization and Validity Date:** see the product packing

**MS Registration:** 80195520033

**Technical Responsible:** Marcos Rojas

**CREA:** 5061214879

## WARRANTY CERTIFICATE

**BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

## **CUSTOMER SERVICE (CS)**

Phone: +55 11 5013 0010

E-mail: [sac@bramsys.com.br](mailto:sac@bramsys.com.br)

[www.bramsys.com.br](http://www.bramsys.com.br)