
Kit Cãnula p.Block

Reg. ANVISA: 80195520024.

BRAMSYS

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AV. DR. ALBERTO SAN JUAN, 350

18530-000 – TIETÊ – SP

CNPJ: 04.729.995/0001-87

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P.BLOCK CANNULA

FINALIDADE:

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P.BLOCK CANNULA

PURPOSE:

The p.Block Cannula are designed to perform the procedures plexus block and peripheral nerves.

INDICATIONS OF USE

The p.Block Cannula is a procedure for blocking the plexus and peripheral nerves by means of anesthetic/medication solutions, using electric percutaneous stimulation for target location. These procedures may be performed jointly or separately at the discretion of the physician.

The cannula is introduced into the patient and its location is determined through the use of fluoroscopy, ultrasound or x-ray and through the percutaneous electrical stimulation generates motor stimulation in order to confirm the location of the plexus or peripheral nerves, allowing the accomplishment of Blockade Through the injection of anesthetic fluids / medications.

The usage instructions that are within the product must be read before opening the sterile package, assuring by that a correct and safe use.

COMPOSITION

Primary packaging: Surgical grade paper envelope containing 1 (PBK2110 e PBK2115) or 2 (PBK2110-2 e PBK2115-2) cannulae of the same business model wrapped by a plastic protective tube, 1 surgical planning scale and 1 surgical marker.

Secondary packaging: Surgical grade paper envelope or individual packing box.

Third packaging: Transport box.

USE INSTRUCTIONS

- The p.Block Cannula is a procedure for blocking the plexus and peripheral nerves by means of anesthetic/medication solutions, using electric percutaneous stimulation for target location. These procedures may be performed jointly or separately at the discretion of the physician.

- The cannula is introduced into the patient and its location is determined through the use of fluoroscopy, ultrasound or x-ray and through the percutaneous electrical stimulation generates motor stimulation in order to confirm the location of the plexus or peripheral nerves, allowing the accomplishment of Blockade Through the injection of anesthetic fluids / medications.

- The usage instructions that are within the product must be read before opening the sterile package, assuring by that a correct and safe use.

CONTRAINDICATIONS

- Infection or hematoma near the site chosen for the puncture;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Stimulation: do not perform during defibrillation or cardioversion;
- Changes in coagulation, for example during therapy with anticoagulants are a relative contraindication;
- Lymphangitis, for example, in case of axillary access;
- Post-lymphadenectomy status, for example, in case of axillary access;
- Acute neuritis of the plexus;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The user, before using the product, must carefully read the instructions for use that will help.

correct and safe use, as well as protecting against possible risks.

- Do not use the product after its expiration date, always check if the sterilization indication is correct.

- The product must only be used by qualified medical professionals, who have the knowledge and experience necessary for the proper

use of the product, knowledge of handling sterile products and, mainly, of the importance of maintaining a totally sterile environment in the operating rooms. , where care and precaution with every detail are essential so that there is no risk of contamination.

- Check the integrity of the packaging before using the product, in case of violation, **DO NOT USE**, as the guarantee of its sterility will be compromised.
- This product is sterilized by Ethylene Oxide (ETO).
- The product must be transported and stored, always protected against all types of moisture, dirt and contact with flammable, corrosive, toxic or explosive substances.
- The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.
- Do not store this product next to products of different categories.
- It is prohibited to reprocess this product.
- The product must not be reused and/or re-sterilized, it is for **SINGLE** use, and must be destroyed and disposed of in **HOSPITAL WASTE**.
- The p.Block Cannula must only be used in hospitals and clinics, under aseptic conditions and handled by physicians duly qualified and trained to perform the procedure.

STORAGE CONDITIONS

Storage temperature: 5°C to 60°C

Storage humidity: 10% to 80%

Lot number, Sterilization and Validity Date: see the product packing.

MS Registration: 80195520024

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance under number UY44WX94W79M, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br