
Microdissector Bramslys

Reg. ANVISA: 80195520026.

BRAMSYS

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BRAMSYS MICRODISSECTOR

PURPOSE

The Bramsys Microdissector is a monopolar electrosurgical device designed to perform precision cutting and dissection on soft tissues in various types of surgical procedures.

COMPOSITION

(1) electrode

USE INSTRUCTIONS

- Check that the packaging is not damaged or wet;
- Check that the sterilization indication is correct;
- Check the sterilization expiration date;
- Open the packaging and remove the packaged and sealed product;
- Following the procedures for handling sterile products, open the second package and remove the product;
- Connect the Microdissector to the electro-surgical pen and this to the generating unit that must be connected to the corresponding electrical network;
- Turn on the generating unit;
- Adjust the cutting power (ablate) to 30W and the coagulation power (coag) to 15W;
- Perform the applicable surgical procedure;
- If necessary, readjust the cutting (ablate) or coagulation (coag) powers to better suit the type of procedure being performed, knowing that the maximum powers that should be used are 90W for cutting and 45W for coagulation ;
- After the procedure, turn off the generating unit and disconnect the Microdissector.

CONDITIONS FOR STORAGE AND TRANSPORTATION

To maintain the integrity of the product, it must:

- Always be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt and contact with flammable, corrosive, toxic or explosive substances;
- The packaging must not have holes or tears, and must be kept away from sharp or sharp instruments;
- Do not place under weight or volumes, as this may cause damage to the product's technical characteristics;
- Do not place them next to or next to different categories, as they could cause damage to their technical characteristics;
- They must be packaged and transported following good practices for the transport of sterilized products;
- This product is sterilized by ethylene oxide (ETO) and the packaging must not be damaged before use;
- The product must be transported in its original packaging;
- Storage and transport temperature: 5°C to 60°C;
- Storage and transport humidity: 10% to 80%.

OPERATION CONDITIONS

The user, before using the product, must carefully read the instructions for use, which will assist in correct and safe use, in addition to protecting against possible risks.

- Do not use the product after its expiration date;
- Check that the sterilization indication is correct;
- The product must be used by qualified medical professionals who have the necessary knowledge and experience for the proper use of the product and handling of sterile products;
- Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised;
- The product must be used in a controlled surgical environment;
- Temperature: 10°C to 40°;
- Humidity: 20% to 90%.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The Microdissector electrosurgical electrode is compatible with electrosurgical pens that follow the characteristics described below, which in turn must be compatible with electrosurgical generators.

Compatible electrosurgical pen technical specifications:

- 3-pin connector for manually operated pen;
- Pen slot depth: 27mm
- Electrode rod fitting diameter: 2.40mm
- The user, before using the product, must carefully read the instructions for use, which will assist in correct and safe use, as well as protect against possible risks.
- When using rigid return plates (steel plates), always use conductive gel between the contact surface of the plate and the patient, positioning it in such a way that there is the largest possible contact area between the plate and the patient .
- Do not handle or position the Microdissector on the patient or user when the electrical current is active.
- Do not use in flammable or explosive atmospheres.
- Do not use the product after its expiration date and always check that the sterilization indication is correct.
- The product should only be used by qualified medical professionals, who have the knowledge and experience necessary to properly use the product, knowledge of handling sterile products and, above all, the importance of maintaining a completely sterile environment in operating rooms , where care and precaution with every detail are essential so that there is no risk of contamination.
- The product must not be reused and/or restyled, being for SINGLE use only, and must be destroyed and discarded in HOSPITAL WASTE.
- The product must be transported and stored always protected against all types of humidity, dirt and contact with flammable, corrosive, toxic or explosive substances.
- The packaging must not have holes or tears, and must be kept away from sharp or sharp instruments.
- Do not store or store this product close to products from different categories.
- Keep the product away from sharp or pointed objects.
- Check the integrity of the packaging before using the product. In case of violation, DO NOT USE, as the guarantee of sterility will be compromised.
- Reprocessing this product is prohibited.
- This product is for single use, sterilized by Ethylene Oxide (ETO).
- Totally safe product if used by a qualified professional and respecting the usage specifications described in this document

Batch number, sterilization date and expiry date: see product packaging.

MS Registration: 80195520026

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillanc, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

www.bramsys.com.br

Item	Compr. / Length	Formato / Format	Ângulo da ponta ativa / Ángulo de punta activo / Active tip angle
BMS-D444	20mm	R	180°
BMS-D445	30mm	R	180°
BMS-D446	40mm	R	180°
BMS-D447	70mm	R	180°
BMS-D448	90mm	R	180°
BMS-D449	120mm	R	180°
BMS-D450	140mm	R	180°
BMS-D451	160mm	R	180°
BMS-D452	30mm	R	45°
BMS-D453	40mm	R	45°
BMS-D454	30mm	R	45°
BMS-D455	30mm	R	90°
BMS-D456	70mm	R	45°
BMS-D457	90mm	R	45°
BMS-D458	120mm	R	45°
BMS-D459	140mm	R	45°
BMS-D460	170mm	R	45°
BMS-D461	170mm	R	45°
BMS-D462	140mm	R	45°
BMS-D463	120mm	R	45°
BMS-D464	90mm	R	45°
BMS-D465	70mm	R	45°
BMS-D466	170mm	R	90°
BMS-D467	140mm	R	90°
BMS-D468	120mm	R	90°
BMS-D469	90mm	R	90°
BMS-D470	70mm	R	90°
BMS-D471	150mm	CV	180°
BMS-D472	130mm	CV	180°
BMS-D473	110mm	CV	180°
BMS-D474	90mm	CV	180°
BMS-D475	160mm	CT	45°
BMS-D476	140mm	CT	45°
BMS-D477	150mm	AG	180°
BMS-D478	120mm	AG	180°
BMS-D479	100mm	AG	180°
BMS-D480	70mm	AG	180°

LEGENDA | LEYENDA | LABEL

R - Reto | Recto | Straight

CV - Curvo | Curved

CT - Contorno | Contour

AG - Angulado | Angular | Angulado