
# Kit Cânula

## Cifoplastia d.Kypho

Reg. ANVISA: 80195520031.

# BRAMSYS

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## USE INSTRUCTIONS

#### D.KYPHO CYPHOPLASTY CANNULA KIT

The d.Kypho Kyphoplasty Cannula Kit is intended for minimally invasive percutaneous access to the fractured or damaged vertebra through the working cannula. After accessing the vertebra, the balloon is used to restore the height of the vertebral body, forming a cavity to be filled through the application of bone cement. The Kit also allows the collection of samples for biopsy.

#### INDICATIONS FOR USE

The d.Kypho Kyphoplasty Cannula Kit is indicated for minimally invasive kyphoplasty procedures in the treatment of spinal fractures or injuries.

#### COMPOSITION

Primary packaging: Blister or surgical grade paper envelope containing all items that make up the kit, as shown below.

Secondary packaging: Individual packaging box, containing primary packaging and Inflator Device.

Tertiary packaging: Transport box.

Item	Ø	Lenght	Quantity. DKU	Quantity. DKB	Capacity
Work expanding cannula	11G	120mm	1	2	-

Biopsy cannula	15G	185mm	1	1	
Syringe with wings	-	-	4	4	2,5ml
Drill	13G	187mm	1	1	-
Disp. Insufflator w/ adapt. 3 way	-	-	1	2	20ml
Balloon catheter	-	10, 15 ou 20mm	1	2	450psi
Bone cement application cannula	13G	187mm	3	6	-
Bone cement mixer	-	-	-	-	40g

#### INFLATION DEVICE SPECIFICATIONS

- Operating range: -0.4 ATM to 30 ATM (1ATM = 1BAR = 14.7 PSI);
- Accuracy:  $\pm 3.0\%$  of normal full scale;
- Fluid administered: 0.45ml ( $\pm 0.07$ ml of fluid administered) for every 360° clockwise of the syringe plunger handle;
- Liquid temperature: 10°C to 40°C;
- Humidity and operating range: 20% to 90%;

#### INSTRUCTIONS FOR USE

1. Define the puncture site for accessing the vertebra
2. Insert working expansive cannula to desired target
3. Remove work cannula mandrel
4. Insert the Drill into the expanding working cannula
5. Make circular movements to drill the vertebrae, in order to prepare it to insert the Balloon Catheter.
6. Remove the Drill from the working expanding cannula
7. Insert the Biopsy Cannula into the working Expansion Cannula
8. Make circular movements to collect the samples
9. Remove the biopsy cannula from the working expanding cannula
10. Remove the biopsy cannula mandrel
11. Insert the Guide Wire into the Biopsy Cannula to remove collected material and deposit it in an appropriate container
12. Collect 10ml of Contrast using a syringe
13. Connect the syringe to the 3-way adapter
14. Position the 3-way adapter lock to inject the contrast
15. Connect the 3-way adapter to the Inflator Device
16. Unlock the inflator device plunger
17. Inject the contrast into the Inflator Device
18. Uncouple the syringe from the 3-way adapter
19. Position the 3-way Adapter lock to use the Balloon Catheter
20. Remove the air between the plunger of the Inflator Device and the outlet of the 3-way adapter
21. Locking the Inflator Device Plunger
22. Insert the Balloon Catheter into the working expanding cannula until the positioning mark is completely inserted into it
23. Perform positioning adjustment, if necessary, until the radiopaque marks of the Balloon Catheter are outside the work expansive Cannula
24. Remove mandrel from Balloon Catheter
25. Connect 3-way Adapter to the Balloon Catheter access
26. Rotate the Plunger of the Inflator Device clockwise to inject the contrast into the Balloon Catheter and inflate it in order to prepare the cavity in the vertebral body
27. Prepare the bone cement using the bone cement mixer
28. Attach the syringe with handles to the bone cement mixer
29. Collect the bone cement in the syringe with handles
30. Uncouple the syringe with wings from the bone cement mixer
31. Connect the syringe with handles to the Cannula for applying bone cement
32. Inject the bone cement into the Cannula for application of bone cement until it is completely filled
33. Disconnect the winged syringe from the Cannula for application of bone cement
34. Rotate the Inflator Device plunger counterclockwise until all contrast is removed from the Balloon Catheter.
35. Unlock the Inflator Device plunger and pull it to perform negative pressure on the Balloon Catheter and ensure that it is completely empty.
36. Locking the Inflator Device Plunger
37. Remove the Balloon Catheter from the Working Cannula
38. Insert the Cannula for application of bone cement in the work expansive Cannula

39. Insert Mandrel into Bone Cement Delivery Cannula to transfer cement to vertebra
40. Remove the Bone Cement Cannula from the Expansive Work Cannula
41. Remove the working expansive cannula
42. Cover the incision with sterile gauze and suture, if necessary
43. Dispose of all material used in the procedure in an appropriate place, following the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.

#### **CONTRAINDICATIONS**

- Infection or hematoma near the chosen puncture site
- Pregnancy;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Changes in coagulation, for example during anticoagulant therapy are a relative contraindication;
- Failure to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast, antibiotics and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical or psychological condition that increases the risk of reduced performance of the procedure.

#### **PRECAUTIONS, RESTRICTIONS, STORAGE AND WARNINGS**

- The user, before using the product, must carefully read the instructions for use that will help him to use it correctly and safely, in addition to protecting against possible risks.
- Do not use the product beyond its expiration date, always check that the sterilization indication is correct.
- The product should only be used by qualified medical professionals, who have the knowledge and experience necessary for the proper use of the product, knowledge of handling sterile products and, mainly, the importance of maintaining a completely sterile environment in the operating rooms, where care and precaution with every detail are essential so that there is no risk of contamination.
- Check the integrity of the package before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
- This product is for single use, sterilized by Ethylene Oxide (ETO).
- It is prohibited to reprocess this product.
- The product must always be transported and stored protected against excessive humidity, high temperatures, exposure to sunlight, dirt and contact with flammable, corrosive, toxic or explosive substances and away from sharp or pointed instruments.
- Do not place next to or together with different categories, as they could damage its technical characteristics and compromise its integrity / quality.
- No modifications are allowed to this product.
- Disposal of the product after use must follow the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.
- Temperature for storage and transport: 5°C to 60°C
- Humidity for storage and transport: 10% to 80%

**No. Batch Number, Sterilization Date and Expiration Date:** see product packaging.

**ANVISA:** 80195520031

**Technical Manager:** Marcos Rojas

**CREA:** 5061214879

#### **WARRANTY CERTIFICATE**

**BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

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