
Kit Cânula

Cifoplastia t.Kypho

Reg. ANVISA: 80195520030.

BRAMSYS

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USE INSTRUCTIONS

T.KYPHO CYPHOPLASTY CANNULA KIT

The t.Kypho Kyphoplasty Cannula Kit is intended for minimally invasive percutaneous access to the fractured or damaged vertebra through the working cannula. After accessing the vertebra, the balloon is used to restore the height of the vertebral body, forming a cavity to be filled through the application of bone cement. The Kit also allows the collection of samples for biopsy.

INDICATIONS FOR USE

The t.Kypho Kyphoplasty Cannula Kit is indicated for minimally invasive kyphoplasty procedures in the treatment of spinal fractures or injuries.

COMPOSITION

Primary packaging: Blister or surgical grade paper envelope containing all items that make up the kit, as shown below.

Secondary packaging: Individual packaging box, containing primary packaging and Inflator Device.

Tertiary packaging: Transport box.

Item	Ø	Lenght	Quantity. TKU	Quantity. TKB	Capacity
Percutaneous access cannula	11G	100mm	1	2	-
Straight and faceted guide wire	13G	280mm	2	4	-
Work expanding cannula	8G	120mm	1	2	-
Biopsy cannula	15G	185mm	1	1	
Syringe with wings	-	-	4	4	2,5ml
Drill	10G	187mm	1	1	-
Disp. Insufflator w/ adapt. 3 way	-	-	1	2	20ml
Balloon catheter	-	10, 15 ou 20mm	1	2	450psi
Bone cement application cannula	9G	187mm	3	6	-
Bone cement mixer	-	-	-	-	40g

INFLATION DEVICE SPECIFICATIONS

- Operating range: -0.4 ATM to 30 ATM (1ATM = 1BAR = 14.7 PSI);
- Accuracy: $\pm 3.0\%$ of normal full scale;
- Fluid administered: 0.45ml (± 0.07 ml of fluid administered) for every 360° clockwise of the syringe plunger handle;
- Liquid temperature: 10°C to 40°C;
- Humidity and operating range: 20% to 90%;

INSTRUCTIONS FOR USE

1. Define the puncture site for accessing the vertebra;
2. Insert Percutaneous Access Cannula to the desired target;
3. Remove mandrel from the Percutaneous Access Cannula;
4. Insert guidewire into percutaneous access cannula;
5. Remove Percutaneous Access Cannula;
6. Insert Expanding Working Cannula through the Guide Wire;
7. Remove guidewire;
8. Remove the working cannula's mandrel;
9. Insert the Drill into the expansive working cannula;
10. Make circular movements to perforate the vertebra, in order to prepare it to insert the Balloon Catheter;
11. Remove Drill from the work expanding cannula;
12. Insert biopsy cannula into the working expanding cannula;
13. Make circular movements to collect the samples;
14. Remove biopsy cannula from working expanding cannula;
15. Remove the biopsy cannula mandrel;
16. Insert the guidewire into the biopsy cannula to remove the collected material and deposit it in an appropriate container;
17. Collect 10ml of Contrast using a Syringe;
18. Connect the Syringe to the 3-way Adapter;
19. Position the 3-way Adapter lock to inject the contrast;
20. Connect the 3-way Adapter to the Inflator Device;
21. Unlock the Inflator Device Plunger;
22. Inject the contrast into the inflator device;
23. Uncouple the Syringe from the 3-way Adapter;
24. Position the 3-way Adapter lock to use the Balloon Catheter;
25. Remove the air between the plunger of the Inflator Device and the outlet of the 3-way Adapter;
26. Lock the Inflator Device plunger;
27. Insert Balloon Catheter in the work expansive Cannula until the positioning mark is completely inserted in it;
28. Perform positioning adjustment, if necessary, until the radiopaque marks of the Balloon Catheter are outside the working expansive Cannula;
29. Remove the mandrel from the Balloon Catheter;
30. Connect the 3-way Adapter to the Balloon Catheter access;
31. Rotate the plunger of the Inflator Device clockwise to inject the contrast into the Balloon Catheter and inflate it in order to prepare the cavity in the vertebral body;
32. Prepare the bone cement using the bone cement mixer;
33. Attach Syringe with handles to Bone Cement Mixer;
34. Collect bone cement in the Syringe with handles;
35. Uncouple Syringe with wings from Bone Cement Mixer;
36. Connect the Syringe with handles to the Cannula for applying bone cement;
37. Inject bone cement into the Cannula for bone cement application until it is completely filled;
38. Disconnect the Syringe with handles from the Cannula for application of bone cement;
39. Rotate the plunger of the Inflator Device counterclockwise until all the contrast is removed from the Balloon Catheter;

40. Unlock the plunger of the inflator device and pull it to perform negative pressure on the balloon catheter and ensure that it is completely empty;
41. Lock the plunger of the inflator device;
42. Remove Balloon Catheter from the Working Cannula;
43. Insert Cannula for application of bone cement in the work expansive Cannula;
44. Insert mandrel into Cannula for application of bone cement to transfer the cement to the vertebra;
45. Remove cannula for application of bone cement from the work expansive cannula;
46. Remove working expansive cannula;
47. Cover the incision with sterile gauze and suture, if necessary;
48. Discard all material used in the procedure in an appropriate place

CONTRAINDICATIONS

- Infection or hematoma near the chosen puncture site
- Pregnancy;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Changes in coagulation, for example during anticoagulant therapy are a relative contraindication;
- Failure to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast, antibiotics and/or local anesthetics used in the procedure;
- Any medical, post-surgical, congenital, anatomical or psychological condition that increases the risk of reduced performance of the procedure.

PRECAUTIONS, RESTRICTIONS, STORAGE AND WARNINGS

- The user, before using the product, must carefully read the instructions for use that will help him to use it correctly and safely, in addition to protecting against possible risks.
- Do not use the product beyond its expiration date, always check that the sterilization indication is correct.
- The product should only be used by qualified medical professionals, who have the knowledge and experience necessary for the proper use of the product, knowledge of handling sterile products and, mainly, the importance of maintaining a completely sterile environment in the operating rooms, where care and precaution with every detail are essential so that there is no risk of contamination.
- Check the integrity of the package before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
- This product is for single use, sterilized by Ethylene Oxide (ETO).
- It is prohibited to reprocess this product.
- The product must always be transported and stored protected against excessive humidity, high temperatures, exposure to sunlight, dirt and contact with flammable, corrosive, toxic or explosive substances and away from sharp or pointed instruments.
- Do not place next to or together with different categories, as they could damage its technical characteristics and compromise its integrity / quality.
- No modifications are allowed to this product.
- Disposal of the product after use must follow the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.
- Temperature for storage and transport: 5°C to 60°C
- Humidity for storage and transport: 10% to 80%

No. Batch Number, Sterilization Date and Expiration Date: see product packaging.

ANVISA: 80195520030

Technical Manager: Marcos Rojas | **CREA:** 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance under number UY44WX94W79M, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br