
# Kit Cãnula HealBone

Reg. ANVISA: 80195529010.

## **BRAMSYS**

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### **HEALBONE CANNULA KIT**

#### **PURPOSE / INDICATION FOR USE**

The HealBone Cannula Kit is a medical device for treating conditions such as osteoarthritis, patellar chondromalacia, osteochondral lesion, osteonecrosis, pseudoarthrosis, tendinopathy, subchondral lesion, bone marrow edema, or insufficiency fracture, among other degenerative or inflammatory conditions that result in chronic pain for the patient. The product is designed for the collection, preparation, and administration of bone graft in a minimally invasive manner, providing greater precision and safety to the procedure. This solution can be used in therapeutic procedures aimed at pain relief and functional improvement, being applied to areas previously identified as sites of inflammation or degeneration, according to clinical indication.

#### **PRINCIPLE OF OPERATION / MECHANISM OF ACTION**

Functionality by Item:

##### **1. Aspiration Cannula**

Diamond tip to facilitate puncture

Fixed tip to prevent contamination by peripheral blood

Alternating side holes to provide greater variety during collection  
 11G gauge for efficient collection and minimally invasive access  
 Centimeter markings to facilitate cannula retraction and collection at different levels  
 Luer lock connection to ensure a secure seal during collection

## 2. Depth Control Device

Assists in retracting the aspiration cannula for collection at different levels

## 3. Filter

Removes bone fragments and other debris.

## 4. 3-Way Adapter

Allows safe connection of devices for handling and collection  
 Three-way connections with luer lock to ensure a secure seal

## 5. Extension Tube

Provides ergonomics during the use of the application cannula  
 Flexible and transparent tube for visualization during application  
 Luer lock connections at both ends to ensure a secure seal

## 6. Injection Cannula

Application of the solution at the lesion site or treatment area  
 Sharp tip to facilitate tissue penetration  
 18G or 20G gauge for application  
 Length of 100 mm or 150 mm  
 Luer lock connection to ensure a secure seal  
 Echogenic tip for ultrasound visualization  
 Centimeter markings

## 7. Planning Scale and Surgical Marker

Assists in marking and planning the desired target during the procedure.

## COMPOSITION

(1) Aspiration cannula, (1) filter, (1) 3-way adapter, (1) depth control device, (1) extension tube, (1) application cannula, (1) surgical marker, (1) surgical planning scale.

## COMMERCIAL PRESENTATION

Primary packaging: Sealed blister with Tyvek paper containing: (1) Aspiration Cannula, (1) Filter, (1) 3-Way Adapter, (1) Depth Control Device, (1) Extension Tube, (1) Application Cannula, (1) Surgical Marker, (1) Surgical Planning Scale.

Secondary packaging: Individual packaging box containing the primary packaging.

Code	Injection Cannula			Aspiration Cannula
	Length (mm)	Diameter (G)	Qty.	Qty.
KCHB-18100	100	18	1	1
KCHB-18150	150	18	1	1
KCHB-20100	100	20	1	1
KCHB-20150	150	20	1	1

## USE INSTRUCTIONS

- Ensure that the packaging is not damaged.
- Remove all items from the sterile packaging in a clean environment using aseptic techniques.
- Ensure that all items are intact.
- Position the patient according to the surgical procedure to be performed.  
 Note: If necessary, use the Surgical Planning Scale and Surgical Marker to assist with puncture planning and locating the desired target.
- Prepare the skin puncture area using aseptic and antiseptic techniques.
- Perform the percutaneous puncture using the Aspiration cannula.  
 Note: If necessary, use fluoroscopy, ultrasound, or X-ray to assist with puncture and target monitoring.
- Reach the deepest point to begin the collection.
- Using the Depth Control Device, retract the cannula to perform additional collections as needed.
- Remove the cannula after completing the collection.  
 Note: If necessary, use the filter to remove bone fragments and other debris, and the 3-way adapter to prepare the collected substance before administration.
- Position the patient according to the desired application area.
- Prepare the skin puncture area using aseptic and antiseptic techniques.
- Perform the percutaneous puncture using the injection cannula.  
 Note: If necessary, use the extension tube to apply the collected material.

- Remove the cannula after application.
- After the procedure is complete, all materials should be decharacterized and disposed of in an appropriate location, following applicable regulations for the disposal of sharp and contaminated materials adopted by the hospital.

#### **STORAGE AND TRANSPORT CONDITIONS**

To maintain the integrity of the product, it should:

Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging should not be punctured or torn and should be kept away from cutting or sharp instruments.

Do not place under weight or pressure, as it may damage the product's technical characteristics.

Avoid placing it next to or with different categories of products, as this could damage its technical characteristics.

Be handled and transported following best practices for transporting sterilized products.

This product is sterilized with ethylene oxide (ETO) and must not have damaged packaging before use.

The product should be transported in its original packaging.

Storage temperature: Between 15°C and 30°C.

#### **CONTRAINDICATIONS**

It is not recommended for use in patients receiving intensive anticoagulant treatment or those with severe hemorrhagic conditions. Medical advice is essential when considering performing a biopsy in patients with bleeding disorders or those on anticoagulant therapy.

#### **ADVERSE EFFECTS**

Adverse effects may include, but are not limited to: discomfort, bruising, infection, bleeding, damage to adjacent tissue, and blood vessel perforation.

#### **PRECAUTIONS, RESTRICTIONS AND WARNINGS**

Avoid using the product if the packaging shows damage.

Do not use the product if it has passed its expiration date.

This product is intended for single use only and must be discarded after use. Reuse is strictly prohibited.

Reuse or re-sterilization may compromise the product's structural integrity and/or cause damage that can result in injury or illness to the patient.

Reuse or re-sterilization may also present a risk of contamination and/or cause infection in the patient or cross-contamination, including but not limited to the transmission of infectious diseases between patients.

It is advisable to avoid using the product if particularly dangerous maneuvers are required, or when subjected to twisting or bending forces that could compromise its integrity and proper functioning.

Post-collection care may vary depending on the technique used and the physiological conditions of each patient.

A thorough clinical judgment is required when considering performing a collection in patients undergoing anticoagulant therapy or with bleeding disorders.

This product should be used by a physician who is fully familiar with the

#### **WARRANTY CERTIFICATE**

The company BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency, guarantees for 24 MONTHS, from the sterilization date, the quality and sterile condition of this product, except in cases of misuse or non-compliance with the guidelines contained in this document.

#### **SERVIÇO DE ATENDIMENTO AO CLIENTE (SAC) SERVICIO DE ATENDIMIENTO AL CLIENTE (SAC) CUSTOMER SERVICE (CS)**

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