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
Site: www.bramsys.com.br

USE INSTRUCTIONS

TRIGEMINAL BALLOON COMPRESSION CANNULA KIT 14G-100

The Trigeminal Balloon Compression Cannula Kit 14G-100 is specifically indicated for the treatment of trigeminal neuralgia through the technique of compressing the trigeminal root with an inflatable balloon.

INDICATIONS FOR USE

The balloon is inserted through a cannula in order to be positioned at the desired target. Through the injection of contrast liquid, the balloon is inflated to compress the trigeminal nerve.

COMPOSITION

Primary packaging: sealed blister with Tyvek containing: 1 cannula, 1 mandrel with blunt tip, 1 mandrel with faceted tip, 1 balloon of 1.5mm/4.0F, 1 syringe of 1 ml, and 1 top limiter.

Secondary packaging: individual packaging box, containing primary packaging.

Tertiary packaging: transport box.

Code	Ø	Length
BMS-D156	14G	100mm

INSTRUCTIONS FOR USE

- Packaging opening:
 - 1.1 Check if the packaging is not damaged or wet.
 - 1.2 Verify that the sterilization indication is correct.
 - 1.3 Check the expiration date of sterilization.
 - 1.4 Open the packaging box and remove the sealed blister.
 - 1.5 Following sterile product handling procedures, open the blister and remove the product.
- Balloon inspection:
 - 2.1 Remove the cap and the guide from the catheter's end.
 - 2.2 Inflate the balloon with 0.75 ml of contrast agent.
 - 2.3 Create a vacuum in the syringe.
 - 2.4 Repeat operations 2 and 3 until all air is removed.
 - 2.5 Remove the air from the syringe.
 - 2.6 Fill the syringe until it reaches 0.75 ml of contrast agent.
 - 2.7 Connect the syringe to the balloon catheter and inflate it with the contrast agent until it reaches 0.75 ml, which is the recommended maximum inflation volume. If the balloon does not inflate, has a leak, or inflates asymmetrically (eccentrically), do not use it.
- Adjustment of the top limiter:
 - 3.1 Insert the top limiter into the catheter.
 - 3.2 Insert the balloon into the cannula until the catheter's tip is in the desired position for surgery.
 - 3.3 Lock the top limiter onto the catheter by tightening the screw.
 - 3.4 Remove the catheter from the cannula.
- Cannula introduction:
 - 4.1 Choose the mandrel to be used (blunt or faceted) and insert the mandrel into the cannula.
 - 4.2 Puncture the cheek with the mandrel inside the cannula until reaching the foramen oval without penetrating it.
 - 4.3 Remove the mandrel from the cannula.
- Introduction of the balloon-catheter for trigeminal nerve compression:
 - 5.1 Remove the mandrel from the catheter.
 - 5.2 Connect the syringe with 0.75 ml of contrast agent to the catheter through the luer connection.
 - 5.3 Inflate and deflate the contrast agent from the balloon several times, with the balloon's tip facing downward, until all air is removed from the balloon.
 - 5.4 Insert the catheter into the cannula until the top limiter fits into the cannula.
 - 5.5 Slowly inflate the balloon to compress the trigeminal nerve. The balloon should take on a pear-like shape.
 - 5.6 Keep the balloon inflated, compressing the nerve for 40 to 60 seconds.
- Completion of the surgical procedure:
 - 6.1 Empty the balloon by withdrawing the contrast agent with the syringe.
 - 6.2 Remove the balloon along with the cannula as a single unit.

6.3 After the procedure, dispose of all parts in HOSPITAL WASTE.

CONTRAINDICATIONS

- Allergic to natural latex rubber.

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.

Do not place under weight or heavy volumes, as this may damage the product's technical characteristics.

Do not place near or together with different categories, as this may compromise its technical characteristics.

It must be stored and transported following best practices for handling sterilized products.

This product is sterilized with ethylene oxide (ETO) and the packaging must not be damaged before use.

The product must be transported in its original packaging.

Storage and transportation temperature: Between 15°C and 30°C.

Storage and transportation humidity: Between 30% and 80%.

OPERATING CONDITIONS

The user, before using the product, must carefully read the instructions for use that will assist in its correct and safe use, in addition to protecting against possible risks.

Do not use the product past its expiration date.

Ensure that the sterilization indication is correct.

The product should be used by qualified medical professionals with the necessary knowledge and experience for the proper use of the product and handling of sterile products.

Verify the integrity of the packaging before using the product; in case of violation, DO NOT USE it, as the guarantee of its sterility will be compromised.

PRECAUTIONS, RESTRICTIONS, STORAGE AND WARNINGS

- The correct and safe use of this product for its intended purpose is governed by the conditions described in this user manual.

- The product must not be reused and/or reprocessed, as it is for SINGLE USE only and must be destroyed and discarded as MEDICAL WASTE.

- This product is sterilized using Ethylene Oxide (ETO).

- No modifications to this product are allowed.

- Disposal must comply with the requirements of RDC No. 222, dated March 28, 2018, and other applicable laws and regulations.

No. Batch Number, Sterilization Date and Expiration Date: see product packaging.

ANVISA: 80195520015

Technical Manager: Marcos Rojas | **CREA:** 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance under number UY44WX94W79M, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsyst.com.br