
Eletrodo Bipolar PLMe

Reg. ANVISA: 80195520025.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) – is in accordance with RDC No. 751/2022. A PDF reader is required, and the program can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the product purchased.

To obtain the instruction for use in printed format, contact us via email at sac@bramsys.com.br. There is no additional cost for sending the product instruction for use in printed or digital format.

This instruction for use is applicable to all commercial models of the Bipolar PLMe family.

USE INSTRUCTIONS

ELECTRODE BIPOLAR PLMe

INDICATION FOR USE / PURPOSE

The Bipolar PLMe is intended to produce ablation and coagulation of soft tissues and for hemostasis of blood vessels in minimally invasive, video-guided procedures in joints, including the spine.

COMPOSITION

(1) Electrode

COMMERCIAL MODELS

Code	Tip model	Drive
PLME-BTW	Bolt	Pedal
PLME-IFW	Infinity	Pedal
PLME-SBW	Super Bolt	Pedal
PLME-STW	Star	Pedal
PLME-XGW	Xgon	Pedal

The commercial models are supplied in the following presentation:

- Primary packaging: PET blister (polyethylene terephthalate) sealed with Tyvek paper, permeable to ethylene oxide for sterilization, containing (1) electrode along with five (5) traceability labels.
- Secondary packaging: Individual packaging box made of duplex paper containing the primary packaging.
- Tertiary packaging: Transport box.

Each package contains only one (1) commercial model of the Bipolar PLMe.

INSTRUCTIONS FOR USE

- Check that the packaging is not damaged or wet;
- Verify that the sterilization indicator is correct;
- Check the sterilization expiration date;
- Open the packaging and remove the sealed product;
- Following sterile product handling procedures, open the secondary packaging and remove the product;
- Prepare the patient prior to the procedure according to standard protocols;
- To use the suction function, attach the suction line adapter to the standard hospital suction equipment. Ensure that the flow regulator is open and that the vacuum is within the range of 200 mmHg (7.87 in Hg) to 400 mmHg (15.75 in Hg);
- Connect the Bipolar PLMe cable end to the generator unit receptacle;
- Insert the Bipolar PLMe through the trocar and place it in contact with the tissue to be ablated;
- Select the ablation function using the foot pedal. The electrodes must always remain in contact with the tissue to be ablated and must be in motion when active. When the Bipolar PLMe is activated, the generator unit emits a single-tone audio signal. Some bubbling and an orange glow on the Bipolar PLMe are normal. If a different setting is desired, adjust the configuration using the foot pedal or the control buttons on the front panel of the generator unit until the desired setting is reached;
- Periodically check the electrodes to ensure they remain intact and that the Bipolar PLMe is functioning correctly. Use the recommended preset point to achieve the desired result;
- For coagulation, select the coagulation function using the foot pedal. If a different setting is desired, adjust it using the control buttons on the front panel of the generator unit until the desired setting is reached;
- After the procedure, turn off the generator unit, disconnect the suction tubing, disconnect the electrode, render it unusable, and dispose of it in HOSPITAL WASTE.

CONTRAINDICATIONS

- Non-arthroscopic surgical procedures;
- Procedures without the use of a conductive irrigant;
- Patients for whom an arthroscopic procedure is contraindicated for any reason;
- Patients with cardiac pacemakers or other implanted electronic devices.

ADVERSE EFFECTS

- Failure to use the suction adapter or provide adequate suction may cause thermal injuries to the physician or patient;
- During the procedure, ensure adequate fluid flow and circulation of saline solution to prevent unnecessary heating, which could cause tissue damage;
- As a consequence of electrosurgery, surrounding tissues may be damaged due to iatrogenic injuries.

STORAGE AND TRANSPORT CONDITIONS

- The product must always be transported and stored protected from any

moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances;

- The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments;
- Do not place weights or other items on top, as this may damage the product's technical characteristics;
- Do not store next to or with different product categories, as this may affect their technical characteristics;
- Products should be stored and transported following good practices for the handling of sterilized items;
- Check the integrity of the packaging before use. If the packaging is compromised, DO NOT USE, as sterility cannot be guaranteed.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- The Bipolar PLMe is designed for use with generator units that are compatible with the following technical specifications:
 - Operating frequency: 100–450 kHz
 - Operating voltage: 0–320 Vrms
 - Operating power: 400 W max
- Before using the product, the user must carefully read the instructions for use, which will guide correct and safe usage and protect against potential risks;
- Do not use the product if the packaging is damaged;
- The product should only be used by qualified medical professionals who have the necessary knowledge and experience to use the product properly, understand sterile product handling, and, most importantly, appreciate the importance of maintaining a completely sterile environment in the operating room, where attention to detail is critical to prevent contamination;
- Do not touch the active tip of the Bipolar PLMe when it is activated.
- Use only conductive media (e.g., normal saline or Ringer's lactate, etc.). Do NOT use non-conductive media (e.g., sterilized water, air, gas, glycine, etc.);
- Do not allow fluids to contact the connection end of the Bipolar PLMe cable that links to the electrosurgical unit;
- The Bipolar PLMe is designed for ablation and/or coagulation purposes and not for mechanically displacing tissue by applying force. Doing so may bend the electrodes or even separate them, causing damage to the device and/or cracks in the spacer;
- Do NOT use the Bipolar PLMe as a lever to enlarge the surgical site or to gain access to tissue, as this could bend or separate the electrodes, causing damage to the device and/or cracks in the spacer;
- The Bipolar PLMe with suction incorporates suction lumens, designed to remove bubbles and/or small particles from the surgical field and not for suctioning or evacuating large volumes from the surgical field;
- Under arthroscopic guidance, ensure that the tip of the Bipolar PLMe (including the return electrode) is completely surrounded by irrigation solution during use;
- Do not allow the patient to come into contact with metal objects connected to ground;
- Excessive use of the electrodes against bony surfaces may result in excessive wear;
- Always observe fire safety precautions. Sparks and heating associated with electrosurgery can be a fire hazard;
- Do NOT use flammable anesthetics or oxidizing gases such as nitrous oxide (N₂O) or oxygen;
- As with any electrosurgical unit, electrodes and cables may provide pathways for electrical current. Position the cable to avoid contact with the patient or other wires. Other electrical equipment may experience interference if positioned near the system;
- Do not use metal cannulas, as they may damage the insulation of the Bipolar PLMe and cause burns.
- Monitoring electrodes should be positioned as far as possible from the surgical electrodes when high-frequency surgical equipment and physiological monitoring equipment are used simultaneously on the same patient. Needle-type monitoring electrodes are not recommended;
- It is recommended to use monitoring equipment that integrates high-frequency limitation devices;
- Ensure adequate flow and circulation of saline solution to prevent unnecessary heating, which could cause tissue damage or thermal injuries;
- Do not use preheated saline solution with the Bipolar PLMe, as this may cause tissue damage or thermal injuries;
- For the Bipolar with suction, the suction adapter must be correctly positioned to prevent thermal injuries to the physician or patient. Failure to provide adequate suction may result in thermal injuries;
- For the Bipolar PLMe with suction lumens, using the lumens as the sole irrigation outlet is not recommended, as it may cause thermal injuries to the physician and patient. The Bipolar PLMe with suction should be used together with a secondary flow source;
- The product must not be reused or resterilized, as it is single-use only and must be destroyed and disposed of in HOSPITAL WASTE;

- The product is completely safe when used by a qualified professional and according to the usage specifications described in this document;
 - A decrease in ablation rate may indicate that the Bipolar PLMe is not functioning properly. Replace the Bipolar PLMe if the problem persists;
 - The Bipolar PLMe is not intended for continuous ablation of soft tissues.
- Electrode wear during use and the rate of wear depend on variables that cannot be replicated in all cases. Factors affecting electrode wear include, but are not limited to, ablation rate, high power settings, prolonged use against bony surfaces, and minimal fluid and suction management;
- Reprocessing this product is prohibited.

ANVISA Registration No.: 80195520025

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br