

Kit Cânula Safe Block

Reg. ANVISA: 80195520028.

BRAMSYS

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CANNULA KIT SAFE BLOCK

FINALIDAD

The Safe Block Cannula Kit is indicated for performing plexus and peripheral nerve block procedures using anesthetic solutions/medications, using percutaneous electrical stimulation to locate the target.

INDICATIONS OF USE:

- The kit is indicated for carrying out procedures to block the plexus and peripheral nerves using anesthetic solutions/medicines, using percutaneous electrical stimulation to locate the target. These procedures can be performed together or separately, at the doctor's discretion.
- The cannula is introduced into the patient and its location is determined through the use of fluoroscopy, ultrasonography or radiography and, through percutaneous electrical stimulation, generates motor stimulation to confirm the location of the plexus or the peripheral nerves, allowing blockage through the injection of anesthetic fluids / medications.
- For a safer location, the radiopaque adhesive grade is used, where it is positioned and adhered to the patient's skin in the area of the surgical incision in the preoperative period and in the asepsis time. The radiopaque adhesive is shown in the fluoroscopy, ultrasound or X-ray image, guiding the site to access the required surgical point, acting as a skin reference point in the radiological image obtained, where reference points are used to identify surgery at the incision site.
- The instructions for use must be read before opening the sterile packaging, thus ensuring its correct and safe use.

COMPOSITION

Primary packaging: surgical grade paper envelope containing 1 to 2 cannulas surrounded by a plastic protective tube.

Secondary packaging: surgical grade paper envelope or individual packaging box, containing 1 radiopaque adhesive grid.

Tertiary packaging: transport box.

Code	Diameter	Length	Quantity
SBK21100P-1	21G	100mm	1
SBK21150P-1	21G	150mm	1
SBK21100P-2	21G	100mm	2
SBK21150P-2	21G	150mm	2

USE INSTRUCTIONS

- Verify that the product packaging is not damaged;
- Remove the Cannula Safe Block and the radiopaque adhesive grid from your sterile containers in a clean environment using aseptic techniques;
- Make sure the cannula and the radiopaque adhesive grid are integral;
- Position the patient according to the surgical procedure to be performed;
- Perform the asepsis at the incision site, position and attach the radiopaque adhesive grid to the patient's skin at the incision site;
- Perform an image of fluoroscopy, ultrasound or x-rays;
- Mark the incision site with a pen through the open windows of the radiopaque adhesive grid;
- Remove the radiopaque adhesive grid completely or at the discretion of the medical professional; the adhesive can be held in the patient so that the procedure is performed by opening the windows, as marked;
- Aspirate the anesthetic solution / drug to be injected into a syringe.
- Insert the cannula of the safe block into the patient, as marked on the radiopaque adhesive grid;
- Connect the stimulation cable to the standard low frequency stimulator connector for nerve location through the stimulus, and then connect the syringe containing the anesthetic solution / syringe into the injector tube with the luer lock connection for nerve block;
- After completion of the procedure, the material should be disposed of and discarded in an appropriate place, following the rules applicable to the disposal of contaminated and punctured materials adopted by the hospital.
- The product should not be reused and / or resterilized, being of SINGLE use;

- Read the instructions for use that accompany the product before opening the sterile package, ensuring its correct and safe use.
- The product should not be reused and or resterilized, being of single use;
- After the procedure, discard them in HOSPITAL WASTE.

CONTRAINDICATIONS

- Infection or hematoma in the vicinity of the site chosen for a puncture;
- Preoperative neurological changes are a relative contraindication;
- Stimulation: Do not perform for defibrillation or cardioversion;
- Changes in coagulation, for example, during anticoagulant therapy are a relative contraindication;
- Lymphangitis (eg in the case of axillary access);
- Post-lymphadenectomy status (for example, in the case of axillary access);
- acute plexus neuritis;
- Inability to obtain patient consent;
- Known hypersensitivity to the resources used.

PRECAUTIONS, RESTRICTIONS AND WARNINGS:

- The user, before using the product, should carefully read the use instructions that will aid him in the correct and safe use and will also protect against possible risks.
- Do not use the product out of the validity date, always verify if the sterilization indication is correct.
- The product should only be used by qualified medical professionals, who have the necessary knowledge and experience for the use of the product, knowledge of the handling of sterile products and, mainly, of the importance of maintaining an atmosphere totally sterile in the surgery rooms, where the care and precaution with each detail are fundamental in order to avoid contamination risks.
- This product is sterilized by Ethylene Oxide (ETO).
- The product must be transported and conditioned against all types of dirt, grime and contact with flammable, corrosive, toxic or explosive substances.
- The package should not be punctured or torn, it should be kept away from sharp or pointed instruments.
- Do not to keep or store this product close or near to products of different categories.
- Verify the integrity of the packing before using the product. In case of violation, DO NOT USE, because the warranty of the product sterility will be committed,
- It is prohibited to reprocess this product.
- This product is appropriate for only one use, sterilized by Ethylene Oxide (EO).

This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

Lot number, Sterilization and Validity Date: see the product packing

MS Registration: 80195520028

Technical Responsible: Marcos Rojas

CREA: 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

www.bramsys.com.br