
Kit Cânula G.E.N Descartável

Reg. ANVISA: 80195529007.

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ATTENTION: Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format — digital medium (www.bramsys.com.br/instrucoes-de-uso) — complies with RDC No. 751/2022. It requires the use of a PDF reader, which can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the purchased product.

To obtain the instruction for use in printed format, contact us by email at sac@bramsys.com.br. There is no additional cost for receiving the product's instruction for use in either printed or digital format.

This instruction for use is compatible with the commercial model DISPOSABLE G.E.N. CANNULA KIT.

USE INSTRUCTIONS

KIT CÁNULA G.E.N DESECHABLE

INDICATION FOR USE/PURPOSE

The Disposable G.E.N. Cannula Kit is indicated for the treatment of chronic pain through minimally invasive surgical radiofrequency procedures with controlled temperature, in order to relieve pain and improve functional capacity. The cannula features a radiopaque, non-stick, insulating coating, a thermocouple stylet, and a sonovisible active tip, allowing the procedure to be guided by ultrasound or fluoroscopy. The Disposable G.E.N. Cannula Kit is a single-use product designed to insulate and conduct the radiofrequency electrode. Its composition includes: a cannulated introducer cannula with a radiopaque, sonovisible, non-stick insulated surface, a thermocouple stylet, and an active tip. The instructions for use accompanying the product must be read before opening the sterile packaging to ensure correct and safe use.

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- **Primary packaging:** Surgical-grade paper envelope containing 1, 2, or 3 cannulas wrapped in a protective plastic tube.
- **Secondary packaging:** Surgical-grade paper envelope or individual packaging box.
- **Tertiary packaging:** Transport box.

COMMERCIAL MODELS

Code	Diameter	Length	Active Tip	Composition
KGEN1-1805-5	18G	50mm	5mm	1 Cannula
KGEN1-1810-10	18G	100mm	10mm	1 Cannula
KGEN1-1810-5	18G	100mm	5mm	1 Cannula
KGEN1-1815-10	18G	150mm	10mm	1 Cannula
KGEN1-1815-15	18G	150mm	15mm	1 Cannula
KGEN1-1815-5	18G	150mm	5mm	1 Cannula
KGEN1-2005-5	20G	50mm	5mm	1 Cannula
KGEN1-2010-10	20G	100mm	10mm	1 Cannula
KGEN1-2010-5	20G	100mm	5mm	1 Cannula
KGEN1-2015-10	20G	150mm	10mm	1 Cannula
KGEN1-2015-15	20G	150mm	15mm	1 Cannula
KGEN1-2015-5	20G	150mm	5mm	1 Cannula
KGEN1-2105-5	21G	50mm	5mm	1 Cannula
KGEN1-2110-10	21G	100mm	10mm	1 Cannula
KGEN1-2110-5	21G	100mm	5mm	1 Cannula
KGEN1-2115-10	21G	150mm	10mm	1 Cannula
KGEN1-2115-15	21G	150mm	15mm	1 Cannula
KGEN1-2115-5	21G	150mm	5mm	1 Cannula
KGEN1-2205-5	22G	50mm	5mm	1 Cannula
KGEN1-2210-10	22G	100mm	10mm	1 Cannula
KGEN1-2210-5	22G	100mm	5mm	1 Cannula
KGEN1-2215-10	22G	150mm	10mm	1 Cannula
KGEN1-2215-15	22G	150mm	15mm	1 Cannula
KGEN1-2215-5	22G	150mm	5mm	1 Cannula
KGEN2-1805-5	18G	50mm	5mm	2 Cannulas
KGEN2-1810-10	18G	100mm	10mm	2 Cannulas
KGEN2-1810-5	18G	100mm	5mm	2 Cannulas
KGEN2-1815-10	18G	150mm	10mm	2 Cannulas
KGEN2-1815-15	18G	150mm	15mm	2 Cannulas
KGEN2-1815-5	18G	150mm	5mm	2 Cannulas
KGEN2-2005-5	20G	50mm	5mm	2 Cannulas
KGEN2-2010-10	20G	100mm	10mm	2 Cannulas
KGEN2-2010-5	20G	100mm	5mm	2 Cannulas
KGEN2-2015-10	20G	150mm	10mm	2 Cannulas
KGEN2-2015-15	20G	150mm	15mm	2 Cannulas
KGEN2-2015-5	20G	150mm	5mm	2 Cannulas
KGEN2-2105-5	21G	50mm	5mm	2 Cannulas

KGEN2-2110-10	21G	100mm	10mm	2 Cannulas
KGEN2-2110-5	21G	100mm	5mm	2 Cannulas
KGEN2-2115-10	21G	150mm	10mm	2 Cannulas
KGEN2-2115-15	21G	150mm	15mm	2 Cannulas
KGEN2-2115-5	21G	150mm	5mm	2 Cannulas
KGEN2-2205-5	22G	50mm	5mm	2 Cannulas
KGEN2-2210-10	22G	100mm	10mm	2 Cannulas
KGEN2-2210-5	22G	100mm	5mm	2 Cannulas
KGEN2-2215-10	22G	150mm	10mm	2 Cannulas
KGEN2-2215-15	22G	150mm	15mm	2 Cannulas
KGEN2-2215-5	22G	150mm	5mm	2 Cannulas
KGEN3-1805-5	18G	50mm	5mm	3 Cannulas
KGEN3-1810-10	18G	100mm	10mm	3 Cannulas
KGEN3-1810-5	18G	100mm	5mm	3 Cannulas
KGEN3-1815-10	18G	150mm	10mm	3 Cannulas
KGEN3-1815-15	18G	150mm	15mm	3 Cannulas
KGEN3-1815-5	18G	150mm	5mm	3 Cannulas
KGEN3-2005-5	20G	50mm	5mm	3 Cannulas
KGEN3-2010-10	20G	100mm	10mm	3 Cannulas
KGEN3-2010-5	20G	100mm	5mm	3 Cannulas
KGEN3-2015-10	20G	150mm	10mm	3 Cannulas
KGEN3-2015-15	20G	150mm	15mm	3 Cannulas
KGEN3-2015-5	20G	150mm	5mm	3 Cannulas
KGEN3-2105-5	21G	50mm	5mm	3 Cannulas
KGEN3-2110-10	21G	100mm	10mm	3 Cannulas
KGEN3-2110-5	21G	100mm	5mm	3 Cannulas
KGEN3-2115-10	21G	150mm	10mm	3 Cannulas
KGEN3-2115-15	21G	150mm	15mm	3 Cannulas
KGEN3-2115-5	21G	150mm	5mm	3 Cannulas
KGEN3-2205-5	22G	50mm	5mm	3 Cannulas
KGEN3-2210-10	22G	100mm	10mm	3 Cannulas
KGEN3-2210-5	22G	100mm	5mm	3 Cannulas
KGEN3-2215-10	22G	150mm	10mm	3 Cannulas
KGEN3-2215-15	22G	150mm	15mm	3 Cannulas
KGEN3-2215-5	22G	150mm	5mm	3 Cannulas

PRODUCT INSTRUCTIONS FOR USE

- 1 - Ensure that the packaging is not damaged.
- 2 - Remove the G.E.N. Cannula from the sterile packaging in a clean environment using aseptic techniques.
- 3 - Verify that the G.E.N. Cannula is intact.
- 4 - Position the patient according to the surgical procedure to be performed.
- 5 - Prepare the skin puncture area using aseptic and antiseptic techniques.
- 6 - Perform percutaneous puncture on the patient using the G.E.N. Cannula, keeping the thermocouple stylet in place.
- 7 - Note: If necessary, use fluoroscopy, ultrasonography, or X-ray to assist with puncture and target monitoring.
- 8 - Connect the cable to the Radiofrequency Generator Unit to perform stimulation, lesioning, and/or neuromodulation procedures with controlled temperature.
- 9 - Remove the G.E.N. Cannula and thermocouple stylet from the patient.
- 10 - After completing the procedure, the material must be rendered unusable and disposed of in an appropriate location, following the regulations applicable to the disposal of sharps and contaminated materials adopted by the hospital.

CONTRAINDICATIONS

The Disposable G.E.N. Cannula Kit is contraindicated in the following situations:

- Infection or hematoma near the selected puncture site;
- Neurological alterations known in the preoperative period (relative contraindication);
- Stimulation: do not perform during defibrillation or cardioversion;
- Coagulation disorders, for example during anticoagulant therapy (relative contraindication);
- Lymphangitis (e.g., in cases of axillary access);
- Post-lymphadenectomy status (e.g., in cases of axillary access);
- Acute plexus neuritis;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used.

COMPATIBILITY WITH OTHER MEDICAL DEVICES

The Disposable G.E.N. Cannula Kit is compatible with and intended for use with the Radiofrequency Generator Unit according to the following technical specifications:

- Power supply voltage range: 90 – 260 VAC
- Power supply frequency: 47 – 63 Hz
- Maximum electrical output power: 50 Wrms
- Wave frequency: 100 – 500 kHz
- Output voltage range: 0 – 150 V
- Output current range: 0 – 999 mA
- Output pulse frequency: 1 – 200 Hz
- Output pulse duration: 0.1 – 50 ms
- Impedance measurement: 0 – 5000 ohms
- Temperature measurement: 25 – 120°C

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will assist in proper and safe operation and help prevent possible risks.
- Do not use the product past its expiration date; always check that the sterilization indicator is correct.
- The product must be used only by qualified medical professionals who have the necessary knowledge and experience for its proper use, are familiar with the handling of sterile products, and understand the importance of maintaining a completely sterile environment in operating rooms, where attention to every detail is essential to avoid risks and contamination.
- The product must not be reused and/or reprocessed; it is for single use only and must be destroyed and disposed of in hospital waste.
- The product must always be transported and stored protected from all types of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.
- Do not store or keep this product near products from different categories.
- Check the integrity of the packaging before using the product; in case of any damage or violation, DO NOT USE, as sterility can no longer be guaranteed.
- Reprocessing of this product is strictly prohibited.
- Single-use product, sterilized by Ethylene Oxide (ETO).

TRANSPORT CONDITIONS

The Disposable G.E.N. Cannula Kit must always be transported and stored protected from all types of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances. The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.

IMPORTANT: The product must be stored and transported in accordance with Good Practices for the Transport of Sterilized Products. The product must be transported in its original packaging, properly stored. Do not place any weight or other volumes on top of it, as this may damage its technical characteristics and compromise its integrity. Do not place it near or together with products from different categories, as this may damage its technical characteristics and compromise its integrity or quality.

Handle with care during transportation, as improper handling may compromise the integrity and quality of the product provided by the manufacturer.

OPERATING CONDITIONS

Before using the product, the user must carefully read the instructions for use, which will assist in its correct and safe operation and help prevent possible risks.

Do not use the product past its expiration date; always verify that the sterilization indicator is correct.

The product must only be used by qualified medical professionals who possess the knowledge and experience necessary for proper use, sterile product handling, and—most importantly—the understanding of maintaining a completely sterile environment in operating rooms, where attention to every detail is essential to prevent contamination risks. Check the integrity of the packaging before use; in case of any damage or violation, DO NOT USE, as the sterility of the product will be compromised. This product is sterilized with Ethylene Oxide (ETO).

It is completely safe when used by qualified professionals and in accordance with the usage specifications described in this document.

Temperature: 10°C to 40°C

Humidity: 20% to 90%

ANVISA Registration No.: 80195529007

WARRANTY CERTIFICATE

The company **BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, duly authorized by the National Health Surveillance Agency, guarantees for **24 MONTHS** from the date of sterilization the quality and sterile condition of this product, except in cases of misuse or failure to follow the instructions contained in this document.

CUSTOMER SERVICE (SAC)

Tel.: (+55 11) 5013-0010

E-mail: sac@bramsys.com.br

Website: www.bramsys.com.br