Reg. ANVISA: 80195529013.

[illegible]

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Read all instructions carefully before use. Follow all warnings and precautions mentioned in these instructions. Failure to do so may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) complies with RDC No. 751/2022. It requires the use of a PDF reader, which can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from that of the purchased product.

To obtain the instruction for use in printed format, contact us via email at sac@bramsys.com.br. There is no additional cost for sending the product's instruction for use in either printed or digital format.

This instruction for use is compatible with all commercial models of the Osteoplasty Cannula Kit family.

USE INSTRUCTIONS

OSTEOPLASTY CANNULA KIT

INDICATION FOR USE

The Osteoplasty Cannula Kit is indicated for osteoplasty procedures in vertebral compression fractures, aiming for stabilization, anatomical repositioning of the vertebral body, and reduction of deformities resulting from bone collapse.

PURPOSE

The Osteoplasty Cannula Kit is intended to provide minimally invasive percutaneous access to the fractured or compromised vertebra, via the use of working cannulas. Following access, an expansion device is used to restore the height of the vertebral body and create a cavity, which must be filled with bone cement for stabilization and structural reinforcement.

COMMERCIAL MODELS AND CONTENT

Code	Working Cannula	Balloon	Composition
OTP11G-10	11G - 120mm	10mm	Unilateral
OTP11G-15	11G - 120mm	15mm	Unilateral
OTP11G-20	11G - 120mm	20mm	Unilateral

The commercial models of the Osteoplasty Cannula Kit product are composed of the following items:

Items	Technical Specifications
(2) Working Cannula	Diameter: 11G Length: 120mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(4) Winged Syringe	Administration Capacity: 2.5 mL Chemical Composition: Polypropylene, Polystyrene, and Rubber
(1) Drill	Diameter: 13G Length: 187 mm Chemical Composition: 440 Stainless Steel and ABS Copolymer
(2) Inflator Device	Operating Range: -0.4 ATM to 30 ATM (1 ATM = 1 BAR = 14.7 PSI) Accuracy: $\pm 3.0\%$ of full scale Administration Capacity: 20 mL Administered Fluid: 0.45 mL (± 0.07 mL) per 360° clockwise rotation of the syringe plunger handle Chemical Composition: Polycarbonate, Polypropylene, and Silicone
(2) 3-Way Luer Lock Adapter	Luer Lock Connection Chemical Composition: Polycarbonate
(2) Balloon Catheter	Length: 10 mm (OTP11G-10), 15 mm (OTP11G-15), and 20 mm (OTP11G-20). Capacity: <400 psi Maximum volume: 10 mm model (3 mL); 15 mm model (4 mL); 20 mm model (6 mL). Chemical composition: Stainless steel 304, polycarbonate, PTFE, and silicone.
(6) Bone Cement Delivery Cannula	Diameter: 13G Length: 187 mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(1) Bone Cement Mixer	Administration Capacity: 40 g Chemical Composition: Polypropylene, Polystyrene, and 304 Stainless Steel

The commercial models are provided in the following presentation:

- Primary packaging: PET (polyethylene terephthalate) blister sealed with Tyvek paper, permeable to ethylene oxide for sterilization, containing all the items that make up the kit, as presented above, along with the 5 traceability labels.
- Secondary packaging: Individual packaging box, made of duplex paper, containing the primary packaging.
- Tertiary packaging: Transport box.

Each package contains only 1 (one) commercial model of the Osteoplasty Cannula Kit.

PRODUCT USE PROCEDURE

- 1- Define the site for access puncture to the vertebra.
- 2- Insert the Working Cannula up to the desired target.
- 3- Remove the mandrel from the Working Cannula.
- 4- Insert the Drill into the Working Cannula.

5- Perform circular motions to perforate the vertebra, in order to prepare it for inserting the Balloon Catheter.
6- Remove the Drill from the Working Cannula.
7- Collect 10ml of Contrast medium with the aid of a syringe.
8- Connect the syringe to the 3-way adapter.
9- Position the lock of the 3-way adapter to inject the contrast.
10- Connect the 3-way adapter to the Insuflator Device.
11- Unlock the plunger of the Insuflator Device.
12- Inject the contrast into the Insuflator Device.
13- Decouple the syringe from the 3-way adapter.
14- Position the lock of the 3-way Adapter for Balloon Catheter use.
15- Remove the air between the plunger of the Insuflator Device up to the outlet of the 3-way adapter.
16- Insert the Balloon Catheter into the Working Cannula.
17- Perform positioning adjustment if necessary until the radiopaque markers of the Balloon Catheter are outside the Working Cannula.
18- Connect the 3-way Adapter to the Balloon Catheter access.
19- Turn the plunger of the Insuflator Device clockwise to inject the contrast into the Balloon Catheter and inflate it, in order to prepare the cavity in the vertebral body.
20- Prepare the bone cement with the aid of the bone cement mixer.
21- Couple the syringe with wings to the bone cement mixer.
22- Collect the bone cement into the syringe with wings.
23- Connect the syringe with wings to the Cannula for bone cement application.
24- Inject the bone cement into the Cannula for bone cement application until it is completely filled.
25- Turn the plunger of the Insuflator Device counter-clockwise until all the contrast is removed from the Balloon Catheter.
26- Unlock the plunger of the Insuflator Device and pull it to create negative pressure in the Balloon Catheter and ensure it is completely empty.
27- Remove the Balloon Catheter from the expansive working Cannula.
28- Insert the Cannula for bone cement application into the expansive working Cannula.
29- Insert the mandrel into the Cannula for bone cement application to transfer the cement to the vertebra.
30- Remove the Cannula for bone cement application from the Working Cannula.
31- Remove the Working Cannula.
32- Dispose of all material used in the procedure in an appropriate location following the regulations applicable to the disposal of sharp/cutting and contaminated material adopted by the hospital.

CONTRAINDICATIONS

- Infection or hematoma in the proximity of the chosen puncture site.
- Pregnancy.
- Neurological conditions known preoperatively are a relative contraindication.
- Coagulation disorders, for example, during anticoagulant therapy, are a relative contraindication.
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to the contrast medium, antibiotics, and/or local anesthetics used in the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedure performance.

CONTRAINDICATIONS

- Infection or hematoma near the selected puncture site.
- Known neurological alterations in the preoperative period are a relative contraindication.
- Stimulation: Do not perform during defibrillation or cardioversion.
- Coagulation disorders, for example, during anticoagulant therapy, are a relative contraindication.
- Lymphangitis, in the case of axillary access.
- Post-lymphadenectomy status, for example, in the case of axillary access.
- Acute plexus neuritis.
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- The user, before using the product, must carefully read the instructions for use, which will assist in the correct and safe use, as well as protect against possible risks.
- Do not use the product after its expiration date; always verify that the sterilization indication is correct.
- The product must only be used by qualified medical professionals who possess the necessary knowledge and experience for the proper use of the product, knowledge of handling sterile products, and especially the

importance of maintaining a completely sterile environment in operating rooms, where care and precaution with every detail are fundamental to avoid risks of contamination.

- Check the integrity of the packaging before using the product; in case of tampering, DO NOT USE, as the guarantee of its sterility will be compromised.
- This product is for single use, sterilized by Ethylene Oxide (ETO).
- Reprocessing this product is prohibited.
- No modification is permitted for this product.
- Disposal of the product after use must follow the regulations applicable to the disposal of sharp/cutting and contaminated material, according to the hospital's Healthcare Waste Management Plan, in agreement with the guidance of the Hospital Infection Control Committee and as foreseen in RDC 222/2018.

STORAGE AND TRANSPORT CONDITIONS

- Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be pierced or torn, and it must be kept away from sharp or pointed instruments.
- Do not place under weight or volumes, as this may cause damage to the product's technical characteristics.
- Do not place near or alongside different categories, as they may cause damage to its technical characteristics.
- Must be packaged and transported following good practices for the transport of sterilized products.
- This product is sterilized by ethylene oxide (ETO) and its packaging must not be damaged before use.
- The product must be transported in its original packaging.
- Storage Temperature: Between 15°C to 30°C.
- Storage Humidity: 30% to 80%.

ANVISA Registration No.: 80195529013

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency (ANVISA) under No. UY44WX94W79M, guarantees the quality and sterile condition of this product for **24 months** from the sterilization date, except in cases of misuse or failure to follow the instructions provided in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br