
Kit Cânula para UnikBlock

Reg. ANVISA: 80195529008.

BRAMSYS

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AV. DR. ALBERTO SAN JUAN, 350

18530-000 – TIETÊ – SP

CNPJ: 04.729.995/0001-87

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KIT CÂNULA UNIKBLOCK

FINALIDADE

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UNIKBLOCK CANNULA KIT

PURPOSE

The cannula is introduced into the patient and its location can be determined anatomically with or without the aid of fluoroscopy, ultrasound (through its echogenic tip) or x-ray. For a safer location, the aid of a radiopaque adhesive grid is used, where it is positioned and adhered to the patient's skin in the area of the surgical incision in a preoperative and aseptic moment. The radiopaque patch is shown on the fluoroscopy or x-ray image, guiding the site to access the required surgical point, acting as a skin reference point in the radiological image obtained, where the reference points are used to identify the incision site. surgical. The UnikBlock Cannula Kit allows the use of percutaneous electrical stimulation together to confirm the location of the selected target, allowing the execution of the desired procedure and the application of temperature-controlled neuromodulation.

INDICATIONS OF USE

The UnikBlock Cannula Kit is indicated for procedures of Plexus Block, Peripheral Nerves and Nerve Blocks of the Musculoskeletal System, guided or not by image method, with or without the use of the neuromodulation resource with controlled temperature and the aid of electrical percutaneous stimulation to target location. These procedures can be performed together, or separately, at the discretion of the physician.

COMPOSITION

1, 2, 3, or 4 cannulas (see number after UNK in commercial code), 1 surgical planning scale, 1 surgical marker, 1 radiopaque adhesive grid.

USE INSTRUCTIONS

Procedure not using the radiopaque adhesive grid:

1. Make sure the packaging is not damaged.
 2. Remove cannula, surgical planning scale, surgical marker and radiopaque adhesive grid from sterile packaging in a clean environment using aseptic techniques.
 3. Make sure the cannula, surgical planning scale, surgical marker and radiopaque adhesive grid are intact.
 4. Position the patient according to the surgical procedure to be performed.
 5. Prepare skin puncture area using aseptic and antiseptic techniques.
 6. Perform the percutaneous puncture with the cannula in the patient.
- NOTE: If necessary, use fluoroscopy, ultrasound, or x-ray to aid puncture and target monitoring.
7. Connect the cable to the neurostimulation equipment or to the Radiofrequency Generating Unit to perform stimulation, injury and/or neuromodulation procedures with motor, sensory and/or controlled temperature stimulation.
 8. Remove the cannula from the patient

9. After the end of the procedure, the material must be de-characterized and disposed of in an appropriate place, following the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.

Procedure using the radiopaque adhesive grid:

1. Make sure the packaging is not damaged
 2. Remove cannula, surgical planning scale, surgical marker and radiopaque adhesive grid from sterile packaging in a clean environment using aseptic techniques.
 3. Make sure the cannula, surgical planning scale, surgical marker and radiopaque adhesive grid are intact.
 4. Position the patient according to the surgical procedure to be performed.
 5. Prepare skin puncture area using aseptic and antiseptic techniques.
 6. Position and secure the radiopaque adhesive grid on the patient's skin at the incision site.
 7. Perform fluoroscopy or x-ray imaging.
 8. Mark the incision site with the surgical planning scale and surgical Marker through the open windows of the radiopaque adhesive grid.
- NOTE: Remove the radiopaque adhesive grid in its entirety or at the discretion of the medical professional, the adhesive can be kept on the patient, so that the procedure is performed by opening the windows, as marked.
9. Perform the percutaneous puncture with the cannula in the patient.
- NOTE: If necessary, use fluoroscopy, ultrasound, or x-ray to aid puncture and target monitoring.
10. Connect the cable to the neurostimulation equipment or to the Radiofrequency Generating Unit to perform stimulation, injury and/or neuromodulation procedures with motor, sensory and/or controlled temperature stimulation.
 11. Remove the cannula from the patient.
 12. After the end of the procedure, the material must be de-characterized and disposed of in an appropriate place, following the regulations applicable to the disposal of sharps and contaminated material adopted by the hospital.

UnikBlock Cannula Kit was duly produced for use and compatibility with the p.Connect Adapter Cable, where it has an individual registration with ANVISA (ANVISA Registration 80195529006).

CONTRAINDICATIONS

- Infection or hematoma near the site chosen for the puncture;
- Neurological alterations known in the preoperative period are a relative contraindication;
- Stimulation: do not perform during defibrillation or cardioversion;
- Changes in coagulation, for example anticoagulant therapy are a relative contraindication;
- Lymphangitis (for example, in case of axillary access);
- Post-Lymphadenectomy status (for example, in case of axillary access);
- Acute plexus neuritis;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used.

STORAGE CONDITIONS

The Disposable UnikBlock Cannula Kit must always be stored protected against all types of moisture, dirt and contact with flammable, corrosive, toxic or explosive substances. The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.

Do not place it under weight or volumes, as they may damage its technical characteristics and compromise its integrity.

Do not store the product where it can be exposed to direct sunlight, high temperatures, humidity or dust as this will compromise its integrity/quality.

Do not place them next to or together with different categories, as they may damage their technical characteristics and compromise their integrity / quality.

Storage temperature: 5°C to 60°C.

Storage humidity: 10% to 80%.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

- The user, before using the product, must carefully read the instructions for use that will help correct and safe use, as well as protecting against possible risks.
- Do not use the product after its expiration date, always check if the sterilization indication is correct.
- The product must only be used by qualified medical professionals, who have the knowledge and experience necessary for the proper use of the product, knowledge of handling sterile products and, mainly, of the importance of maintaining a totally sterile environment in the operating rooms. , where care and precaution with every detail

are essential so that there is no risk of contamination.

- Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
 - This product is sterilized by Ethylene Oxide (ETO).
 - The product must be transported and stored, always protected against all types of moisture, dirt and contact with flammable, corrosive, toxic or explosive substances.
 - The packaging must not be punctured or torn, it must be kept away from sharp or pointed instruments.
 - Do not store or store this product next to products of different categories.
 - Check the integrity of the packaging before using the product, in case of violation, DO NOT USE, as the guarantee of its sterility will be compromised.
 - It is prohibited to reprocess this product.
 - The product must not be reused and/or re-sterilized, it is for SINGLE use, and must be destroyed and disposed of in HOSPITAL WASTE.
- This product is totally safe if used by qualified professional and respecting the usage specifications described in this document.

Lot number, Sterilization and Validity Date: see the product packing.

MS Registration: 80195529008

Technical Responsible: Marcos Rojas | **CREA:** 5061214879

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

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