
Kit Cânula b.Pack

Reg. ANVISA: 80195529015.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) – is in accordance with RDC No. 751/2022. A PDF reader is required, and the program can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the product purchased.

To obtain the instruction for use in printed format, contact us via email at sac@bramsys.com.br. There is no additional cost for sending the product instruction for use in printed or digital format.

This instruction for use is compatible with all commercial models of the b.Pack Cannula Kit family.

USE INSTRUCTIONS

B.PACK CANNULA KIT

INDICATION FOR USE / PURPOSE

The b.Pack Cannula Kit is indicated for anesthetic or neurolytic blocks of peripheral nerves, nerve plexuses, sympathetic ganglia, and joint blocks, as well as for diagnostic or therapeutic infiltrations and punctures (articular or extra-articular), with or without imaging guidance. It provides safety and precision during procedures and can be used with local anesthetics, neurolytic agents, and/or botulinum toxin, as well as with electrical stimulation (sensory and motor) to locate the target. The choice of technique and agents is the responsibility of the physician.

PRINCIPLE OF OPERATION / MECHANISM OF ACTION

The cannula is inserted into the patient, and its location can be determined anatomically, with or without the aid of fluoroscopy, ultrasound (via its echogenic tip), or X-ray. The All One Block Cannula Kit allows the combined use of percutaneous electrical stimulation to confirm the location of the selected target, enabling the execution of the desired procedure and the application of controlled temperature neuromodulation.

COMMERCIAL MODELS AND CONTENT

Code	Ø	Longitud	Active Tip
BPK-18050-5	18G	50mm	5
BPK-18070-5	18G	70mm	5
BPK-18070-10	18G	70mm	10
BPK-18100-5	18G	100mm	5
BPK-18100-10	18G	100mm	10
BPK-18150-5	18G	150mm	5
BPK-18150-10	18G	150mm	10
BPK-20050-5	20G	50mm	5
BPK-20070-5	20G	70mm	5
BPK-20100-5	20G	100mm	5
BPK-20100-10	20G	100mm	10
BPK-20150-5	20G	100mm	5
BPK-20150-10	20G	150mm	10
BPK-21050-5	21G	50mm	5
BPK-21070-10	21G	70mm	10
BPK-21100-5	21G	100mm	5
BPK-21100-10	21G	100mm	10
BPK-21150-5	21G	150mm	5
BPK-21150-10	21G	150mm	10

The commercial models of the b.Pack Cannula Kit consist of the following items:

Items	Technical Specifications
(1) Cannula	Diameter: 18G, 20G, and 21G Length: 50 mm, 70 mm, 100 mm, and 150 mm Active tip: 5 mm and 10 mm Chemical composition: Stainless steel 304L/316L, polyester, polyacetal, low-density polyethylene
(1) Surgical Planning Scale	Chemical composition: PVC
(1) Surgical Marker	Chemical composition: Polypropylene, polyester, Gentian Violet 3%, and isopropanol 30%

The commercial models are supplied in the following packaging presentation:

- Primary packaging: Surgical-grade envelope containing all items that make up the kit, as described above, along with five traceability labels.
- Secondary packaging: Individual packaging box made of duplex paper containing the primary packaging.
- Tertiary packaging: Transport box.

Each package contains only one (1) commercial model of the b.Pack Cannula Kit.

PRODUCT INSTRUCTIONS FOR USE

- 1- Ensure that the packaging is not damaged.

2- Remove the items from the package in a clean environment using aseptic techniques.

3- Position the patient according to the surgical procedure to be performed.

Note: If necessary, use the Surgical Planning Scale and Surgical Marker to assist in planning the puncture and locating the desired target.

4- Prepare the puncture site using aseptic and antiseptic techniques.

5- Perform the percutaneous puncture of the cannula.

Note: If necessary, use fluoroscopy, ultrasonography, or X-ray to assist in targeting the puncture site.

6- Connect the stimulation cable to the stimulator device to assist in nerve localization through sensory and motor stimulation.

7- Perform the procedure as intended.

8- After completing the procedure, all materials must be rendered unusable and disposed of in an appropriate location, following the applicable regulations for the disposal of sharp and contaminated materials adopted by the hospital.

CONTRAINDICATIONS

- Infection or hematoma near the procedure site.
- Pregnancy.
- Preoperative neurological alterations are considered a relative contraindication.
- Coagulation disorders, for example during anticoagulant therapy, are considered relative contraindications.
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to contrast agents, antibiotics, and/or local anesthetics used during the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedural performance.
- The use of general anesthesia is not recommended. During the procedure, a sensory response from the patient may be required; therefore, it is advised that the treatment be performed under local anesthesia.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will assist in its correct and safe handling, as well as protect against possible risks.
- Do not use the product past its expiration date. Always verify that the sterilization indicator is correct.
- The product must be used only by qualified medical professionals who possess the necessary knowledge and experience for its proper use, including the handling of sterile products and, most importantly, the importance of maintaining a completely sterile environment in operating rooms, where attention and care with every detail are essential to prevent contamination risks.
- Check the integrity of the packaging before use. In case of any damage or violation, DO NOT USE the product, as its sterility will be compromised.
- This product is for single use only and is sterilized with Ethylene Oxide (ETO).
- Reprocessing of this product is strictly prohibited. It must be destroyed and disposed of as hospital waste.
- No modification to this product is allowed.
- After use, disposal of the product must follow applicable regulations for sharp and contaminated materials, in accordance with the hospital's Health Service Waste Management Plan, the guidance of the Hospital Infection Control Committee, and the provisions of RDC 222/2018.

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

- Be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- Be kept in packaging that is not punctured or torn, and away from sharp or pointed instruments.
- Not be placed under heavy loads or volumes, as this may damage the product's technical characteristics.
- Not be stored or transported together with different product categories, as this may compromise its technical characteristics.
- Be stored and transported following good practices for the handling of sterilized products.
- This product is sterilized with Ethylene Oxide (ETO) and must not have its packaging damaged before use.
- Be transported in its original packaging.
- Storage temperature: Between 15°C and 30°C.
- Storage humidity: Between 30% and 80%.

ANVISA Registration No.: 80195529015

WARRANTY CERTIFICATE

The company **BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, duly authorized by the National Health Surveillance Agency (ANVISA) under registration No. UY44WX94W79M, guarantees for a period of **24 months** from the date of sterilization the quality and sterile condition of this product, except in cases of misuse or failure to comply with the instructions

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br