
Kit Cântula Cifoplastia d.Kypho

Req. ANVISA: 80195520031.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Read all instructions carefully before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) is in accordance with RDC No. 751/2022. It requires the use of a PDF reader, which can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the acquired product.

To obtain the instruction for use in printed format, contact us via email at sac@bramsys.com.br. There is no additional cost for sending the product's instruction for use in printed or digital format.

This instruction for use is compatible with the commercial model Kit Cannula Kyphoplasty d.Kypho.

USE INSTRUCTIONS

KIT CANNULA KYPHOPLASTY D.KYPHO

INDICATION FOR USE / PURPOSE

The d.Kypho Kyphoplasty Cannula Kit is intended to provide minimally invasive percutaneous access to a fractured or damaged vertebra through the working cannulas. Once the vertebra has been accessed, the balloon is used to restore the height of the vertebral body, creating a cavity that should be filled through the application of bone cement. The kit also allows for the collection of samples for biopsy procedures. The d.Kypho Kyphoplasty Cannula Kit is indicated for minimally invasive kyphoplasty procedures in the treatment of vertebral fractures or spinal damage.

COMMERCIAL PRESENTATION

Code	Composition
DKU11G-10	(1) Working Expansive Cannula (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (1) Inflator Device (1) 3-Way Luer Lock Adapter (1) Balloon Catheter (3) Cannulas for Bone Cement Application (1) Bone Cement Mixer
DKU11G-15	(1) Working Expansive Cannula (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (1) Inflation Device (1) 3-Way Luer Lock Adapter (1) Balloon Catheter (3) Bone Cement Application Cannulas (1) Bone Cement Mixer
DKU11G-20	(1) Working Expansive Cannula (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (1) Inflation Device (1) 3-Way Luer Lock Adapter (1) Balloon Catheter (3) Bone Cement Application Cannulas (1) Bone Cement Mixer
DKB11G-10	(2) Working Expansive Cannulas (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (2) Inflation Devices (2) 3-Way Luer Lock Adapters (2) Balloon Catheters (6) Bone Cement Application Cannulas (1) Bone Cement Mixer
DKB11G-15	(2) Working Expansive Cannulas (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (2) Inflation Devices (2) 3-Way Luer Lock Adapters (2) Balloon Catheters (6) Bone Cement Application Cannulas (1) Bone Cement Mixer
DKB11G-20	(2) Working Expansive Cannulas (1) Biopsy Cannula (1) Guide Wire (4) Winged Syringes (1) Drill (2) Inflation Devices (2) 3-Way Luer Lock Adapters (2) Balloon Catheters (6) Bone Cement Application Cannulas (1) Bone Cement Mixer

TECHNICAL SPECIFICATIONS

- Expansive Working Cannula
 - Diameter: 11G
 - Length: 120 mm

- Biopsy Cannula
 - Diameter: 15G
 - Length: 185 mm
- Guide Wire
 - Diameter: 16G
 - Length: 240 mm
- Winged Syringe
 - Delivery capacity: 2.5 mL
- Drill
 - Diameter: 13G
 - Length: 187 mm
- Inflation Device
 - Operating range: -0.4 ATM to 30 ATM (1 ATM = 1 BAR = 14.7 PSI)
 - Accuracy: $\pm 3.0\%$ of full scale
 - Delivery capacity: 20 mL
 - Delivered fluid: 0.45 mL (± 0.07 mL delivered) per each 360° clockwise rotation of the syringe plunger knob
 - Fluid temperature: 10°C to 40°C
 - Humidity and operating range: 20% to 90%
- 3-Way Adapter
 - Luer lock connection
- Balloon Catheter
 - Length: 10, 15, and 20 mm
 - Capacity: <400 psi
 - Maximum volume: 10 mm model (3 mL); 15 mm model (4 mL); 20 mm model (6 mL)
- Bone Cement Mixer
 - Delivery capacity: 40 g
 - Cannula for bone cement application
 - Diameter: 13G
 - Length: 187 mm

PRODUCT INSTRUCTIONS FOR USE

1. Define the puncture site for vertebral access;
2. Insert the working expandable cannula until reaching the desired target;
3. Remove the stylet from the working expandable cannula;
4. Insert the drill into the working expandable cannula;
5. Perform circular movements to perforate the vertebra in order to prepare it for balloon catheter insertion;
6. Remove the drill from the working expandable cannula;
7. Insert the biopsy cannula into the working expandable cannula;
8. Perform circular movements to collect the samples;
9. Remove the biopsy cannula from the working expandable cannula;
10. Remove the stylet from the biopsy cannula;
11. Insert the guide wire into the biopsy cannula to remove the collected material and deposit it in an appropriate container;
12. Collect 10 mL of contrast medium using a syringe;
13. Connect the syringe to the 3-way adapter;
14. Position the 3-way adapter lock to inject the contrast medium;
15. Connect the 3-way adapter to the inflator device;
16. Unlock the plunger of the inflator device;
17. Inject the contrast medium into the inflator device;
18. Disconnect the syringe from the 3-way adapter;
19. Position the 3-way adapter lock for balloon catheter use;
20. Remove the air between the inflator device plunger and the outlet of the 3-way adapter;
21. Lock the plunger of the inflator device;
22. Insert the balloon catheter into the working expandable cannula until the positioning mark is fully inserted.
23. Adjust the positioning if necessary until the radiopaque markers of the balloon catheter are outside the working expandable cannula;
24. Remove the stylet from the balloon catheter;
25. Connect the 3-way adapter to the access port of the balloon catheter;
26. Turn the plunger knob of the inflator device clockwise to inject contrast medium into the balloon catheter and inflate it, in order to prepare the cavity in the vertebral body;
27. Prepare the bone cement using the bone cement mixer;
28. Attach the winged syringe to the bone cement mixer;
29. Collect the bone cement in the winged syringe;
30. Detach the winged syringe from the bone cement mixer;
31. Connect the winged syringe to the cannula for bone cement application;
32. Inject the bone cement into the cannula for bone cement application until it is completely filled;
33. Disconnect the winged syringe from the cannula for bone cement application;
34. Turn the plunger knob of the inflator device counterclockwise until all contrast medium is removed from the balloon catheter;
34. Unlock the plunger of the inflator device and pull it back to create

negative pressure in the balloon catheter, ensuring it is completely emptied;
36. Lock the plunger of the inflator device;
37. Remove the balloon catheter from the working expandable cannula;
38. Insert the cannula for bone cement application into the working expandable cannula;
39. Insert the stylet into the cannula for bone cement application to transfer the cement into the vertebra;
40. Remove the cannula for bone cement application from the working expandable cannula;
41. Remove the working expandable cannula;
42. Cover the incision with sterile gauze and suture if necessary.
43. Dispose of all materials used in the procedure in an appropriate location, following the regulations applicable to the disposal of sharp and contaminated materials as adopted by the hospital.

CONTRAINDICATIONS

- Infection or hematoma near the chosen puncture site.
- Pregnancy.
- Neurological alterations identified during the preoperative assessment (considered a relative contraindication).
- Coagulation disorders, for example, during anticoagulant therapy (considered a relative contraindication).
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to contrast agents, antibiotics, and/or local anesthetics used in the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedural performance.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will assist in its correct and safe use, as well as protect against possible risks.
 - Do not use the product after its expiration date. Always check that the sterilization indicator is correct.
 - The product must only be used by qualified medical professionals who possess the necessary knowledge and experience for its proper use, including handling sterile products and, most importantly, maintaining a completely sterile environment in operating rooms, where care and attention to every detail are essential to prevent contamination risks.
 - Check the integrity of the packaging before using the product. In case of any damage, DO NOT USE it, as the guarantee of its sterility will be compromised.
 - This product is for single use only and is sterilized by Ethylene Oxide (ETO).
 - Reprocessing of this product is strictly prohibited.
 - The product must be transported and stored protected from excessive humidity, high temperatures, direct sunlight, dirt, and contact with flammable, corrosive, toxic, or explosive substances, and kept away from sharp or pointed instruments.
 - Do not store or place near different product categories, as this may cause damage to its technical characteristics and compromise its integrity/quality.
 - No modification to this product is permitted.
 - Product disposal after use must comply with the hospital's applicable regulations for the disposal of sharp and contaminated materials.
- Storage and transport temperature: 5°C to 60°C
Storage and transport humidity: 10% to 80%

ANVISA Registration No. 80195520031

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency (ANVISA) under No. UY44WX94W79M, guarantees the quality and sterile condition of this product for **24 months** from the sterilization date, except in cases of misuse or failure to follow the instructions provided in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br