





Kit Cânula easyATM

Reg. ANVISA: 80195529014.

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Sitio web: www.bramsys.com.br

Revisión: 0 – 10/10/2025

ATTENTION: Read all instructions carefully before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instructions for Use in a non-printed digital format (www.bramsys.com.br/instrucoes-de-uso) is in accordance with RDC No. 751/2022. A PDF reader is required; to obtain the program, download it for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the instructions corresponding to that date or earlier. Do not use instructions with a revision/date (Revision: XX - day/month/year) different from the purchased product.

To obtain the Instructions for Use in printed format, contact us at sac@bramsys.com.br. There is no additional cost for sending the product instructions in printed or digital format.

USE INSTRUCTIONS

EASYATM CANNULA KIT

INDICATION FOR USE / PURPOSE

The easyATM Cannula Kit is a set of disposable surgical instruments intended for minimally invasive procedures in the temporomandibular joint (TMJ). Its use is aimed at the diagnosis and treatment of various intra-articular pathologies, such as: internal TMJ disorders; inflammatory and degenerative pathologies; intra-articular fibrosis and adhesions; mechanical joint locking; among others.

COMMERCIAL MODEL AND CONTENTS

EATM-1252 – easyATM Arthroscopy Cannula Kit

The commercial model of the easyATM Cannula Kit consists of the following items

Items	Technical Specifications
(1) Trocar Cannula	Shaft diameter: 2.6 mm Shaft length: 52 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Trocar Cap	Length: 23 mm Chemical composition: Polyacetal
(1) Hook Needle Tip	Shaft diameter: 1.9 mm Shaft length: 80 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Curette Needle Tip	Shaft diameter: 1.9 mm Shaft length: 80 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Banana Knife	Shaft diameter: 1.9 mm Shaft length: 82 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Meniscotome Needle Tip	Shaft diameter: 1.9 mm Shaft length: 80 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Blunt Obturator Tip	Shaft diameter: 1.9 mm Shaft length: 75 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Sharp Obturator Tip	Shaft diameter: 1.9 mm Shaft length: 75 mm Chemical composition: Stainless steel 304/316L and polyacetal
(1) Surgical Marker	Chemical composition: Polypropylene, polyester, 3% Gentian Violet, and 30% isopropanol
(1) Surgical Planning Scale	Chemical composition: PVC

The commercial model is provided in the following presentation:

Primary packaging: PET blister (polyethylene terephthalate) sealed with Tyvek paper, permeable to ethylene oxide for sterilization, containing all the items included in the kit as listed above, along with 5 traceability labels.

Secondary packaging: Individual packaging box, made of duplex paper, containing the primary packaging.

Tertiary packaging: Transport box.

Each primary package contains only 1 (one) commercial model of the easyATM Cannula Kit.

PRODUCT INSTRUCTIONS FOR USE

- 1 - Prepare the patient before the procedure according to standard aseptic techniques.
- 2 - Following sterile product handling procedures, open the packaging and remove the product.
- 3 - If necessary, use the surgical marker and surgical planning scale to define the insertion points for the trocar cannula.
- 4 - Remove the trocar cannula cap and insert the sharp obturator tip to initiate the puncture.
- 5 - Once the outer part of the joint is reached, replace the sharp obturator tip with the blunt obturator tip and advance to the desired target location.
- 6 - Use the instruments included in the easyATM Cannula Kit according to the procedure to be performed, as well as the surgeon's technique and assessment.
- 7 - Perform the procedure according to the guidance of a qualified

professional.

8 - Remove the trocar cannula.

9 - After completing the procedure, deactivate and dispose of all materials in the appropriate location, following hospital regulations for the disposal of sharp and contaminated materials.

CONTRAINDICATIONS

- Non-arthroscopic surgical procedures.
- Patients for whom an arthroscopic procedure is contraindicated for any reason.
- Patients with allergies to contrast agents, antibiotics, or local anesthetics used in the procedure.
- Pregnant women.
- Patients with neurological disorders; findings during preoperative assessment are considered a relative contraindication.
- Coagulation disorders, for example, during anticoagulant therapy, are considered a relative contraindication.
- Known hypersensitivity to the materials used.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will guide proper and safe usage and protect against potential risks.
- Do not use the product past its expiration date; always verify that the sterilization indication is correct.
- The product should only be used by qualified medical professionals who possess the knowledge and experience necessary for proper use, including handling of sterile products and, importantly, the ability to maintain a completely sterile environment in the operating room, where careful attention to every detail is essential to avoid contamination risks.
- Check the integrity of the packaging before using the product; if the packaging is compromised, DO NOT USE, as sterility cannot be guaranteed.
- This product is for single use only and sterilized with ethylene oxide (ETO).
- Reprocessing of this product is prohibited.
- No modifications to this product are allowed.
- Disposal of the product after use must follow applicable regulations for the disposal of sharp and contaminated materials, in accordance with the hospital's Health Care Waste Management Plan and the guidance of the Hospital Infection Control Committee, as provided by RDC 222/2018.

STORAGE AND TRANSPORTATION CONDITIONS

To maintain the integrity of the product, it must:

Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.

Have packaging that is not punctured or torn and keep it away from sharp or pointed instruments.

Not be placed under weight or heavy loads, as this may damage the product's technical characteristics.

Not be placed near or together with different product categories, as this may affect its technical characteristics.

Be packaged and transported following good practices for the transport of sterilized products.

Be sterilized with ethylene oxide (ETO) and not have damaged packaging prior to use.

Be transported in its original packaging.

Storage temperature: Between 15°C and 30°C.

Storage humidity: 30% to 80%.

ANVISA Registration No.: 80195529014

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br