
Kit Cânula Cifoplastia t.Kypho

Reg. ANVISA: 80195520030.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) – is in accordance with RDC No. 751/2022. A PDF reader is required, and the program can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the product purchased.

To obtain the instruction for use in printed format, contact us via email at sac@bramsys.com.br. There is no additional cost for sending the product instruction for use in printed or digital format.

This instruction for use is compatible with all commercial models of the Kit Cannula Kyphoplasty t.Kypho family.

USE INSTRUCTIONS

KIT CÂNULA CIFOPLASTIA T.KYPHO

INDICATION FOR USE

The Kit Cannula Kyphoplasty t.Kypho is indicated for minimally invasive kyphoplasty procedures in the treatment of spinal fractures or damage.

PURPOSE

The Kit Cannula Kyphoplasty t.Kypho is intended for minimally invasive percutaneous access to the fractured or damaged vertebra through the working cannulas. Once the vertebra is accessed, a balloon is used to restore the vertebral body height, creating a cavity that is subsequently filled by the application of bone cement.

COMMERCIAL MODELS AND CONTENTS

Code	Working Cannula	Balloon	Composition
TKU11G-10	11G - 120mm	10mm	Unilateral
TKU11G-15	11G - 120mm	15mm	Unilateral
TKU11G-20	11G - 120mm	20mm	Unilateral

The commercial models of the Kit Cannula Kyphoplasty t.Kypho consist of the following items:

Items	Technical Specifications
(1) Working Cannula	Diameter: 11G Length: 120mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(4) Winged Syringe	Administration Capacity: 2.5 mL Chemical Composition: Polypropylene, Polystyrene, and Rubber
(1) Drill	Diameter: 13G Length: 187 mm Chemical Composition: 440 Stainless Steel and ABS Copolymer
(1) Inflator Device	Operating Range: -0.4 ATM to 30 ATM (1 ATM = 1 BAR = 14.7 PSI) Accuracy: $\pm 3.0\%$ of full scale Administration Capacity: 20 mL Administered Fluid: 0.45 mL (± 0.07 mL) per 360° clockwise rotation of the syringe plunger handle Chemical Composition: Polycarbonate, Polypropylene, and Silicone
(1) 3-Way Luer Lock Adapter	Luer Lock Connection Chemical Composition: Polycarbonate
(1) Balloon Catheter	Length: 10mm (TKU11G-10), 15mm (TKU11G-15) and 20mm (TKU11G-20) Capacity: < 400psi Diameter: 11G Maximum volume: 10mm model (3mL); 15mm model (4mL); 20mm model (6mL) Chemical composition: 304 stainless steel, polycarbonate, PTFE and silicone
(3) Bone Cement Delivery Cannula	Diameter: 13G Length: 187 mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(1) Bone Cement Mixer	Administration Capacity: 40 g Chemical Composition: Polypropylene, Polystyrene, and 304 Stainless Steel

The commercial models are supplied in the following presentation:

- **Primary Packaging:** PET (polyethylene terephthalate) blister sealed with Tyvek paper, permeable to ethylene oxide for sterilization, containing all items that make up the kit, as described above, along with 5 traceability labels.
- **Secondary Packaging:** Individual packaging box made of duplex paper, containing the primary packaging.
- **Tertiary Packaging:** Transport box.

Each package contains only 1 (one) commercial model of the Kit Cannula Kyphoplasty t.Kypho.

PRODUCT INSTRUCTIONS FOR USE

- 1- Identify the site for vertebral access puncture.
- 2- Insert the Working Cannula to the desired target.
- 3- Remove the mandrel from the Working Cannula.
- 4- Insert the Drill into the Working Cannula.
- 5- Perform circular movements to drill the vertebra, preparing it for Balloon

Catheter insertion.

- 6- Remove the Drill from the Working Cannula.
- 7- Collect 10 mL of contrast using a syringe.
- 8- Connect the syringe to the 3-Way Adapter.
- 9- Set the 3-Way Adapter to inject the contrast.
- 10- Connect the 3-Way Adapter to the Inflator Device.
- 11- Unlock the plunger of the Inflator Device.
- 12- Inject the contrast into the Inflator Device.
- 13- Disconnect the syringe from the 3-Way Adapter.
- 14- Set the 3-Way Adapter to use the Balloon Catheter.
- 15- Remove air between the Inflator Device plunger and the 3-Way Adapter outlet.
- 16- Lock the plunger of the Inflator Device.
- 17- Insert the Balloon Catheter into the Working Cannula until the positioning mark is completely inside.
- 18- Adjust the position if necessary until the radiopaque markers of the Balloon Catheter are outside the Working Cannula.
- 19- Connect the 3-Way Adapter to the Balloon Catheter access.
- 20- Rotate the Inflator Device plunger clockwise to inject contrast into the Balloon Catheter and inflate it, creating the vertebral body cavity.

ATTENTION: Maximum pressure and volume are described in the technical specifications of the Balloon Catheter.

- 21- Prepare the bone cement using the Bone Cement Mixer.
- 22- Attach the Winged Syringe to the Bone Cement Mixer.
- 23- Collect the bone cement in the Winged Syringe.
- 24- Disconnect the Winged Syringe from the Bone Cement Mixer.
- 25- Connect the Winged Syringe to the Bone Cement Delivery Cannula.
- 26- Inject the bone cement into the Bone Cement Delivery Cannula until fully filled.
- 27- Disconnect the Winged Syringe from the Bone Cement Delivery Cannula.
- 28- Rotate the Inflator Device plunger counterclockwise until all contrast is removed from the Balloon Catheter.
- 29- Unlock the Inflator Device plunger and pull to create negative pressure in the Balloon Catheter, ensuring it is completely empty.
- 30- Lock the Inflator Device plunger.
- 31- Remove the Balloon Catheter from the Working Cannula.
- 32- Insert the Bone Cement Delivery Cannula into the Working Cannula.
- 33- Insert the mandrel into the Bone Cement Delivery Cannula to transfer cement to the vertebra.
- 34- Remove the Bone Cement Delivery Cannula from the Working Cannula.
- 35- Remove the Working Cannula.
- 36- Cover the incision with sterile gauze and suture if necessary.
- 37- Dispose of all materials used in the procedure in an appropriate location, following hospital regulations for disposal of sharps and contaminated materials.

CONTRAINDICATIONS

- Infection or hematoma near the chosen puncture site.
- Pregnancy.
- Neurological alterations identified during the preoperative assessment (considered a relative contraindication).
- Coagulation disorders, for example, during anticoagulant therapy (considered a relative contraindication).
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to contrast agents, antibiotics, and/or local anesthetics used in the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedural performance.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will guide correct and safe use and help protect against potential risks.
- Do not use the product past its expiration date and always verify that the sterilization indicator is correct.
- The product should only be used by qualified medical professionals who possess the necessary knowledge and experience for proper product use, handling of sterile products, and, most importantly, understanding the importance of maintaining a completely sterile environment in the operating room, where attention to detail is critical to avoid contamination risks.
- Check the integrity of the packaging before use. If the packaging is compromised, DO NOT USE, as sterility cannot be guaranteed.
- This product is for single use only and sterilized by Ethylene Oxide (ETO).
- Reprocessing of this product is strictly prohibited.
- No modifications to the product are allowed.
- Disposal of the product after use must comply with regulations for the

disposal of sharps and contaminated materials, according to the hospital's Healthcare Waste Management Plan, in accordance with the guidance of the Hospital Infection Control Committee and as established in RDC 222/2018.

STORAGE AND TRANSPORTATION CONDITIONS

To maintain the integrity of the product, it must:

- Always be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- Remain in packaging that is not punctured or torn and kept away from sharp or pointed instruments.
- Not be placed under weight or heavy loads, as this may damage the product's technical characteristics.
- Not be stored or transported together with different product categories, as this may compromise its technical characteristics.
- Be packaged and transported following good practices for sterilized products.
- This product is sterilized with Ethylene Oxide (ETO) and the packaging must remain intact until use.
- Transport the product in its original packaging.
- Storage temperature: Between 15°C and 30°C.
- Storage humidity: 30% to 80%.

ANVISA Registration No.: 80195520030

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the National Health Surveillance Agency (ANVISA) under No. UY44WX94W79M, guarantees the quality and sterile condition of this product for **24 months** from the sterilization date, except in cases of misuse or failure to follow the instructions provided in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br