
Kit Cãnula Axon Release

Reg. ANVISA: 80195529018.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instructions for Use document in a non-printed format – digital medium (www.bramsys.com.br/instrucoes-de-uso) complies with RDC No. 751/2022. The use of a PDF reader is required. To obtain the software, download it free of charge at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the Instructions for Use corresponding to that date or an earlier one. Do not use Instructions for Use with a revision/date (Revision: XX - day/month/year) different from the purchased product.

To obtain the Instructions for Use in printed format, contact us via e-mail at sac@bramsys.com.br. There is no additional cost for sending the product Instructions for Use in printed or digital format.

These Instructions for Use are compatible with all commercial models of the Axon Release Cannula Kit family.

USE INSTRUCTIONS

AXON RELEASE CANNULA KIT

INDICATION FOR USE / PURPOSE

Purpose: The Axon Release Cannula Kit is a set developed to map, locate, and assist in the anatomical release of peripheral nerve structures, according to the adopted surgical technique, in open surgical procedures.

Its primary purpose is to enable the planning, mapping, and precise localization of the neural target through impedancemetry, sensory and motor neurostimulation, as well as the possibility of applying controlled thermal modulation, according to medical decision.

Indications for Use: The Axon Release Cannula Kit is indicated for use in open surgeries involving:

- Surgical planning, mapping, and localization of peripheral nerves;
- Controlled neuromodulation, when applicable;
- Support for performing neurotomy and neurolysis, according to medical judgment.

The kit components may be used together or separately, according to the adopted surgical technique and the patient's clinical needs.

COMMERCIAL MODELS AND CONTENTS

Code	Length	Active Tip
AXR-5515	55mm	1.5mm
AXR-7515	75mm	1.5mm
AXR-11515	115mm	1.5mm
AXR-5525	55mm	2.5mm
AXR-7525	75mm	2.5mm
AXR-11525	115mm	2.5mm
AXR-5535	55mm	3.5mm
AXR-7535	75mm	3.5mm
AXR-11535	115mm	3.5mm
AXR-5545	55mm	4.5mm
AXR-7545	75mm	4.5mm
AXR-11545	115mm	4.5mm
AXR-5555	55mm	5.5mm
AXR-7555	75mm	5.5mm
AXR-11555	115mm	5.5mm

The commercial models of the Axon Release Cannula Kit are composed of the following items:

Items	Technical Specifications
(1) Probe with ergonomic handle	Diameter: 1.27 mm (18G); Length: 55 mm, 75 mm, or 115 mm; Active tip: 1.5 mm, 2.5 mm, 3.5 mm, 4.5 mm, or 5.5 mm; Cable length: 2.5 meters;
(1) Cannula	Diameter: 20G Length: 70 mm
(1) Surgical Planning Marker	Gentian violet skin marking ink
(1) Surgical Planning Scale	Length: 15 cm

The commercial models are supplied in the following presentation:

- Primary packaging: Surgical-grade pouch containing all items that compose the kit, as presented above, together with the 5 traceability labels.
- Secondary packaging: Individual packaging box made of duplex paperboard, containing the primary packaging.

Each package contains only 1 commercial model of the Axon Release Cannula Kit.

PRODUCT INSTRUCTIONS FOR USE

1. Perform the surgical planning using the anatomical marker and the planning scale to define the incision point and the pathway to the target nerve structure.

2. Connect the probe to the generator unit, when applicable, to perform: Impedancemetry to verify electrical contact; Sensory and motor neurostimulation for target mapping; Controlled thermal modulation, according to medical decision.

Note: If necessary, use the graduated cannula to inject conductive substance (saline solution).

3. Perform the surgical procedure with the support of the kit, according to the guidance of a qualified healthcare professional.

4. Remove the components used at the end of the procedure.

5. Decharacterize and dispose of all materials in an appropriate location, according to regulations for the disposal of sharps and contaminated materials.

CONTRAINDICATIONS

- Infection or hematoma near the procedure site;
- Pregnancy;
- Known preoperative neurological disorders are a relative contraindication;
- Coagulation disorders, for example during anticoagulant therapy, are a relative contraindication;
- Inability to obtain patient consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast agents, antibiotics, and/or local anesthetics used during the procedure;
- Any medical, post-surgical, congenital, anatomical, or psychological condition that increases the risk of reduced procedure performance.

ADVERSE EFFECTS

Transient local pain, hematoma, bleeding, edema, local inflammation, inadvertent injury to adjacent structures (blood vessels, muscles, tendons, or other nerves), infection at the surgical site, and inadequate healing.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Do not use the product beyond its expiration date. Always verify that the sterilization indicator is correct.
- The product must only be used by qualified medical professionals who possess the necessary knowledge and experience for proper use of the product, knowledge in handling sterile products, and, above all, awareness of the importance of maintaining a completely sterile environment in operating rooms, where care and attention to every detail are essential to avoid risks and contamination.
- The product must not be reused and/or reprocessed, as it is intended for SINGLE USE only, and must be destroyed and disposed of as HOSPITAL WASTE.
- The product must always be transported and stored protected from all types of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and must be kept away from cutting or sharp instruments.
- Do not store this product near products of different categories.
- Check the integrity of the packaging before using the product. In case of package damage or violation, DO NOT USE, as the sterility guarantee will be compromised.
- Reprocessing of this product is prohibited.
- This product is single-use and sterilized with Ethylene Oxide (ETO).

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

- Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and must be kept away from cutting or sharp instruments.
- Do not place under weight or stacked volumes, as this may damage the product's technical characteristics.
- Do not place near or together with products from different categories, as they may damage its technical characteristics.
- Be packaged and transported following good practices for the transportation of sterilized products.
- This product is sterilized with Ethylene Oxide (ETO) and the packaging must not be damaged before use.
- The product must be transported in its original packaging.
- Storage and transportation temperature: Between 15°C and 30°C.
- Storage and transportation humidity: 30% to 80%.

ANVISA Registration No.: 80195529018

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the Brazilian Health Regulatory Agency (ANVISA), guarantees for **24 months** from the sterilization date the quality and sterile condition of this product, except in cases of misuse or failure to comply with the guidelines contained in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br