
Kit Cãnula DiscGuide

Reg. ANVISA: 80195529019.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Read all instructions carefully before use. Follow all warnings and precautions provided in these instructions. Failure to comply with these instructions may result in complications.

The availability of these Instructions for Use in a non-printed format – digital media ([www.bramsys.com.br/instrucoes-de-uso])(<http://www.bramsys.com.br/instrucoes-de-uso>) – is in compliance with RDC No. 751/2022. A PDF reader is required. To obtain the software, download it free of charge at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the Instructions for Use corresponding to that date or an earlier version. Do not use Instructions for Use with a revision/date (Revision: XX - day/month/year) different from that of the purchased product.

To obtain the Instructions for Use in printed format, please contact us via e-mail at [sac@bramsys.com.br](<mailto:sac@bramsys.com.br>). There is no additional cost for providing the product's Instructions for Use in either printed or digital format.

These Instructions for Use are applicable to all DCG-165195 models.

USE INSTRUCTIONS

DISCGUIDE CANNULA KIT

INDICATION FOR USE / PURPOSE

The DiscGuide Cannula Kit was developed to provide the instruments required for intradiscal access and decompression procedures of the spinal column, intended for the treatment of intervertebral disc herniation. The kit enables safe and precise access to the affected disc, allowing for the controlled removal of disc material and the reduction of pressure on adjacent neural structures. It can be used in both percutaneous and endoscopic techniques.

COMMERCIAL MODELS AND CONTENTS

Code	Description
DCG-165195	DiscGuide Cannula Kit

The commercial models of the Axon Release Cannula Kit are composed of the following items:

Items	Technical Specifications
(1) 18G Guide Wire	Diameter: 18 Gauge Length: 330 mm
(1) 20G Guide Wire	Diameter: 20 Gauge Length: 330mm
(1) Introducer Cannula	Diameter: 16 Gauge Length: 200mm
(1) Trocar	Diameter: 10 Gauge Length: 165mm
(1) Faceted Cannula	Diameter: 10 Gauge Length: 165mm
(1) Trepine	Diameter: 12 Gauge Length: 195mm
(1) Obturator	Diameter: 12 Gauge Length: 195 mm
(1) Limiting Stop	Internal Diameter: Compatible with 10 Gauge
(1) Surgical Planning Marker	Gentian violet skin-marking ink
(1) Surgical Planning Scale	Length: 15 cm

The commercial models are supplied in the following configuration:

- **Primary Packaging:** Blister pack containing all components of the kit, as described above, together with the 5 traceability labels.
- **Secondary Packaging:** Individual packaging box made of duplex paperboard, containing the primary packaging.

Each package contains only one commercial model of the DiscGuide Cannula Kit.

PRODUCT INSTRUCTIONS FOR USE

1. Define the puncture site for access using fluoroscopy equipment.
2. Insert the introducer cannula to the desired target.
3. Remove the stylet from the introducer cannula.
4. Inject contrast medium to visualize the disc.
5. Insert the guide wire to the desired target.
6. Carefully remove the introducer cannula while maintaining the guide wire at the desired target.
7. Select either the trocar or the faceted cannula to perform the procedure.
8. Attach the obturator to the selected option.
9. Position the limiting stop on the selected option.
10. Introduce the assembly consisting of the selected option, obturator, and limiting stop into the patient over the guide wire.
11. Position the limiting stop on the selected option until it rests against the patient's skin.
12. Lock the limiting stop to stabilize the selected option within the patient.
13. Remove the obturator from the selected option.
14. Insert the trephine into the selected option.
15. Press the trephine against the disc wall while performing rotational movements until the disc is penetrated.
16. Remove the trephine from the selected option.

17. Perform the procedure on the disc through the selected option.
18. After completing the procedure, remove the instruments from the selected option.
19. Remove the selected option from the patient.
20. Upon completion of the procedure, the material must be destroyed and disposed of in an appropriate location, in accordance with applicable regulations for the disposal of sharps and contaminated materials.

CONTRAINDICATIONS

- Infection or hematoma in the vicinity of the procedure site;
- Pregnancy;
- Pre-existing neurological disorders identified during the preoperative assessment are considered a relative contraindication;
- Coagulation disorders, such as those occurring during anticoagulant therapy, are considered a relative contraindication;
- Inability to obtain the patient's informed consent;
- Known hypersensitivity to the materials used;
- Allergy to contrast media, antibiotics, and/or local anesthetics used during the procedure;
- Any medical, post-surgical, congenital, anatomical, or psychological condition that may increase the risk of reduced procedural performance.

ADVERSE EFFECTS

Pain at the procedure site, bleeding, hematoma, local or deep infection (such as discitis), and injury to adjacent structures.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Do not use the product beyond its expiration date. Always verify that the sterilization indicator is valid and correct.
- The product should only be used by qualified medical professionals who possess the knowledge and experience required for its proper use, including expertise in handling sterile products and, most importantly, an understanding of the importance of maintaining a completely sterile environment in operating rooms, where attention to detail and strict precautions are essential to prevent risks and contamination.
- The product must not be reused or reprocessed. It is intended for SINGLE USE ONLY and must be destroyed and disposed of as HOSPITAL WASTE.
- The product must be transported and stored in a manner that protects it from all forms of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.
- Do not store this product near products of different categories.
- Verify the integrity of the packaging before use. If the packaging has been compromised, DO NOT USE the product, as its sterility can no longer be guaranteed.
- Reprocessing of this product is strictly prohibited.
- This product is for single use only and has been sterilized using Ethylene Oxide (ETO).

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

- Be stored and transported at all times protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- Be kept in packaging that is not punctured or torn and away from sharp or pointed instruments.
- Not be placed under weight or heavy loads, as this may damage the product's technical characteristics.
- Not be stored near or together with products from different categories, as they may adversely affect its technical characteristics.
- Be packaged and transported in accordance with good practices for the handling and transportation of sterile products.
- Be kept in its original packaging, as the product is sterilized using Ethylene Oxide (ETO) and the packaging must not be damaged before use.
- Be transported in its original packaging.
- Storage and transportation temperature: Between 15°C and 30°C.
- Storage and transportation humidity: 30% to 80%.

ANVISA Registration No.: 80195529019

WARRANTY CERTIFICATE

BRAMSYS INDÚSTRIA E COMÉRCIO LTDA, duly authorized by the Brazilian Health Regulatory Agency (ANVISA), guarantees for **24 months** from the sterilization date the quality and sterile condition of this product, except in cases of misuse or failure to comply with the guidelines contained in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br