
# Kit Cãnula Low Vert

Reg. ANVISA: 80195529012.

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E-mail: [sac@bramsys.com.br](mailto:sac@bramsys.com.br)

Sitio web: [www.bramsys.com.br](http://www.bramsys.com.br)

**ATTENTION:** Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium ([www.bramsys.com.br/instrucoes-de-uso](http://www.bramsys.com.br/instrucoes-de-uso)) – is in accordance with RDC No. 751/2022. A PDF reader is required, and the program can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the product purchased.

To obtain the instruction for use in printed format, contact us via email at [sac@bramsys.com.br](mailto:sac@bramsys.com.br). There is no additional cost for sending the product instruction for use in printed or digital format.

These instructions for use are compatible with all commercial models of the Kit Cânula Low Vert family.

# USE INSTRUCTIONS

## LOW VERT CANNULA KIT

### INDICATION FOR USE/ PURPOSE

The Low Vert Cannula Kit is indicated for access in minimally invasive vertebroplasty procedures for the treatment of vertebral fractures.

### COMMERCIAL MODELS AND CONTENTS

Code	Description
LVU12011	Low Vert Cannula Kit 11G 120mm
LVB12011	Low Vert Cannula Kit 11G 120mm/2

The Low Vert Cannula Kit presents the following technical specifications for each respective commercial model:

Code	Beveled Cannula			Drill			Application Cannula		
	Length	Ø	Qty.	Length	Ø	Qty.	Length	Ø	Qty.
LVU12011	120	11	1	187	13	1	187	11	3
LVB12011	120	11	2	187	1	1	187	11	6

The commercial models are supplied in the following presentation:

- Primary packaging: PET (polyethylene terephthalate) blister sealed with Tyvek paper permeable to ethylene oxide for sterilization, containing all kit components, according to the presentation described above, together with the 5 traceability labels.
- Secondary packaging: Individual packaging box, made of duplex paper, containing the primary packaging.

Each package contains only 1 commercial model of the Low Vert Cannula Kit.

### PRODUCT USE INSTRUCTIONS

1. Ensure that the packaging is not damaged.
  2. Remove all items from the sterile packaging in a clean environment using aseptic techniques.
  3. Ensure that all items are intact.
  - 4 - Position the patient according to the surgical procedure to be performed.
  5. Define the puncture site for vertebral access.
  6. Insert the expandable working cannula until the desired target is reached.
  7. Remove the stylet from the expandable working cannula.
  8. Insert the Drill into the expandable working cannula.
  9. Perform clockwise rotational movements with the Drill until bone resistance is no longer felt.
  10. Remove the Drill from the expandable working cannula.
  11. Remove the stylet from the cannulas for bone cement application and fill them.
  12. Insert the bone cement application cannula into the expandable working cannula.
  13. Insert the stylet into the bone cement application cannula to transfer the cement into the vertebra.
- NOTE:** Repeat steps 11 to 13 until the vertebra is filled as desired.
14. Remove the bone cement application cannula from the expandable working cannula.
  15. Remove the expandable working cannula.
  16. Dispose of all material used in the procedure in an appropriate location.

### CONTRAINDICATIONS

The use is not recommended in patients presenting a serious hemorrhagic condition, bleeding disorders, or who are undergoing anticoagulant treatment. Medical opinion is essential when considering the performance of surgery in patients under these conditions.

### ADVERSE EFFECTS

Adverse effects include, but are not limited to: discomfort, hematoma formation, infection, bleeding, damage to adjacent tissue, and blood vessel perforation.

## **PRECAUTIONS, RESTRICTIONS AND WARNINGS**

- Avoid using the product if the packaging is damaged;
- Do not use the product if the expiration date has passed;
- The product is intended for single use only and must be discarded after use;
- Reuse is strictly prohibited;
- Reuse or resterilization may compromise the structural integrity of the product and/or cause damage that may result in injury or illness to the patient;
- Reuse or resterilization may also represent a risk of product contamination and/or cause infection in the patient or cross-infection, including, but not limited to, the transmission of infectious diseases between patients;
- It is advisable to avoid using the product if particularly hazardous maneuvers are required, when it is subjected to torsional or bending forces that may compromise its integrity and proper functioning;
- Post-surgical care may vary depending on the technique used and the physiological conditions of each patient;
- Careful clinical judgment is required when considering surgery in patients undergoing anticoagulant therapy or with bleeding disorders;
- This product must be used by a physician who is fully familiar with the indications, contraindications, limitations, typical outcomes, and possible adverse effects related to the use of vertebral puncture.

## **STORAGE AND TRANSPORTATION CONDITIONS**

To maintain the integrity of the product, it must:

- Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and must be kept away from cutting or sharp instruments.
- Do not place under weight or volumes, as they may cause damage to the technical characteristics of the product.
- Do not place near or together with different categories, as they may cause damage to its technical characteristics.
- They must be stored and transported following good practices for the transportation of sterilized products.
- This product is sterilized by ethylene oxide (ETO) and the packaging must not be damaged before use.
- The product must be transported in its original packaging.
- Storage temperature: Between 15 °C and 30 °C.

**ANVISA Registration No. 80195529012**

### **WARRANTY CERTIFICATE**

**BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, properly authorized by the National Agency of Sanitary Surveillance, guarantees for **24 MONTHS**, starting from the sterilization date, the quality and the sterile condition of this product, except in the cases of bad usage or noncompliance of the orientations written in this document.

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