
Kit Cãnula b.Clinic Lite

Reg. ANVISA: 80195529020.

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E-mail: sac@bramsys.com.br

Sitio web: www.bramsys.com.br

ATTENTION: Read all instructions carefully before use. Follow all warnings and precautions described in these instructions. Failure to observe these points may result in complications.

The availability of this Instructions for Use document in a non-printed format – digital media ([\[www.bramsys.com.br/instrucoes-de-uso\]](http://www.bramsys.com.br/instrucoes-de-uso)(<http://www.bramsys.com.br/instrucoes-de-uso>)) is in compliance with RDC No. 751/2022. The use of a PDF reader is required. To obtain the software, download it free of charge at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date and download the version of the Instructions for Use corresponding to that date or an earlier version. Do not use Instructions for Use with a revision/date (Revision: XX - day/month/year) different from that of the purchased product.

To obtain the Instructions for Use in printed format, please contact us via e-mail at [\[sac@bramsys.com.br\]](mailto:sac@bramsys.com.br)(<mailto:sac@bramsys.com.br>). There is no additional cost for receiving the product Instructions for Use in either printed or digital format.

These Instructions for Use are compatible with all commercial models of the b.Clinic Lite Cannula Kit family.

USE INSTRUCTIONS

B.CLINIC LITE CANNULA KIT

INDICATION FOR USE / PURPOSE

The b.Clinic Lite Cannula Kit is intended for interventional procedures involving percutaneous access and precise positioning in specific anatomical structures, assisting physicians in performing therapeutic and diagnostic techniques related to pain management, without temperature measurement. The device is intended to enable access and the performance of interventional stimulation and neuromodulation procedures, allowing precise positioning in target regions for medical pain management techniques.

COMMERCIAL MODELS AND CONTENTS

Code	Ø	Length	Active Tip
BCL-18050-5	18G	50mm	5mm
BCL-18100-5	18G	100mm	5mm
BCL-18150-5	18G	150mm	5mm
BCL-18100-10	18G	100mm	10mm
BCL-18150-10	18G	150mm	10mm
BCL-18100-15	18G	100mm	15mm
BCL-18150-15	18G	150mm	15mm
BCL-20050-5	20G	50mm	5mm
BCL-20100-5	20G	100mm	5mm
BCL-20150-5	20G	150mm	5mm
BCL-20100-10	20G	100mm	10mm
BCL-20150-10	20G	150mm	10mm
BCL-20100-15	20G	100mm	15mm
BCL-20150-15	20G	150mm	15mm
BCL-21050-5	21G	50mm	5mm
BCL-21100-5	21G	100mm	5mm
BCL-21150-5	21G	150mm	5mm
BCL-21100-10	21G	100mm	10mm
BCL-21150-10	21G	150mm	10mm
BCL-21100-15	21G	100mm	15mm
BCL-21150-15	21G	150mm	15mm

The commercial models of the b.Clinic Lite Cannula Kit product are composed of the following items:

Items	Technical Specifications
(1) b.Clinic Lite Cannula	Diameter: 18G, 20G, and 21G Working Length: 50 mm, 100 mm, 150 mm Active Tip: 5 mm, 10 mm, 15 mm Shape: Straight Bevel: Sharp Echogenic Tip Injection tube for substance delivery with Luer Lock connection Coating: Transparent polyester Centimeter markings Cable Length: 60 cm
(1) Surgical Planning Scale	Chemical composition: PVC
(1) Surgical Marker	Chemical composition: Polypropylene, polyester, 3% gentian violet, and 30% isopropanol alcohol

The commercial models are supplied in the following configuration:

- Primary packaging: a surgical-grade pouch containing all items that make up the kit, as described above, together with the 5 traceability labels.
- Secondary packaging: an individual packaging box made of duplex paperboard, containing the primary packaging.

Each package contains only 1 (one) commercial model of the b.Clinic Lite Cannula Kit.

PRODUCT INSTRUCTIONS FOR USE

1. Remove the protective tube.
2. Position the cannula on the patient.

3. Connect the cannula to the generator unit.
4. Perform the procedure according to the guidance of a qualified professional.

NOTE: If necessary, use the injection tube to administer substances.

5. Remove the cannula.

6. After completing the procedure, all materials must be decharacterized and disposed of in an appropriate location, following the applicable regulations for the disposal of sharp and contaminated materials.

CONTRAINDICATIONS

- Infection or hematoma near the procedure site.
- Pregnancy.
- Neurological changes identified in the preoperative evaluation are a relative contraindication.
- Coagulation disorders, for example during anticoagulant therapy, are a relative contraindication.
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to contrast agents, antibiotics, and/or local anesthetics used in the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedural performance.
- The use of general anesthesia is not recommended. During the procedure, a sensory response from the patient may be necessary; therefore, it is recommended that the treatment be performed under local anesthesia.

ADVERSE EFFECTS

Pain at the application site; hematoma; bleeding; ecchymosis; local edema; inflammation; infection at the puncture site; abscess.

PRECAUTIONS, RESTRICTIONS, AND WARNINGS

- Before using the product, the user must carefully read the instructions for use, which will assist in its correct and safe use and help protect against possible risks.
- Do not use the product past its expiration date; always check that the sterilization indicator is correct.
- The product must only be used by qualified medical professionals who possess the knowledge and experience required for its proper use, including handling sterile products and, most importantly, maintaining a fully sterile environment in operating rooms, where care and attention to every detail are essential to avoid contamination risks.
- Check the integrity of the packaging before using the product. In case of violation, **DO NOT USE** it, as sterility can no longer be guaranteed.
- This product is for single use only and is sterilized with Ethylene Oxide (ETO).
- Reprocessing this product is prohibited; it must be destroyed and discarded as **HOSPITAL WASTE**.
- No modification to this product is permitted.
- Product disposal after use must follow applicable regulations for the disposal of sharp and contaminated materials, in accordance with the hospital's Healthcare Waste Management Plan and the guidance of the Hospital Infection Control Committee, as established in RDC 222/2018.

STORAGE AND TRANSPORT CONDITIONS

To maintain the integrity of the product, it must:

- Be stored and transported protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and should be kept away from sharp or pointed instruments.
- Not be placed under heavy weights or volumes, as these may damage the product's technical characteristics.
- Not be placed near or together with products of different categories, as they may damage its technical characteristics.
- Be stored and transported following good practices for handling sterilized products.
- This product is sterilized by ethylene oxide (ETO) and must not have its packaging damaged before use.
- Be transported in its original packaging.
- Storage temperature: Between 15°C and 30°C.
- Storage humidity: 30% to 80%.

ANVISA Registration No.: 80195529020

WARRANTY CERTIFICATE

The company **BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, duly authorized by the National Health Surveillance Agency (ANVISA), guarantees for a period of **24 months** from the date of sterilization the quality and sterile condition of this product, except in cases of misuse or failure to comply with the instructions contained in this document.

CUSTOMER SERVICE (CS)

Phone: +55 11 5013 0010

E-mail: sac@bramsys.com.br

Site: www.bramsys.com.br