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**E-mail:** sac@bramsys.com.br

**Sitio web:** www.bramsys.com.br

**ATTENTION:** Carefully read all instructions before use. Follow all warnings and precautions mentioned in these instructions. Failure to observe these points may result in complications.

The availability of this Instruction for Use in a non-printed format – digital medium ([www.bramsys.com.br/instrucoes-de-uso](http://www.bramsys.com.br/instrucoes-de-uso)) – is in accordance with RDC No. 751/2022. A PDF reader is required, and the program can be downloaded for free at: <https://www.adobe.com/br/acrobat/pdf-reader.html>.

Check the product label for its manufacturing date, and download the version of the instruction for use corresponding to that date or an earlier one. Do not use an instruction for use with a revision/date (Revision: XX - day/month/year) different from the product purchased.

To obtain the instruction for use in printed format, contact us via email at [sac@bramsys.com.br](mailto:sac@bramsys.com.br). There is no additional cost for sending the product instruction for use in printed or digital format.

This instruction for use is compatible with all commercial models of the Kit Cannula Kyphoplasty t.Kypho family.

# USE INSTRUCTIONS

## KIT CÂNULA CIFOPLASTIA T.KYPHO

### INDICATION FOR USE

The t.Kypho Kyphoplasty Cannula Kit is indicated for minimally invasive kyphoplasty procedures in the treatment of fractures or damage to the spine. The product is intended for minimally invasive percutaneous access to the fractured or damaged vertebra through working cannulas. After accessing the vertebra, a balloon is used to restore the height of the vertebral body, creating a cavity that must be filled with bone cement.

### COMMERCIAL MODELS AND CONTENTS

Code	Working Cannula	Balloon	Composition
TKU11G-10	11G - 120mm	10mm	Unilateral
TKU11G-15	11G - 120mm	15mm	Unilateral
TKU11G-20	11G - 120mm	20mm	Unilateral

The commercial models of the Kit Cannula Kyphoplasty t.Kypho consist of the following items:

Items	Technical Specifications
(1) Working Cannula	Diameter: 11G Length: 120mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(4) Winged Syringe	Administration Capacity: 2.5 mL Chemical Composition: Polypropylene, Polystyrene, and Rubber
(1) Drill	Diameter: 13G Length: 187 mm Chemical Composition: 440 Stainless Steel and ABS Copolymer
(1) Inflator Device	Operating Range: -0.4 ATM to 30 ATM (1 ATM = 1 BAR = 14.7 PSI) Accuracy: $\pm 3.0\%$ of full scale Administration Capacity: 20 mL Administered Fluid: 0.45 mL ( $\pm 0.07$ mL) per 360° clockwise rotation of the syringe plunger handle Chemical Composition: Polycarbonate, Polypropylene, and Silicone
(1) 3-Way Luer Lock Adapter	Luer Lock Connection Chemical Composition: Polycarbonate
(1) Balloon Catheter	Length: 10mm (TKU11G-10), 15mm (TKU11G 15) and 20mm (TKU11G-20) Capacity: < 400psi Diameter: 11G Maximum volume: 10mm model (3mL); 15mm model (4mL); 20mm model (6mL) Chemical composition: 304 stainless steel, polycarbonate, PTFE and silicone
(3) Bone Cement Delivery Cannula	Diameter: 13G Length: 187 mm Chemical Composition: 304 Stainless Steel and ABS Copolymer
(1) Bone Cement Mixer	Administration Capacity: 40 g Chemical Composition: Polypropylene, Polystyrene, and 304 Stainless Steel

The commercial models are supplied in the following presentation:

- **Primary Packaging:** PET (polyethylene terephthalate) blister sealed with Tyvek paper, permeable to ethylene oxide for sterilization, containing all items that make up the kit, as described above, along with 5 traceability labels.
- **Secondary Packaging:** Individual packaging box made of duplex paper, containing the primary packaging.
- **Tertiary Packaging:** Transport box.

Each package contains only 1 (one) commercial model of the Kit Cannula Kyphoplasty t.Kypho.

### PRODUCT INSTRUCTIONS FOR USE

- 1- Identify the site for vertebral access puncture.
- 2- Insert the Working Cannula to the desired target.
- 3- Remove the mandrel from the Working Cannula.
- 4- Insert the Drill into the Working Cannula.
- 5- Perform circular movements to drill the vertebra, preparing it for Balloon Catheter insertion.
- 6- Remove the Drill from the Working Cannula.

- 7- Collect 10 mL of contrast using a syringe.
- 8- Connect the syringe to the 3-Way Adapter.
- 9- Set the 3-Way Adapter to inject the contrast.
- 10- Connect the 3-Way Adapter to the Inflator Device.
- 11- Unlock the plunger of the Inflator Device.
- 12- Inject the contrast into the Inflator Device.
- 13- Disconnect the syringe from the 3-Way Adapter.
- 14- Set the 3-Way Adapter to use the Balloon Catheter.
- 15- Remove air between the Inflator Device plunger and the 3-Way Adapter outlet.
- 16- Lock the plunger of the Inflator Device.
- 17- Insert the Balloon Catheter into the Working Cannula until the positioning mark is completely inside.
- 18- Adjust the position if necessary until the radiopaque markers of the Balloon Catheter are outside the Working Cannula.
- 19- Connect the 3-Way Adapter to the Balloon Catheter access.
- 20- Rotate the Inflator Device plunger clockwise to inject contrast into the Balloon Catheter and inflate it, creating the vertebral body cavity.

**ATTENTION: Maximum pressure and volume are described in the technical specifications of the Balloon Catheter.**

- 21- Prepare the bone cement using the Bone Cement Mixer.
- 22- Attach the Winged Syringe to the Bone Cement Mixer.
- 23- Collect the bone cement in the Winged Syringe.
- 24- Disconnect the Winged Syringe from the Bone Cement Mixer.
- 25- Connect the Winged Syringe to the Bone Cement Delivery Cannula.
- 26- Inject the bone cement into the Bone Cement Delivery Cannula until fully filled.
- 27- Disconnect the Winged Syringe from the Bone Cement Delivery Cannula.
- 28- Rotate the Inflator Device plunger counterclockwise until all contrast is removed from the Balloon Catheter.
- 29- Unlock the Inflator Device plunger and pull to create negative pressure in the Balloon Catheter, ensuring it is completely empty.
- 30- Lock the Inflator Device plunger.
- 31- Remove the Balloon Catheter from the Working Cannula.
- 32- Insert the Bone Cement Delivery Cannula into the Working Cannula.
- 33- Insert the mandrel into the Bone Cement Delivery Cannula to transfer cement to the vertebra.
- 34- Remove the Bone Cement Delivery Cannula from the Working Cannula.
- 35- Remove the Working Cannula.
- 36- Cover the incision with sterile gauze and suture if necessary.
- 37- Dispose of all materials used in the procedure in an appropriate location, following hospital regulations for disposal of sharps and contaminated materials.

#### **CONTRAINDICATIONS**

- Infection or hematoma near the chosen puncture site.
- Pregnancy.
- Neurological alterations identified during the preoperative assessment (considered a relative contraindication).
- Coagulation disorders, for example, during anticoagulant therapy (considered a relative contraindication).
- Inability to obtain patient consent.
- Known hypersensitivity to the materials used.
- Allergy to contrast agents, antibiotics, and/or local anesthetics used in the procedure.
- Any medical, postoperative, congenital, anatomical, or psychological condition that increases the risk of reduced procedural performance.

#### **PRECAUTIONS, RESTRICTIONS, AND WARNINGS**

- Before using the product, the user must carefully read the instructions for use, which will guide correct and safe use and help protect against potential risks.
- Do not use the product past its expiration date and always verify that the sterilization indicator is correct.
- The product should only be used by qualified medical professionals who possess the necessary knowledge and experience for proper product use, handling of sterile products, and, most importantly, understanding the importance of maintaining a completely sterile environment in the operating room, where attention to detail is critical to avoid contamination risks.
- The product must not be reused and/or reprocessed, being for SINGLE USE only; it must be destroyed and disposed of as HOSPITAL WASTE.
- The product must be transported and stored always protected from all types of moisture, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and must be kept away from cutting or sharp instruments. Do not store this product near products from different categories.

- Check the integrity of the packaging before using the product; if compromised, DO NOT USE, as sterility will be affected.
- Reprocessing of this product is prohibited.
- This product is for single use only, sterilized by Ethylene Oxide (ETO).
- Reprocessing of this product is prohibited.
- The product must be transported and stored always protected from excessive humidity, high temperatures, sunlight exposure, dirt, and contact with flammable, corrosive, toxic, or explosive substances, and kept away from cutting or sharp instruments.
- Do not place near or together with different categories, as they may cause damage to its technical characteristics and compromise its integrity/quality.
- No modifications are allowed to this product.
- Disposal of the product after use must follow applicable regulations for the disposal of sharp and contaminated materials adopted by the hospital.

#### **STORAGE AND TRANSPORTATION CONDITIONS**

To maintain the integrity of the product, it must:

- Be stored and transported always protected from sunlight, high temperatures, excessive humidity, dirt, and contact with flammable, corrosive, toxic, or explosive substances.
- The packaging must not be punctured or torn and must be kept away from cutting or sharp instruments.
- Do not place under weight or volumes, as they may cause damage to the technical characteristics of the product.
- Do not place near or together with different categories, as they may cause damage to its technical characteristics.
- They must be stored and transported following good practices for the transportation of sterilized products.
- This product is sterilized by ethylene oxide (ETO) and the packaging must not be damaged before use.
- The product must be transported in its original packaging.
- Storage and transportation temperature: Between 15 °C and 30 °C.
- Storage and transportation humidity: Between 30% and 80%.

**ANVISA Registration No.: 80195520030**

#### **WARRANTY CERTIFICATE**

**BRAMSYS INDÚSTRIA E COMÉRCIO LTDA**, duly authorized by the National Health Surveillance Agency (ANVISA), guarantees the quality and sterile condition of this product for **24 months** from the sterilization date, except in cases of misuse or failure to follow the instructions provided in this document.

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